



CERTIFICATE OF CONFORMANCE

Date: 24-Aug-20

This certification is provided as full assurance that the following product code and lot number were manufactured in accordance via a Quality System to ISO 9001 with prescribed procedures and specifications.

Product Code: 21109

Lot Number: 081620VX06

A handwritten signature in black ink, appearing to read "C. Inthavongsa".

Champagne Inthavongsa
Quality Supervisor
cinthavongsa@contecinc.com



Certificate of Analysis

Product Name: **PREempt RTU (US) 5G**
 EPA Reg. No.: 74559-1
 Product Code: **21109**
 Lot Number: L081620VX06
 Date of Testing: 08/13/2020
 Date of Expiry: 08/2023

Test Item	Specification	Result	Test Method
Appearance	Clear, colourless solution	Clear, colourless solution	Visual Inspection
Hydrogen Peroxide Content	0.50 - 0.55%	0.54	Virox No 1FP
pH	2.50 - 3.50	2.85	Virox No 2FP
Specific Gravity	1.005 - 1.015	1.01	Virox No 3FP

This is to certify that the above product met Finished Product Specifications at the time of manufacture

Prepared By: CMiller

Date: 08/13/2020