



HYPERSENSE® X11+

Latex Exam Gloves

PRODUCT DESCRIPTION

SPECIFICATIONS

PRODUCT:	Hypersense X11+		
GRADE:	Powder-Free Exam		
MATERIAL:	Latex		
THICKNESS	Finger:	14.0mil	0.36mm
	Palm:	12.0mil	0.31mm
	Cuff:	8.7mil	0.22mm
COLOR	Outer:	Blue	
	Inner:	Blue	
SURFACE	Outer:	Fully Textured	
	Inner:	Smooth	
LENGTH:	12.0"	305mm	
WEIGHT*:	19.9g		
FIT TYPE:	Ambidextrous		
CUFF TYPE:	Beaded		
WATERTIGHT AQL:	1.5		
ORIGIN COUNTRY:	China		

*Specification based on size large gloves

PHYSICAL PROPERTIES

	PRE-AGING	POST-AGING
TENSILE STRENGTH:	27 MPa	24 MPa
ELONGATION:	800%	750%
FORCE AT BREAK:		

INTEGRATED TECHNOLOGIES

- Advanced tactile sensitivity of natural latex.
- Fully textured surface delivers excellent grip.
- Extended cuff provides extra protection.



STANDARDS & CERTIFICATIONS

EN 420:2003	CAT III COMPLEX DESIGN
EN 374-3:2003	FDA 21 CFR 177.2600
EN 455:2000-2015	Component materials
ASTM D3578	comply with FDA
ASTM F1671-97	regulations for food
ASTM F739-12	contact.
REACH EU 1907/2006	
DIRECTIVE 2011/65/EU	
EU Reg 528/2012	



Please visit swsafety.com/symbols for information on standards used on this document.

ORDERING

ITEM NO.	SIZE	PALM WIDTH	PACKAGING
76051-020	S 6½ - 7	3.3" 85mm	50 gloves per box 10 boxes per case
76051-022	M 7½ - 8	3.8" 95mm	50 gloves per box 10 boxes per case
76051-024	L 8½ - 9	4.2" 105mm	50 gloves per box 10 boxes per case
76051-026	XL 9½ - 10	4.4" 113mm	50 gloves per box 10 boxes per case

APPLICATIONS

- Chemical Handling
- Epoxy and Fiberglass

STORAGE

Keep out of sunlight. Store in a cool, dry place. Keep away from sources of ozone or ignition.

