



Item Number ED150
Item Edetate Disodium, Dihydrate, USP
CAS Number 6381-92-6
Molecular Formula $C_{10}H_{14}N_2Na_2O_8 \cdot 2H_2O$
Molecular Weight 372.24

Assay (Dried Basis)..... 99.0-101.0%
pH of a 5% Solution @ 25°C..... 4.0-6.0
Loss on Drying..... 8.7-11.4%

MAXIMUM LIMITS

Calcium..... To pass test
Heavy Metals..... 0.005%
Nitrilotriacetic Acid..... To pass test
Residual Solvents..... To pass test

Spectrum Chemicals & Laboratory Products

Corporate Headquarters: **East Coast Facility:**
14422 S. San Pedro St. 755 Jersey Ave.
Gardena, CA 90248 New Brunswick, NJ 08901
(800) 813-1514 • (310) 516-8000 • Fax (310) 516-9843