

MPR-1412-PA | 76355-510  
MPR-1412R-PA | 76355-512

# Pharmaceutical Refrigerator

Accurate Temperature Control  
with Fast Recovery to Setpoint

48.2 cu.ft. | 1364 L  
48.0 cu.ft. | 1359 L



MPR-1412-PA | 76355-510

The MPR-1412-PA | 76355-510 (shelves) and MPR-1412R-PA | 76355-512 (drawers) large capacity pharmaceutical refrigerators deliver a stable, reliable temperature control for clinical, pharmaceutical, biomedical storage and processing. Fast recovery following door openings assures safe storage during frequent access and high ambient conditions.

## Positive Airflow Assures Interior Uniformity, Quick Recovery

Forced air circulation with multiple blowers assures temperature uniformity at all shelf or drawer levels. Internal shelf or drawer clearances permit the positive airflow to restore temperature to setpoint following door openings. This cooling power and airflow management compensates for heat generated by instrumentation if used inside the cabinet.



## Integrated Circulation Avoids Temperature Stratification

The top-mounted air circulation blower directs air down the back, side and front of the load in drawers or shelves. Regardless of location, all products are subjected to the same safe temperature to assure reproducibility and repeatability of stored product viability.

## Precision Temperature Control Protects Inventory

Many pharmaceuticals are degraded or destroyed if accidentally frozen. It is important that temperature be fixed at the desired setpoint in accordance with drug manufacturers' insert specifications. MPR Series pharmaceutical refrigerators operate over a temperature range of 2°C to 23°C, adjustable in 0.1°C increments. Actual interior air temperature is displayed to 0.1°C.



## Microprocessor Controls Simplify Operation

A microprocessor controller monitors all functions. Setpoint is factory set at 5°C. A soft key interface allows setpoint adjustment throughout the operating range. Audible and visual warnings indicate door ajar or temperature deviation from setpoint. The large easy-to-read LED display confirms actual temperature to 0.1°C resolution.

## Choose Your Inventory Management Preference

The MPR Series offers two choices in inventory control. The MPR-1412-PA | 76355-510 includes 8 open-wire shelves strong enough to hold bulk loads. The MPR-1412R-PA | 76355-512 is equipped with 8 heavy-duty, fully-extendable painted steel pull-out drawers for more convenient access to stored product such as large bottles or reagent kits.



## Choose Your Inventory Management Preference

Open wire shelves (shown) are adjustable to allow for loads of varying heights such as bulk storage, or process equipment such as fraction collectors. The Model MPR-1412R-PA | 76355-512 with pull-out drawers (not shown) offers a practical solution with additional ergonomic advantages for handling bulk items.

