

spectrum

chemicals & laboratory products

A Division of Spectrum Chemical Mfg. Corp.

Dear Customer,

This File Contains Both The ANSI Material Safety Data Sheet and The GHS Safety Data Sheet For The Same Product

Spectrum is currently transitioning all chemical product labeling from the ANSI¹ format to the GHS² format (see note below). In order to ensure that you receive complete labeling during the transition, we have included both the ANSI MSDS and the GHS SDS in a single file. The ANSI MSDS is given first, followed by the GHS SDS. Please use whichever matches the container label.

Why It Matters:

The complete precautionary labeling for this chemical consists of BOTH the label on the container AND the matching Material Safety Data Sheet (for ANSI labels) or Safety Data Sheet (for GHS labels). Both elements of the labeling [Label + (M)SDS] are written to be read and understood together, so as to provide complete precautionary information. It is intended for you to read and understood BOTH before handling or using the chemical.

Picking the Right One: 2 Easy Ways To Tell Whether Your Container Has an ANSI Label or a GHS Label

- 1) GHS labels: any pictogram displayed in the upper left-hand corner will be inside a red diamond.
ANSI labels: pictograms, if present, will be inside individual black boxes.
- 2) GHS labels: on the bottom of the right-hand panel of the label, locate the Lot Number. Directly to the left will be a string of control characters, followed by a single letter.
For GHS labels, the string of characters will end in "GHS:"

Label in ANSI Format

<p>CAUTION! MAY BE HARMFUL IF SWALLOWED MAY CAUSE EYE AND SKIN IRRITATION MAY AFFECT BEHAVIOR AND METABOLISM</p> <p>Do not taste or swallow. Avoid contact with eyes, skin and clothing. Avoid breathing mist or vapor. Avoid prolonged or repeated exposure. Use with adequate ventilation. Wash thoroughly after handling.</p> <p>FIRST AID: In case of contact, flush affected area with plenty of water for at least 15 minutes. Remove if worn. If inhaled, remove to fresh air. If not breathing, give artificial respiration. If breathing is difficult, give oxygen. If irritation persists, call a physician.</p> <p>KEEP FROM CHILDREN</p>	<p>SPECTRUM CHEMICALS & LABORATORY PRODUCTS</p> <p>BE159 SIZ SY</p> <p>Benzyl Benzoate (Benzoic Acid Phenylmethyl Ester)</p> <p>U.S.P. CAS 120-51-4</p> <p>CAUTION: For manufacturing, processing or repacking. Read and understand the label and Material Safety Data Sheet (MSDS) prior to use.</p> <p>For chemical emergency, call (800)424-9300</p> <p>www.SpectrumChemical.com</p>	<p>$C_{11}H_{12}O_2$ F.W. 212.24</p> <p>Assay 99.0-100.5% Specific Gravity @ 25°C 1.116-1.120 Congealing Temperature Min: 18.0°C Refractive Index @ 20°C 1.565-1.570 Acidity To pass test</p> <p>MAXIMUM LIMITS</p> <p>Aldehyde 0.05% Residual Solvents To pass test</p> <p>LIGHT SENSITIVE. Keep tightly closed in light-resistant containers.</p> <p>FLUSHED WITH NITROGEN</p>
---	--	---

Lot No. XQ###

SPECTRUM CHEMICAL MFG. CORP. Gardena, CA 90248 • New Brunswick, NJ 08901

CORPORATE OFFICES
14422 South San Pedro Street
Gardena, California 90248
PHONE 310.516.8000
FAX 310.516.9843

Label in GHS Format

WARNING!

- May irritate or sensitize • May cause central nervous system effects based on animal data
- Do not use or breathe • Wear protective gloves
- After handling
- WASH AND IF SWALLOWED: Call a POISON CENTER or doctor/physician if you feel unwell. Avoid reuse!

KEEP FROM CHILDREN

SPECTRUM™

BE159 SIZ SY

Benzyl Benzoate
(Benzoic Acid Phenylmethyl Ester)

U.S.P.
CAS 129-51-4

CAUTION: For manufacturing processing or shipping, read and understand the label and Safety Data Sheet (SDS) prior to use.

Chemical Emergency: (800)424-9088

www.SpectrumChemical.com

$C_{15}H_{14}O_2$ F.W. 212.24

Assay	99.0-100.5%
Specific Gravity @ 25°C	1.116-1.120
Freezing Temperature	Min. 18.0°C
Refractive Index @ 20°C	1.568-1.570
Acidity	To pass test

MAXIMUM LIMITS

Aldehyde	0.05%
Residual Solvents	To pass test

LIGHT SENSITIVE: Keep tightly closed in light-resistant containers.

FLUSHED WITH NITROGEN

Lot No. XQ####

¹ American National Standards Institute

² Globally Harmonized System for Hazard Communication

Sincerely,

Regulatory Affairs

SAFETY DATA SHEET

Preparation Date: 3/02/2015

Revision Date: 3/02/2015

Revision Number: G1

Product identifier

Product code: CH123
Product Name: CHLOROBUTANOL, ANHYDROUS, NF

Other means of identification

Synonyms: 1,1,1-Trichloro-2-methyl-2-propanol
Trichloro-tert-butanol
Trichloro-t-butyl alcohol
1,1,1-Trichloro-tert-butyl alcohol

CAS #: 57-15-8
RTECS # UC0175000
CI#: Not available

Recommended use of the chemical and restrictions on use

Recommended use: Plasticizer for cellulose esters and ethers. Antimicrobial. Preservative.
Uses advised against No information available

Supplier: Spectrum Chemicals and Laboratory Products, Inc.
14422 South San Pedro St.
Gardena, CA 90248
(310) 516-8000

Order Online At: <https://www.spectrumchemical.com>

Emergency telephone number Chemtrec 1-800-424-9300
Contact Person: Martin LaBenz (West Coast)
Contact Person: Ibad Tirmiz (East Coast)

2. HAZARDS IDENTIFICATION

Classification

This chemical is considered hazardous by the 2012 OSHA Hazard Communication Standard (29 CFR 1910.1200)

Acute toxicity - Oral	Category 4
Skin corrosion/irritation	Category 2
Serious eye damage/eye irritation	Category 2A
Specific target organ toxicity (single exposure)	Category 3

Label elements

Warning**Hazard statements**

Harmful if swallowed
Causes skin irritation
Causes serious eye irritation
May cause respiratory irritation

**Hazards not otherwise classified (HNOC)**

Not Applicable

Other hazards

Not available

Precautionary Statements - Prevention

Wash face, hands and any exposed skin thoroughly after handling
Do not eat, drink or smoke when using this product
Avoid breathing dust/fume/gas/mist/vapors/spray
Wear protective gloves
Wear eye/face protection
Use only outdoors or in a well-ventilated area

Precautionary Statements - Response

Specific treatment (see .? on this label)
IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. If eye irritation persists: Get medical advice/attention.
IF ON SKIN: Wash with plenty of soap and water
If skin irritation occurs: Get medical advice/attention
Take off contaminated clothing and wash before reuse
IF INHALED: Remove victim to fresh air and keep at rest in a position comfortable for breathing. Call a POISON CENTER or doctor/physician if you feel unwell.
IF SWALLOWED: Call a POISON CENTER or doctor/physician if you feel unwell
Rinse mouth

Precautionary Statements - Storage

Store in a well-ventilated place. Keep container tightly closed
Store locked up

Precautionary Statements - Disposal

Dispose of contents/container to an approved waste disposal plant

3. COMPOSITION/INFORMATION ON INGREDIENTS

Components	CAS-No.	Weight %	Trade Secret
Chlorobutanol, Anhydrous 57-15-8	57-15-8	100	*

4. FIRST AID MEASURES

4. FIRST AID MEASURES

First aid measures

General Advice:

Poison information centres in each State capital city can provide additional assistance for scheduled poisons (13 1126).

Skin Contact:

Wash off immediately with soap and plenty of water removing all contaminated clothes and shoes. Get medical attention. If skin irritation persists, call a physician.

Eye Contact:

Flush eye with water for 15 minutes. Get medical attention.

Inhalation:

Move to fresh air. If not breathing, give artificial respiration. If breathing is difficult, give oxygen. Get medical attention.

Ingestion:

Do not induce vomiting without medical advice. Never give anything by mouth to an unconscious person. Obtain medical attention.

Most important symptoms and effects, both acute and delayed

Symptoms

Irritating to eyes and skin. Irritating to respiratory system. Central nervous system effects. Drowsiness. Sedation. Sleepiness.

Indication of any immediate medical attention and special treatment needed

Notes to Physician:

Treat symptomatically

Protection of first-aiders

First-Aid Providers: Avoid exposure to blood or body fluids. Wear gloves and other necessary protective clothing. Dispose of contaminated clothing and equipment as bio-hazardous waste

5. FIRE-FIGHTING MEASURES

Extinguishing Media

Suitable Extinguishing Media:

Dry chemical. Carbon dioxide (CO₂). Water spray mist or foam.

Unsuitable Extinguishing Media:

No information available.

Specific hazards arising from the chemical

Hazardous Combustion Products:

Carbon monoxide; Carbon dioxide; hydrogen chloride

Specific hazards:

May be combustible at high temperatures.

Special Protective Actions for Firefighters

Specific Methods:

No information available.

Special Protective Equipment for Firefighters:

As in any fire, wear self-contained breathing apparatus pressure-demand, MSHA/NIOSH (approved or equivalent) and full protective gear

6. ACCIDENTAL RELEASE MEASURES

Personal precautions, protective equipment and emergency procedures

Personal Precautions: Keep people away from and upwind of spill/leak. Ensure adequate ventilation. Avoid contact with skin, eyes and clothing. Do not touch damaged containers or spilled material unless wearing appropriate protective clothing. Use personal protective equipment. Avoid dust formation. Remove all sources of ignition.

Environmental precautions Prevent further leakage or spillage if safe to do so. Prevent product from entering drains. Prevent from entering into soil, ditches, sewers, waterways, and/or ground water.

Methods and material for containment and cleaning up

Methods for containment Stop leak if you can do it without risk. Cover with plastic sheet to prevent spreading.

Methods for cleaning up Sweep up and shovel into suitable containers for disposal. Clean contaminated surface thoroughly.

7. HANDLING AND STORAGE

Precautions for safe handling

Technical Measures/Precautions:

Provide sufficient air exchange and/or exhaust in work rooms. Avoid dust formation. Keep away from incompatible materials.

Safe Handling Advice:

Wear personal protective equipment. Avoid contact with skin, eyes and clothing. Avoid dust formation. Do not ingest. Do not breathe vapours/dust. Keep away from heat and sources of ignition. Handle in accordance with good industrial hygiene and safety practice.

Conditions for safe storage, including any incompatibilities

Technical Measures/Storage Conditions:

Keep container tightly closed in a dry and well-ventilated place. Store at room temperature in the original container. Store away from incompatible materials.

Incompatible Materials:

Oxidizing agents. Acids. Acid anhydrides. Acid chlorides.

8. EXPOSURE CONTROLS/PERSONAL PROTECTION

Control parameters

National occupational exposure limits

United States

Components	OSHA	NIOSH	ACGIH	AIHA WHEEL
Chlorobutanol, Anhydrous - 57-15-8	None	None	None	None

Canada

Components	Alberta	British Columbia	Ontario	Quebec
Chlorobutanol, Anhydrous - 57-15-8	None	None	None	None

Australia and Mexico

Components	Australia	Mexico
Chlorobutanol, Anhydrous 57-15-8	None	None

Appropriate engineering controls

Engineering measures to reduce exposure:

Ensure adequate ventilation. Use process enclosures, local exhaust ventilation, or other engineering controls to keep airborne levels below recommended exposure limits. If user operations generate dust, fume or mist, use ventilation to keep exposure to airborne contaminants below the exposure limit.

Individual protection measures, such as personal protective equipment

Personal Protective Equipment

Eye protection: Goggles.

Skin and body protection: Long sleeved clothing. Chemical resistant apron. Gloves.

Respiratory protection: Effective dust mask. Wear respirator with dust filter..

Hygiene measures: Avoid contact with skin, eyes and clothing. Wash hands before breaks and immediately after handling the product. When using, do not eat, drink or smoke.

9. PHYSICAL AND CHEMICAL PROPERTIES

9. PHYSICAL AND CHEMICAL PROPERTIES

Physical state: Solid.	Appearance: Crystals. Crystalline.	Color: Colorless. White.
Odor: No information available	Taste No information available	Formula: C4H7Cl3O
Molecular/Formula weight: 177.46	Flash point (°C): No data available	Flashpoint (°C/°F): No information available.
Flash Point Tested according to: Not available	Lower Explosion Limit (%): No information available	Upper Explosion Limit (%): No information available
Autoignition Temperature (°C/°F): No information available	pH: No information available	Melting point/range(°C/°F): 97-99°C/206.6-210.2°C
Boiling point/range(°C/°F): 167°C/332.6°F	Decomposition temperature(°C/°F): No information available	Bulk density: No information available
Specific gravity: No information available	Vapor pressure @ 20°C (kPa): No information available	Density (g/cm3): No information available
Evaporation rate: No information available	Vapor density: No information available	VOC content (g/L): No information available
Odor threshold (ppm): No information available	Partition coefficient (n-octanol/water): 2.03	Viscosity: No information available
Miscibility: No information available	Solubility: Insoluble in cold water Soluble in hot water Soluble in Alcohol Soluble in Ether Soluble in Glycerol Soluble in Chloroform Soluble in Acetone	

10. STABILITY AND REACTIVITY

Reactivity

Reactive with oxidizing agents
Reactive with strong acids

Chemical stability

Stability:

Stable at normal conditions

Possibility of Hazardous Reactions:

Hazardous polymerization does not occur

Conditions to avoid:

Heat. Avoid dust formation. Incompatible materials.

Incompatible Materials:

Oxidizing agents. Acids. Acid anhydrides. Acid chlorides.

Hazardous decomposition products:

Carbon monoxide. Carbon dioxide. Hydrogen chloride.

Other Information

Corrosivity:

No information available

Special Remarks on Corrosivity: No information available

11. TOXICOLOGICAL INFORMATION

Information on likely routes of exposure

Principal Routes of Exposure:

Ingestion. Inhalation.

Acute Toxicity

Component Information

Chlorobutanol, Anhydrous - 57-15-8

LD50/oral/rat = 510 mg/kg

LD50/oral/mouse = No information available

LD50/dermal/rabbit = No information available

LD50/dermal/rat = No information available

LC50/inhalation/rat = No information available

LC50/inhalation/mouse = No information available

Other LD50 or LC50 information = No information available

Product Information

LD50/oral/rat =

VALUE- Acute Tox Oral = 510mg/kg

LD50/oral/mouse =

Value - Acute Tox Oral = No information available

LD50/dermal/rabbit

VALUE-Acute Tox Dermal = No information available

LD50/dermal/rat

VALUE -Acute Tox Dermal = No information available

LC50/inhalation/rat

VALUE-Vapor = No information available

VALUE-Gas = No information available

VALUE-Dust/Mist = No information available

LC50/Inhalation/mouse

VALUE-Vapor = No information available

VALUE - Gas = No information available

VALUE - Dust/Mist = No information available

Symptoms

Skin Contact: Causes skin irritation.

Eye Contact: Causes eye irritation.

Inhalation Irritating to respiratory system
Ingestion Harmful if swallowed. May affect behavior/central nervous system (somnolence, sedation, sleepiness).

Aspiration hazard No information available

Product code: CH123

Product name: CHLOROBUTANOL,
ANHYDROUS, NF

7 / 12

Delayed and immediate effects as well as chronic effects from short and long-term exposure

Chronic Toxicity May be habit forming if ingested for a prolonged period of time

Sensitization: No information available

Mutagenic Effects: May affect genetic material
Mutations in microorganisms
Experiments with bacteria and/or yeast have shown mutagenic effects

Carcinogenic effects: Not considered carcinogenic

Components	ACGIH - Carcinogens	IARC	NTP	OSHA HCS - Carcinogens	Australia - Prohibited Carcinogenic Substances	Australia - Notifiable Carcinogenic Substances
Chlorobutanol, Anhydrous	Not listed	Not listed	Not listed	Not listed	Not listed	Not listed

Reproductive toxicity No data is available

Reproductive Effects: No information available
Developmental Effects: No information available
Teratogenic Effects: No information available

Specific Target Organ Toxicity

STOT - single exposure No information available
STOT - repeated exposure No information available
Target Organs: Central nervous system.

12. ECOLOGICAL INFORMATION

Ecotoxicity

Ecotoxicity effects: No data available.
Persistence and degradability: No information available
Bioaccumulative potential: No information available
Mobility: No information available

13. DISPOSAL CONSIDERATIONS

Disposal Methods

Waste from residues / unused products:
Waste must be disposed of in accordance with Federal, State and Local regulation.

Contaminated packaging:
Empty containers should be taken for local recycling, recovery or waste disposal

Components	RCRA - F Series Wastes	RCRA - K Series Wastes	RCRA - P Series Wastes	RCRA - U Series Wastes
Chlorobutanol, Anhydrous	None	None	None	None

14. TRANSPORT INFORMATION

DOT

UN-No: Not Regulated
Proper Shipping Name: No information available
Hazard Class: No information available
Subsidiary Risk: Not applicable
Packing Group: None
ERG No: No information available
Marine Pollutant: No data available
DOT RQ (lbs): No information available

TDG (Canada)

UN-No: Not Regulated
Proper Shipping Name: No information available
Hazard Class: No information available
Subsidiary Risk: No information available
Packing Group: No information available
Description: No information available

ADR

UN-No: Not Regulated
Proper Shipping Name: No information available
Hazard Class: No information available
Packing Group: No information available
Subsidiary Risk: No information available
Classification Code: No information available
Description: No information available
CEFIC Tremcard No: No information available

IMO / IMDG

UN-No: Not Regulated
Proper Shipping Name: No information available
Hazard Class: No information available
Subsidiary Risk: No information available
Packing Group: No information available
Description: No information available
IMDG Page: No information available
Marine Pollutant: No information available
MFAG: No information available
Maximum Quantity: No information available

RID

UN-No: Not Regulated
Proper Shipping Name: No information available
Hazard Class: No information available
Subsidiary Risk: No information available
Packing Group: No information available
Classification Code: No information available
Description: No information available

14. TRANSPORT INFORMATION

ICAO

UN-No: Not Regulated
Proper Shipping Name: No information available
Hazard Class: No information available
Subsidiary Risk: No information available
Packing Group: No information available
Description: No information available

IATA

UN-No: Not Regulated
Proper Shipping Name: No information available
Hazard Class: No information available
Subsidiary Risk: No information available
Packing Group: No information available
Description: No information available

15. REGULATORY INFORMATION

International Inventories

Components	U.S. TSCA	KOREA KECL	Philippines (PICCS)	Japan ENCS	CHINA	Australia (AICS)	EINECS-No.
<i>Chlorobutanol, Anhydrous</i>	Present	Present KE-34081	Present	Present (2)-2002	Present [29435]	Present	Present 200-317-6

U.S. Regulations

California Prop. 65: Safe Drinking Water and Toxic Enforcement Act of 1986.

Chemicals Known to the State of California to Cause Cancer:

This product does not contain a chemical requiring a warning under California Prop. 65. (See table below)

Chemicals Known to the State of California to Cause Reproductive Toxicity:

This product does not contain a chemical requiring a warning under California Prop. 65. (See table below)

Components	Carcinogen	Developmental Toxicity	Male Reproductive Toxicity	Female Reproductive Toxicity:
<i>Chlorobutanol, Anhydrous</i>	Not Listed	Not Listed	Not Listed	Not Listed

CERCLA/SARA

Components	CERCLA - Hazardous Substances and their Reportable Quantities	Section 302 Extremely Hazardous Substances and TPQs	Section 302 Extremely Hazardous Substances and RQs	Section 313 - Chemical Category	Section 313 - Reporting <i>de minimis</i>
<i>Chlorobutanol, Anhydrous</i>	None	None	None	None	None

U.S. TSCA

Components	TSCA Section 5(a)2 - Chemicals With Significant New Use Rules (SNURS)	TSCA 8(d) -Health and Safety Reporting
<i>Chlorobutanol, Anhydrous</i>	Not Applicable	Not Applicable

Canada

WHMIS hazard class:

D2B Toxic materials

Product code: CH123

Product name: CHLOROBUTANOL,
ANHYDROUS, NF

10 / 12

Canada Controlled Products Regulation:

This product has been classified according to the hazard criteria of the CPR (Controlled Products Regulation) and the MSDS contains all of the information required by the CPR.

Inventory

Components	Canada (DSL)	Canada (NDSL)
Chlorobutanol, Anhydrous	Present	Not Listed

Components	CEPA Schedule I - Toxic Substances	CEPA - 2010 Greenhouse Gases Subject to Mandatory Reporting
Chlorobutanol, Anhydrous	Not listed	Not listed

EU Classification**R-phrase(s)**

R22 - Harmful if swallowed.

R36/37/38 - Irritating to eyes, respiratory system and skin.

S -phrase(s)

S26 - In case of contact with eyes, rinse immediately with plenty of water and seek medical advice.

Components	Classification	Concentration Limits:	Safety Phrases
Chlorobutanol, Anhydrous		No information	

The product is classified in accordance with Annex VI to Directive 67/548/EEC

Indication of danger:

Xi - Irritant.

Xn - Harmful.

**16. OTHER INFORMATION**

16. OTHER INFORMATION



Preparation Date: 3/02/2015
Revision Date: 3/02/2015
Prepared by: Sonia Owen

Disclaimer:

All chemicals may pose unknown hazards and should be used with caution. This Safety Data Sheet (SDS) applies only to the material as packaged. If this product is combined with other materials, deteriorates, or becomes contaminated, it may pose hazards not mentioned in this SDS. The physical properties reported in this SDS are obtained from the literature and do not constitute product specifications. Information contained herein does not constitute a warranty, whether expressed or implied, as to the safety, merchantability or fitness of the goods for a particular purpose. Spectrum Chemicals & Laboratory Products, Inc. assumes no responsibility for results obtained or for incidental or consequential damages, including lost profits, arising from the use of these data. No warranty against infringement of any patent, copyright or trademark is made or implied. It shall be the user's responsibility to develop proper methods of handling and personal protection based on the actual conditions of use. While this SDS is based on technical data judged to be reliable, Spectrum assumes no responsibility for the completeness or accuracy of the information contained herein.

End of Safety Data Sheet

Material Safety Data Sheet

<p>NFPA</p> 	<p>HMIS</p> <table border="1" style="margin-left: auto; margin-right: auto;"> <tr> <td style="background-color: #00FFFF;">Health Hazard</td> <td style="text-align: center; font-weight: bold;">2</td> </tr> <tr> <td style="background-color: #FFCCCC;">Fire Hazard</td> <td style="text-align: center; font-weight: bold;">1</td> </tr> <tr> <td style="background-color: #FFFF00;">Reactivity</td> <td style="text-align: center; font-weight: bold;">0</td> </tr> </table>	Health Hazard	2	Fire Hazard	1	Reactivity	0	<p>Personal Protective Equipment</p>  <p>See Section 15.</p>
Health Hazard	2							
Fire Hazard	1							
Reactivity	0							

Section 1. Chemical Product and Company Identification		Page Number: 1
Common Name/ Trade Name	Chlorobutanol	Catalog Number(s). YY1080, CH123
Manufacturer	SPECTRUM CHEMICAL MFG. CORP. 14422 S. SAN PEDRO STREET GARDENA, CA 90248	CAS# 57-15-8
Commercial Name(s)	Not available.	RTECS UC0175000
Synonym	1,1,1-Trichloro-tert-butyl alcohol	TSCA TSCA 8(b) inventory: Chlorobutanol
Chemical Name	Not available.	CI# Not available.
Chemical Family	Not available.	<p style="color: blue; margin: 0;">IN CASE OF EMERGENCY</p> <p style="color: blue; margin: 0;">CHEMTREC (24hr) 800-424-9300</p> <p style="margin: 0;">CALL (310) 516-8000</p>
Chemical Formula	C4H6CL3OH	
Supplier	SPECTRUM CHEMICAL MFG. CORP. 14422 S. SAN PEDRO STREET GARDENA, CA 90248	

Section 2. Composition and Information on Ingredients					
Name	CAS #	Exposure Limits			% by Weight
		TWA (mg/m ³)	STEL (mg/m ³)	CEIL (mg/m ³)	
1) Chlorobutanol	57-15-8				100
Toxicological Data on Ingredients	<p>Chlorobutanol LD50: Not available. LC50: Not available.</p>				

Section 3. Hazards Identification	
Potential Acute Health Effects	Extremely hazardous in case of ingestion. Very hazardous in case of skin contact (irritant), of eye contact (irritant), of inhalation. Inflammation of the eye is characterized by redness, watering, and itching. Skin inflammation is characterized by itching, scaling, reddening, or, occasionally, blistering.
Potential Chronic Health Effects	<p>Extremely hazardous in case of ingestion. Very hazardous in case of skin contact (irritant), of eye contact (irritant), of inhalation.</p> <p>CARCINOGENIC EFFECTS: Not available. MUTAGENIC EFFECTS: Not available. TERATOGENIC EFFECTS: Not available. DEVELOPMENTAL TOXICITY: Not available.</p> <p>The substance is toxic to liver. Repeated or prolonged exposure to the substance can produce target organs damage. Repeated or prolonged inhalation of dust may lead to chronic respiratory irritation.</p>

Continued on Next Page

Section 4. First Aid Measures

Eye Contact	Check for and remove any contact lenses. Immediately flush eyes with running water for at least 15 minutes, keeping eyelids open. Cold water may be used. Do not use an eye ointment. Seek medical attention.
Skin Contact	After contact with skin, wash immediately with plenty of water. Gently and thoroughly wash the contaminated skin with running water and non-abrasive soap. Be particularly careful to clean folds, crevices, creases and groin. Cover the irritated skin with an emollient. If irritation persists, seek medical attention. Wash contaminated clothing before reusing.
Serious Skin Contact	Wash with a disinfectant soap and cover the contaminated skin with an anti-bacterial cream. Seek medical attention.
Inhalation	Allow the victim to rest in a well ventilated area. Seek immediate medical attention.
Serious Inhalation	Not available.
Ingestion	Do not induce vomiting. Loosen tight clothing such as a collar, tie, belt or waistband. If the victim is not breathing, perform mouth-to-mouth resuscitation. Seek immediate medical attention.
Serious Ingestion	Not available.

Section 5. Fire and Explosion Data

Flammability of the Product	May be combustible at high temperature.
Auto-Ignition Temperature	Not available.
Flash Points	Not available.
Flammable Limits	Not available.
Products of Combustion	These products are carbon oxides (CO, CO ₂).
Fire Hazards in Presence of Various Substances	Not available.
Explosion Hazards in Presence of Various Substances	Risks of explosion of the product in presence of mechanical impact: Not available. Risks of explosion of the product in presence of static discharge: Not available.
Fire Fighting Media and Instructions	SMALL FIRE: Use DRY chemical powder. LARGE FIRE: Use water spray, fog or foam. Do not use water jet.
Special Remarks on Fire Hazards	Not available.
Special Remarks on Explosion Hazards	Not available.

Section 6. Accidental Release Measures

Small Spill	Use appropriate tools to put the spilled solid in a convenient waste disposal container. Finish cleaning by spreading water on the contaminated surface and dispose of according to local and regional authority requirements.
Large Spill	Use a shovel to put the material into a convenient waste disposal container. Finish cleaning by spreading water on the contaminated surface and allow to evacuate through the sanitary system.

Section 7. Handling and Storage

Precautions	Keep away from heat. Keep away from sources of ignition. Empty containers pose a fire risk, evaporate the residue under a fume hood. Ground all equipment containing material. Do not breathe dust. In case of insufficient ventilation, wear suitable respiratory equipment. If you feel unwell, seek medical attention and show the label when possible. Avoid contact with skin and eyes.
Storage	Keep container dry. Keep in a cool place. Ground all equipment containing material. Keep container tightly closed. Keep in a cool, well-ventilated place. Combustible materials should be stored away from extreme heat and away from strong oxidizing agents.

Section 8. Exposure Controls/Personal Protection

Engineering Controls	Use process enclosures, local exhaust ventilation, or other engineering controls to keep airborne levels below recommended exposure limits. If user operations generate dust, fume or mist, use ventilation to keep exposure to airborne contaminants below the exposure limit.
Personal Protection	Splash goggles. Lab coat. Dust respirator. Be sure to use an approved/certified respirator or equivalent. Gloves.
Personal Protection in Case of a Large Spill	Splash goggles. Full suit. Dust respirator. Boots. Gloves. A self contained breathing apparatus should be used to avoid inhalation of the product. Suggested protective clothing might not be sufficient; consult a specialist BEFORE handling this product.
Exposure Limits	Not available.

Section 9. Physical and Chemical Properties

Physical state and appearance	Solid.	Odor	Not available.
Molecular Weight	177.46 g/mole	Taste	Not available.
pH (1% soln/water)	Not available.	Color	Not available.
Boiling Point	167°C (332.6°F)		
Melting Point	99°C (210.2°F)		
Critical Temperature	Not available.		
Specific Gravity	Not available.		
Vapor Pressure	Not applicable.		
Vapor Density	Not available.		
Volatility	Not available.		
Odor Threshold	Not available.		
Water/Oil Dist. Coeff.	Not available.		
Ionicity (in Water)	Not available.		
Dispersion Properties	Not available.		
Solubility	Very slightly soluble in cold water.		

Section 10. Stability and Reactivity Data

Stability	The product is stable.
Instability Temperature	Not available.
Conditions of Instability	Not available.
Incompatibility with various substances	Not available.

Continued on Next Page

Corrosivity Non-corrosive in presence of glass.

Special Remarks on Reactivity Not available.

Special Remarks on Corrosivity Not available.

Polymerization No.

Section 11. Toxicological Information

Routes of Entry Eye contact. Inhalation. Ingestion.

Toxicity to Animals LD50: Not available.
LC50: Not available.

Chronic Effects on Humans The substance is toxic to liver.

Other Toxic Effects on Humans Extremely hazardous in case of ingestion.
Very hazardous in case of skin contact (irritant), of inhalation.

Special Remarks on Toxicity to Animals Not available.

Special Remarks on Chronic Effects on Humans Not available.

Special Remarks on other Toxic Effects on Humans Not available.

Section 12. Ecological Information

Ecotoxicity Not available.

BOD5 and COD Not available.

Products of Biodegradation Possibly hazardous short term degradation products are not likely. However, long term degradation products may arise.

Toxicity of the Products of Biodegradation The products of degradation are more toxic.

Special Remarks on the Products of Biodegradation Not available.

Section 13. Disposal Considerations

Waste Disposal

Section 14. Transport Information


DOT Classification Not a DOT controlled material (United States).

Identification Not applicable.

Special Provisions for Transport Not applicable.

Continued on Next Page

DOT (Pictograms)



Section 15. Other Regulatory Information and Pictograms

Federal and State Regulations

TSCA 8(b) inventory: Chlorobutanol

California Proposition 65 Warnings

Other Regulations

OSHA: Hazardous by definition of Hazard Communication Standard (29 CFR 1910.1200).

Other Classifications

WHMIS (Canada) CLASS D-2B: Material causing other toxic effects (TOXIC).

DSCL (EEC) R38- Irritating to skin.
R41- Risk of serious damage to eyes.

HMIS (U.S.A.)

Health Hazard	2
Fire Hazard	1
Reactivity	0
Personal Protection	E


National Fire Protection Association (U.S.A.)




WHMIS (Canada) (Pictograms)




DSCL (Europe) (Pictograms)



TDG (Canada) (Pictograms)



ADR (Europe) (Pictograms)



Protective Equipment



Gloves.



Lab coat.



Dust respirator. Be sure to use an approved/certified respirator or equivalent. Wear appropriate respirator when ventilation is inadequate.



Splash goggles.

Section 16. Other Information

MSDS Code C4046

References Not available.

Other Special Considerations Not available.

Validated by Sonia Owen on 8/11/2006.

Verified by Sonia Owen.

Printed 9/11/2006.

CALL (310) 516-8000

Notice to Reader

All chemicals may pose unknown hazards and should be used with caution. This Material Safety Data Sheet (MSDS) applies only to the material as packaged. If this product is combined with other materials, deteriorates, or becomes contaminated, it may pose hazards not mentioned in this MSDS. It shall be the user's responsibility to develop proper methods of handling and personal protection based on the actual conditions of use. While this MSDS is based on technical data judged to be reliable, Spectrum Quality Products, Inc. assumes no responsibility for the completeness or accuracy of the information contained herein.