

27 / Nov / 2014

Manufacturer:

AJ Innuscreen GmbH
Robert-Rössle-Straße 10
13125 Berlin/ Germany
Phone: +49 (0) 36 41 / 77-94 00
Fax: +49 (0) 36 41 / 77-76 77 76

Distribution:

Analytik Jena AG
Konrad-Zuse-Straße 1
07745 Jena/ Germany
www.bio.analytik-jena.com
lifescience@analytik-jena.com

Summary - Material Safety Data Sheet

Product: innuPREP Proteinase K

Summary may be subject to change – please inquire for update!

Safety data sheet

according to 1907/2006/EC, Article 31

Printing date 09.05.2014

Revision: 09.05.2014

1 Identification of the substance/mixture and of the company/undertaking

Product identifier

Trade name: Proteinase K

Relevant identified uses of the substance or mixture and uses advised against:

No further relevant information available.

Application of the substance / the preparation: Enzyme

Uses advised against: No further information available.

Details of the supplier of the safety data sheet

Manufacturer/Supplier:

AJ Innuscreen GmbH
 Biotechnologie Campus Berlin-Buch
 Robert-Rössle-Straße 10
 13125 Berlin
 Germany
 Tel.: +49 (0) 30 9489-3380
 Fax: +49 (0) 30 9489-3381
 E-Mail: info@aj-innuscreen.de
 Internet: www.aj-innuscreen.de

Further information obtainable from: AJ Innuscreen GmbH

Emergency telephone number:

Poison Information Center / Vergiftungs-Informationen-Zentrale
 Mathildenstraße 1
 79106 Freiburg
 Germany
 Tel.: +49 (0) 761 19 240
 Fax: +49 (0) 761 270 44 570
 E-Mail: giftinfo@uniklinik-freiburg.de
 Internet: www.giftberatung.de

2 Hazards identification

Classification of the substance or mixture:

Classification according to Regulation (EC) No 1272/2008:



GHS08 health hazard

Resp. Sens. 1 H334 May cause allergy or asthma symptoms or breathing difficulties if inhaled.



GHS07

Skin Irrit. 2 H315 Causes skin irritation.
 Eye Irrit. 2 H319 Causes serious eye irritation.
 Skin Sens. 1 H317 May cause an allergic skin reaction.
 STOT SE 3 H335 May cause respiratory irritation.

Classification according to Directive 67/548/EEC or Directive 1999/45/EC:



Xi; Irritant

R36/37/38: Irritating to eyes, respiratory system and skin.



Xn; Sensitising

R42: May cause sensitisation by inhalation.

(Contd. on page 2)

Safety data sheet

according to 1907/2006/EC, Article 31

Printing date 09.05.2014

Revision: 09.05.2014

Trade name: Proteinase K

(Contd. of page 1)

Information concerning particular hazards for human and environment:

The product has to be labelled due to the calculation procedure of the "General Classification guideline for preparations of the EU" in the latest valid version.

Classification system:

The classification is according to the latest editions of the EU-lists, and extended by company and literature data.

Label elements:**Labelling according to Regulation (EC) No 1272/2008:**

The product is classified and labelled according to the CLP regulation.

Hazard pictograms:

GHS07 GHS08

Signal word: Danger**Hazard-determining components of labelling:**

Proteinase, Tritirachium album serine

Hazard statements:

H315 Causes skin irritation.

H319 Causes serious eye irritation.

H334 May cause allergy or asthma symptoms or breathing difficulties if inhaled.

H317 May cause an allergic skin reaction.

H335 May cause respiratory irritation.

Precautionary statements:

P101 If medical advice is needed, have product container or label at hand.

P102 Keep out of reach of children.

P103 Read label before use.

P261 Avoid breathing dust/fume/gas/mist/vapours/spray.

P280 Wear protective gloves/protective clothing/eye protection/face protection.

P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

P362 Take off contaminated clothing and wash before reuse.

P405 Store locked up.

P501 Dispose of contents/container in accordance with local/regional/national/international regulations.

Other hazards:**Results of PBT and vPvB assessment:**

PBT: Not applicable.

vPvB: Not applicable.

3 Composition/information on ingredients

Chemical characterization: Mixtures

Description: Mixture of substances listed below with nonhazardous additions.

Dangerous components:

CAS: 39450-01-6 Proteinase, Tritirachium album serine 50-100%

EINECS: 254-457-8 Xi R36/37/38; Xn R42

Resp. Sens. 1, H334; Skin Irrit. 2, H315; Eye Irrit. 2, H319

Additional information: For the wording of the listed risk phrases refer to section 16.

(Contd. on page 3)

Safety data sheet

according to 1907/2006/EC, Article 31

Printing date 09.05.2014

Revision: 09.05.2014

Trade name: Proteinase K

(Contd. of page 2)

4 First aid measures

Description of first aid measures:**General information:**

Take affected persons out into the fresh air.

Do not leave affected persons unattended.

After inhalation: Supply fresh air; consult doctor in case of complaints.**After skin contact:** Immediately wash with water and soap and rinse thoroughly.**After eye contact:**

Rinse opened eye for several minutes under running water.

Call a doctor immediately.

After swallowing: Call a doctor immediately.**Information for doctor:****Most important symptoms and effects, both acute and delayed:**

No further relevant information available.

Indication of any immediate medical attention and special treatment needed:

No further relevant information available.

5 Firefighting measures

Extinguishing media:**Suitable extinguishing agents:**CO₂, powder or water spray. Fight larger fires with water spray or alcohol resistant foam.**For safety reasons unsuitable extinguishing agents:** Water with full jet**Special hazards arising from the substance or mixture:** No further relevant information available.**Advice for firefighters:****Protective equipment:** Do not inhale explosion gases or combustion gases.

6 Accidental release measures

Personal precautions, protective equipment and emergency procedures: Wear protective clothing.**Environmental precautions:** Do not allow product to reach sewage system or any water course.**Methods and material for containment and cleaning up:** Pick up mechanically.**Reference to other sections:**

See Section 7 for information on safe handling.

See Section 8 for information on personal protection equipment.

See Section 13 for disposal information.

7 Handling and storage

Handling:**Precautions for safe handling:** Store in cool, dry place in tightly closed receptacles.**Information about fire - and explosion protection:** No special measures required.**Conditions for safe storage, including any incompatibilities:****Storage:****Requirements to be met by storerooms and receptacles:** Store only in the original receptacle.**Information about storage in one common storage facility:** Not required.**Further information about storage conditions:** None.**Specific end use(s):** No further relevant information available.

(Contd. on page 4)

Safety data sheet

according to 1907/2006/EC, Article 31

Printing date 09.05.2014

Revision: 09.05.2014

Trade name: Proteinase K

(Contd. of page 3)

8 Exposure controls/personal protection

Additional information about design of technical facilities: No further data; see item 7.

Control parameters:

Ingredients with limit values that require monitoring at the workplace:

The product does not contain any relevant quantities of materials with critical values that have to be monitored at the workplace.

Exposure controls:

Personal protective equipment:

General protective and hygienic measures:

The usual precautionary measures are to be adhered to when handling chemicals.

Immediately remove all soiled and contaminated clothing

Wash hands before breaks and at the end of work.

Avoid contact with the eyes and skin.

Respiratory protection: Use suitable respiratory protective device in case of insufficient ventilation.

Protection of hands: Protective gloves

Material of gloves:

The selection of the suitable gloves does not only depend on the material, but also on further marks of quality and varies from manufacturer to manufacturer. As the product is a preparation of several substances, the resistance of the glove material can not be calculated in advance and has therefore to be checked prior to the application.

Penetration time of glove material:

The exact break through time has to be found out by the manufacturer of the protective gloves and has to be observed.

Eye protection: Tightly sealed goggles

9 Physical and chemical properties

Information on basic physical and chemical properties:

General Information:

Appearance:

Form: Solid

Colour: White

Odour: Nearly odourless

Odour threshold: Not determined.

pH-value: 6.5 (as aqueous solution)

Change in condition:

Melting point/Melting range: Undetermined.

Boiling point/Boiling range: Undetermined.

Flash point: Not applicable.

Flammability (solid, gaseous): Not applicable.

Ignition temperature:

Decomposition temperature: Not determined.

Self-igniting: Product is not selfigniting.

Danger of explosion: Product does not present an explosion hazard.

Explosion limits:

Lower: Not determined.

Upper: Not determined.

(Contd. on page 5)

Safety data sheet

according to 1907/2006/EC, Article 31

Printing date 09.05.2014

Revision: 09.05.2014

Trade name: Proteinase K

(Contd. of page 4)

Vapour pressure:	Not applicable.
Density:	Not determined.
Relative density:	Not determined.
Vapour density:	Not determined.
Evaporation rate:	Not determined.
Solubility in / Miscibility with water:	Soluble.
Partition coefficient (n-octanol/water):	Not determined.
Viscosity:	
Dynamic:	Not determined.
Kinematic:	Not determined.
Other information:	No further relevant information available.

10 Stability and reactivity

Reactivity:
Chemical stability:
Thermal decomposition / conditions to be avoided:
No decomposition if used according to specifications.
Possibility of hazardous reactions: No dangerous reactions known.
Conditions to avoid: No further relevant information available.
Incompatible materials: No further relevant information available.
Hazardous decomposition products: No dangerous decomposition products known.

11 Toxicological information

Information on toxicological effects:
Acute toxicity:
Primary irritant effect:
on the skin: Irritant to skin and mucous membranes.
on the eye: Strong irritant with the danger of severe eye injury.
Sensitization:
Sensitization possible through skin contact.
Sensitization possible through inhalation.

12 Ecological information

Toxicity:
Aquatic toxicity: No further relevant information available.
Persistence and degradability: No further relevant information available.
Behaviour in environmental systems:
Bioaccumulative potential: No further relevant information available.
Mobility in soil: No further relevant information available.
Additional ecological information:
General notes:
Water hazard class 1 (German Regulation) (Self-assessment): slightly hazardous for water
Do not allow undiluted product or large quantities of it to reach ground water, water course or sewage system.
Results of PBT and vPvB assessment:
PBT: Not applicable.
vPvB: Not applicable.

(Contd. on page 6)

Safety data sheet

according to 1907/2006/EC, Article 31

Printing date 09.05.2014

Revision: 09.05.2014

Trade name: Proteinase K

(Contd. of page 5)

Other adverse effects: No further relevant information available.

13 Disposal considerations

Waste treatment methods:**Recommendation:**

Must not be disposed together with household garbage. Do not allow product to reach sewage system.

Uncleaned packaging:**Recommendation:** Disposal must be made according to official regulations.

14 Transport information

UN-Number:

ADR, ADN, IMDG, IATA Void

UN proper shipping name:

ADR, ADN, IMDG, IATA Void

Transport hazard class(es):

ADR, ADN, IMDG, IATA

Class: Void**Packing group:**

ADR, IMDG, IATA Void

Environmental hazards:**Marine pollutant:** No**Special precautions for user:** Not applicable.**Transport in bulk according to Annex II of MARPOL73/78 and the IBC Code** Not applicable.**UN "Model Regulation":** -

15 Regulatory information

Safety, health and environmental regulations/legislation specific for the substance or mixture:

No further relevant information available.

Chemical safety assessment: A Chemical Safety Assessment has not been carried out.

16 Other information

This information is based on our present knowledge. However, this shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.

Relevant phrases:

H315 Causes skin irritation.

H319 Causes serious eye irritation.

H334 May cause allergy or asthma symptoms or breathing difficulties if inhaled.

R36/37/38 Irritating to eyes, respiratory system and skin.

R42 May cause sensitisation by inhalation.

Department issuing MSDS:

ADEEGO GmbH

Wildsteig 4A

42113 Wuppertal

Tel. +49 (0) 202 2575700

Fax. +49 (0) 202 2575701

Email: mail@adeego.de

www.adeego.de

(Contd. on page 7)

Safety data sheet

according to 1907/2006/EC, Article 31

Printing date 09.05.2014

Revision: 09.05.2014

Trade name: Proteinase K

(Contd. of page 6)

Contact: Dr. Axel Deeg**Abbreviations and acronyms:**

ADR: Accord européen sur le transport des marchandises dangereuses par Route (European Agreement concerning the International Carriage of Dangerous Goods by Road)

IMDG: International Maritime Code for Dangerous Goods

IATA: International Air Transport Association

GHS: Globally Harmonized System of Classification and Labelling of Chemicals

EINECS: European Inventory of Existing Commercial Chemical Substances

ELINCS: European List of Notified Chemical Substances

CAS: Chemical Abstracts Service (division of the American Chemical Society)

Skin Irrit. 2: Skin corrosion/irritation, Hazard Category 2

Eye Irrit. 2: Serious eye damage/eye irritation, Hazard Category 2

Resp. Sens. 1: Sensitisation - Respirat., Hazard Category 1

Skin Sens. 1: Sensitisation - Skin, Hazard Category 1

STOT SE 3: Specific target organ toxicity - Single exposure, Hazard Category 3
