

Commerical Kitting

Specimen shipping systems

Complete shipping systems for category A, category B, exempt human and exempt animal specimen transport via air or ground



Refrigerated

Temperature control specimen shipping systems

Therapak has taken the worry out of shipping biological substances and complying with shipping and labeling requirements. Therapak's Specimen Shipping Systems have been validated to meet requirements for shipping biological samples, whether they are Category B specimens, Exempt Animal or Exempt Human specimens, or Category A, Infections samples.

Components for each shipping system can be purchased separately by visiting vwr.com, contacting your VWR sales rep or calling 800 932 5000.

When your sample shipment requires that your payload is kept refrigerated, our compliant systems are validated to maintain necessary temperatures, and are available with refrigerant packs or foam bricks.

Each shipping system includes a shipping box, an Expanded Polystyrene (EPS) cooler, leakproof bag, absorbent material, and necessary labeling and inserts.



REFRIGERATED SHIPPING SYSTEMS

Cat. No. 10029-294

ID, cm (in.)	OD, cm (in.)	Capacity	Compliance	Leakproof Feature	Therapak No.	VWR Cat. No.	Case of
17.78 x 17.78 x 17.78 (7 x 7 x 7")	30.48 x 30.48 x 30.48 (12 x 12 x 12»)	6 Tubes	Cat. B	95kPa Bag	38028	10029-326	4
20.32 x 15.24 x 20.32 (8 x 6 x 8")	27.94 x 22.86 x 27.94 (11 x 9 x 11»)	6 Tubes	Cat. B	95kPa Bag	37813	10029-294	3
27.94 x 25.4 x 22.86 (11 x 10 x 9")	35.56 x 33.02 x 30.48 (14 x 13 x 12»)	6 Tubes	Cat. B	95kPa Bags (3)	37816	10029-300	4
*27.94 x 25.4 x 22.86 (11 x 10 x 9")	35.56 x 33.02 x 30.48 (14 x 13 x 12»)	6 Tubes	Cat. B	95kPa Canister	37814	10029-296	4
20.32 x 15.24 x 20.32 (8 x 6 x 8")	27.94 x 22.86 x 27.94 (11 x 9 x 11»)	6 Tubes	Cat. B	95kPa Canister	37815	10029-298	3
15.24 x 11.43 x 6.35 (6 x 4½ x 2½")	20.32 x 15.24 x 11.43 (8 x 6 x 4½»)	6 Tubes	Cat. B/EHS	95kPa Bag	380125	76203-686	12
20.32 x 15.24 x 10.16 (8 x 6 x 4")	27.94 x 22.86 x 18.42 (11 x 9 x 7¼»)	6 Tubes	Cat. B	95kPa Cansiter	380135	76203-890	12
20.32 x 15.24 x 20.32 (8 x 6 x 8")	27.94 x 22.86 x 27.94 (11 x 9 x 11")	6 Tubes	EAS	Biohazard Bag	56884	10029-418	3
20.32 x 15.24 x 20.32 (8 x 6 x 8")	27.94 x 22.86 x 27.94 (11 x 9 x 11")	6 Tubes	EHS	Biohazard Bag	56876	10029-402	3
25.4 x 25.4 x 25.4 (10 x 10 x 10")	30.48 x 30.48 x 30.48 (12 x 12 x 12")	1 Cat A Infectious Shipper (76203-876)	Overpack only	—	53212	—	2



Cat. No. 10029-296



Cat. No. 10029-298



Cat. No. 10029-294



Cat. No. 10029-300

COMPONENTS USED IN REFRIGERATED SYSTEMS

VWR Cat. No.	Therapak No.	Description	Details	Unit										
89170-936	10316G	Aqui-Pak Absorbent Pouch	6 Tubes	Case/ 900	●	●	●	●	●	●	●	●	●	●
89170-946	15650G	95kPa Bag	17.78 x 27.94 cm (7 x 11")	Pack/ 100		●	●			●				
10497-176	16022G	Biohazard Bag	15.24 x 22.86 cm (6 x 9")	Case/ 1,000	●								●	●
—	16022G	Zip Style Bag	15.24 x 15.24 cm (6 x 6")	Case/ 1,000					●					
—	16032G	Zip Style Bag	15.24 x 20.32 cm (6 x 8")	Case/ 1,000								●		
89170-958	23072G	95kPa Medium Canister	7.366 x 15.24 cm (2.9 dia. x 6")	Case/ 15					●			●		
89170-960	23080G	95kPa Large Canister	10.16 cm dia. x 17.78 cm (4" dia. x 7")	Case/ 6				●		●				
—	233175G-P	White Box	20.32 x 16.51 x 11.43 cm (8 x 6½ x 4½")	Each										
—	23681G	Insert for #56514	1 Piece	Each	●									
—	44912G	Bubble Wrap Sheet	30.48 x 30.48 cm (12 x 12")	Roll/ 750	●	●	●	●	●			●		
10029-350	54530G	Dry Ice Label	13.97 x 13.97 cm (5½ x 5½")	Roll/ 500	●									
10029-352	54607G	Biohazard Symbol Label	7.62 x 5.08 cm (3 x 2")	Roll/ 1,000	●	●	●	●	●			●	●	●
10029-354	54609G	Frozen Specimen Label	5.08 x 2.22 cm (2 x 7⁄8")	Roll/ 1,000										
10029-356	54635G	Double Up Arrow Label	7.62 x 10.16 cm (3 x 4")	Roll/ 500	●				●			●		
10029-358	54751G	EAS Label	8.89 x 2.54 cm (3½ x 1")	Roll/ 500									●	
10029-360	54752G	UN3373 Category B Label	8.89 x 10.16 cm (3½ x 4")	Roll/ 500		●			●	●				
10029-362	54765G	List of Contents Card	10.8 x 6.8 (4¼ x 2¾")	Pack/ 100		●	●	●	●	●		●		
10029-364	54782G	UN3373 Category B Label	20.32 x 10.16 cm (8 x 4")	Roll/ 500			●	●	●	●				
10029-368	54814G	EHS Label	6.35 x 5.08 cm (2½ x 2")	Roll/ 500	●						●			●
—	56398G	Foam Brick	2.13 kg (7.5 oz.)	Case/ 72						●				
10029-372	56400G	Refrigerant Pack	2.27 kg (8 oz.)	Case/ 72		●	●	●	●	●		●	●	●
10497-252	56514G	Small EPS Cooler	27.94 x 25.4 x 27.94 cm (11 x 10 x 11")	Case/ 12	●	●			●				●	●
—	56524G	Box for Small EPS Cooler	27.94 x 25.4 x 27.94 cm (11 x 10 x 11")	Pack/ 25	●	●			●				●	●
—	56535G	Medium EPS Cooler	35.56 x 33.02 x 30.48 cm (14 x 13 x 12")	Case/ 4				●	●					
—	56536G	Box for Medium EPS Cooler	38.1 x 34.29 x 53.8 cm (15 x 13½ x 12¼")	Pack/ 25			●	●						
—	56908G	Small EPS Cooler	20.32 x 16.51 x 11.43 cm (8 x 6½ x 4½")	Case/ 12						●				
—	56947G	Extra Small EPS Cooler	27.94 x 22.86 x 18.42 cm (11 x 9 x 7¼")	Case/ 12								●		
—	56948G	Box for Extra Small EPS Cooler	28.58 x 23.5 x 19.05 cm (11¼ x 9¼ x 7½")	Pack/ 25								●		

Frozen



When your sample shipment requires that your payload is kept frozen, Therapak's compliant systems are designed to be used with dry ice to keep samples within the desired temperature range.

Each shipping system includes a shipping box, an EPS cooler, a leakproof secondary container, absorbent material, and necessary labeling and inserts.

Cat. No. 10029-304

FROZEN SHIPPING SYSTEMS

ID, cm (in.)	OD, cm (in.)	Capacity	Compliance	Leakproof Feature	Therapak No.	VWR Cat. No.	Case of
30.48 x 25.4 x 22.86 (12 x 10 x 9")	38.1 x 34.29 x 31.12 (15 x 13.5 x 12.25")	81 2mL Vials	Cat. B	Biohazard Bag	37818	10029-304	4
30.48 x 25.4 x 22.86 (12 x 10 x 9")	38.1 x 34.29 x 31.12 (15 x 13.5 x 12.25")	81 2mL Vials	Cat. B	Biohazard Bag	37808	10029-272	4
40.64 x 33.02 x 27.94 (16 x 13 x 11")	46.99 x 40.64 x 35.56 (18.5 x 16 x 14")	72 Tubes	Cat. B	95 kPa Canisters (4)	37809	10029-274	2
35.56 x 33.02 x 30.48 (14 x 13 x 12")	38.1 x 34.29 x 31.12 (15 x 13.5 x 12.25")	18 Tubes	Cat. B	95 kPa Canisters	37817	10029-302	2
20.32 x 17.78 x 20.32 (8 x 7 x 8")	27.94 x 25.4 x 27.94 (11 x 10 x 11")	6 Tubes	Cat. B	95 kPa Canisters	37806	10029-262	3
20.32 x 15.24 x 10.16 (8 x 6 x 4")	27.94 x 22.86 x 18.42 (11 x 9 x 7.25")	6 Tubes	Cat. B	95 kPa Canisters	380134	76203-888	12
20.32 x 17.78 x 20.32 (8 x 7 x 8")	27.94 x 25.4 x 27.94 (11 x 10 x 11")	6 Tubes	EAS	Biohazard Bag	56885	10029-420	3
20.32 x 17.78 x 20.32 (8 x 7 x 8")	27.94 x 25.4 x 27.94 (11 x 10 x 11")	6 Tubes	EHS	Biohazard Bag	56877	10029-404	3



Cat. No. 10029-272



Cat. No. 10029-274



Cat. No. 10029-404



Cat. No. 10029-262

COMPONENTS USED IN FROZEN SYSTEMS

VWR Cat. No.	Therapak No.	Description	Details	Unit	VWR Cat. No. and Therapak No.							
					10029-304 37818	10029-272 37808	10029-274 37809	10029-302 37817	10029-262 37806	76203-888 380134	10029-420 56885	10029-404 56877
89170-930	10307G	Absorbent Sheet	15.24 x 15.24 cm (6 x 6")	Pack/ 500	•							
76203-694	10312G	Absorbent Sheet	30.48 x 30.48 cm (12 x 12")	Pack/ 500		•						
89170-934	10314G	Aqui-Pak Absorbent Pouch	4 Tubes	Case/ 900								
89170-936	10316G	Aqui-Pak Absorbent Pouch	6 Tubes	Case/ 900			•	•	•	•	•	•
10497-176	16019G	Biohazard Bag	15.24 x 22.86 cm (6 x 9")	Case/ 1,000							•	•
—	16022G	Zip Style Bag	15.24 x 15.24 cm (6 x 6")	Case/ 1,000			•	•				
—	16032G	Zip Style Bag	20.32 x 20.32 cm (8 x 8")	Case/ 1,000						•		
—	16072G	Zip Style Bag	30.48 x 30.48 cm (12 x 12"), 2 mil	Case/ 1,000		•						
—	16079G	Zip Style Bag	20.32 x 25.4 cm (8 x 10"), 2 mil	Case/ 1,000	•				•			
10497-178	16532G	Biohazard Bag, 2-pouch	25.4 x 25.4 cm (10 x 10")	Case/ 1,000		•						
10497-180	16604G	Biohazard Bag, 2-pouch	20.32 x 20.32 cm (8 x 8")	Case/ 1,000	•							
89170-958	23072G	95kPa Medium Canister	7.37 cm dia x 15.24 cm (3" dia. x 6")	Case/ 15					•	•		
89170-960	23080G	95kPa Large Canister	10.16 cm dia x 17.78 cm (4" dia. x 7")	Case/ 6			•	•				
—	23677G	Insert for #56535	4 Pieces	Each	•	•						
—	23678G	Insert for #56517	4 Pieces	Each			•					
—	23681G	Insert for #56514	1 Pieces	Each					•		•	•
—	23801G	Medium "X" Insert	2 Pieces	Each				•				
—	34940G	Cryobox	7.62 cm (3"), 81 Slots	Case/ 10		•						
—	44912G	Bubble Wrap Sheet	30.48 x 30.48 cm (12 x 12")	Roll/ 750	•	•		•	•	•	•	•
10029-350	54530G	Dry Ice Label	13.97 x 13.97 cm (5½ x 5½")	Roll/ 500	•	•	•	•	•	•	•	•
10029-352	54607G	Biohazard Symbol Label	7.62 x 5.08 cm (3 x 2")	Roll/ 1,000	•	•	•	•	•	•	•	•
10029-354	54609G	Frozen Specimen Label	5.08 x 2.22 cm (2 x 7/8")	Roll/ 1,000			•	•				
10029-364	54635G	Double Up Arrow Label	7.62 x 10.16 cm (3 x 4")	Roll/ 500			•	•		•	•	•
10029-358	54751G	EAS Label	8.89 x 2.54 cm (3½ x 1")	Roll/ 500							•	
10029-360	54752G	UN3373 Category B Label	8.89 x 10.16 cm (3½ x 4")	Roll/ 500					•	•		
10029-362	54765G	Lists of Contents Card	10.8 x 6.8 cm (4¼ x 2¾")	Pack/ 100	•	•	•	•	•	•		
10029-364	54782G	UN3373 Category B Label	20.32 x 10.16 cm (8 x 4")	Roll/ 500	•	•	•	•				
10029-368	54814G	EHS Label	6.35 x 5.08 cm (2½ x 2")	Roll/ 500								•
10497-252	56514G	Small EPS Cooler	27.94 x 25.4 x 27.94 cm (11 x 10 x 11")	Case/ 12					•		•	•
—	56517G	Large EPS Cooler	45.72 x 40.64 x 35.56 cm (18 x 16 x 14")	Case/ 2			•					
—	56518G	Box for Large EPS Cooler	46.99 x 40.64 x 35.56 cm (18½ x 16 x 14")	Pack/ 25			•					
—	56524G	Box for Small EPS Cooler	27.94 x 25.4 x 27.94 cm (11 x 10 x 11")	Pack/ 25					•		•	•
—	56535G	Medium EPS Cooler	35.56 x 33.02 x 30.48 cm (14 x 13 x 12")	Case/ 4	•	•		•				
—	56536G	Box for Medium EPS Cooler	38.1 x 34.29 x 31.12 cm (15 x 13½ x 12¼")	Pack/ 25	•	•		•				
—	34932G	Cryobox	5.08 cm (2"), 81 Slots	Each	•							
—	56947G	Extra Small EPS Cooler	27.94 x 22.86 x 18.42 cm (11 x 9 x 7¼")	Pack/ 25						•		
—	56948G-P	Box for Extra Small EPS Cooler	28.58 x 23.5 x 19.05 cm (11¼ x 9¼ x 7½")	Case/ 12						•		

Ambient



Therapak's Ambient Shipping Systems protect sample shipments that are not temperature sensitive.

Each Ambient Shipping System includes a tuck-style corrugated box for two-way shipping, leakproof secondary sample containment, absorbent material, and necessary labeling and inserts.

Cat. No. 10029-264

AMBIENT SHIPPING SYSTEMS

ID, cm (in.)	OD, cm (in.)	Capacity	Compliance	Leakproof Feature	Therapak No.	VWR Cat. No.	Case of
17.78 x 10.16 x 10.16 (7 x 4 x 4")	20.32 x 10.8 x 10.8 (8 x 4.25 x 4.25")	6 Tubes	Cat. B	95kPa Bag	37801	10029-264	16
28.58 x 22.86 x 5.33 (11.25 x 9 x 2.1")	31.12 x 23.5 x 5.72 (12.25 x 9.25 x 2.25")	6 Tubes	Cat. B	95kPa Bag	37980	10029-322	15
17.78 x 10.16 x 10.16 (7 x 4 x 4")	20.32 x 10.8 x 10.8 (8 x 4.25 x 4.25")	6 Tubes	Cat. B	95kPa Bag	37802	10029-266	16
17.78 x 10.16 x 10.16 (7 x 4 x 4")	20.32 x 10.8 x 10.8 (8 x 4.25 x 4.25")	6 Tubes	Cat. B	95kPa Bag	37812	10029-280	16
17.78 x 10.16 x 10.16 (7 x 4 x 4")	20.32 x 10.8 x 10.8 (8 x 4.25 x 4.25")	6 Tubes	Cat. B	95kPa Canister	37804	10029-258	16
17.78 x 10.16 x 10.16 (7 x 4 x 4")	20.32 x 10.8 x 10.8 (8 x 4.25 x 4.25")	6 Tubes	Cat. B	95kPa Canister	37805	10029-260	16
13.335 x 10.16 x 10.16 (5.25 x 4 x 4")	15.24 x 10.16 x 10.16 (6 x 4 x 4")	4 Tubes	Cat. B	95kPa Canister	37803	—	36
10.16 x 10.16 x 22.86 (4 x 4 x 9")	10.16 x 10.16 x 23.11 (4 x 4 x 9.1")	6 Swabs	Cat. B	95kPa Canister	38054	10029-328	12
27.94 x 26.67 x 10.16 (11 x 10.5 x 4")	28.58 x 27.31 x 10.8 (11.25 x 10.75 x 4.25")	18 Tubes	Cat. B	95kPa Bag	37925	10029-312	6
16.51 x 12.19 x 5.33 (6.5 x 4.8 x 2.1")	21.34 x 17.27 x 9.4 (8.4 x 6.8 x 3.7")	4 Tubes	Cat. B	95kPa Bag	380126	76203-668	12
15.24 x 7.62 x 6.35 (6 x 3 x 2.5")	18.03 x 10.67 x 9.14 (7.1 x 4.2 x 3.6")	6 Tubes	Cat. B	95kPa Bag	380127	76203-690	16
12.7 x 15.21 x 4.57 (5 x 6 x 1.8")	13.21 x 15.24 x 4.57 (5.2 x 6 x 1.8")	4 Tubes	Cat. B	95kPa Bag	380129	76203-882	20
27.94 x 26.67 x 10.16 (11 x 10.5 x 4")	28.58 x 27.31 x 10.8 (11.25 x 10.75 x 4.25")	6 x 9" Bag (x10)	EAS	Biohazard Bag	56883	10029-416	6
17.78 x 10.16 x 10.16 (7 x 4 x 4")	20.32 x 10.8 x 10.8 (8 x 4.25 x 4.25")	6 Tubes	EAS	Biohazard Bag	56882	10029-414	12
13.34 x 10.16 x 9.53 (5.25 x 4 x 3.75")	15.24 x 10.8 x 10.16 (6 x 4.25 x 4")	50 mL	EAS	Biohazard Bag	56880	10029-410	36
13.34 x 10.16 x 9.53 (5.25 x 4 x 3.75")	15.24 x 10.8 x 10.16 (6 x 4.25 x 4")	4 Tubes	EAS	Biohazard Bag	56881	10029-412	36
23.37 x 7.6 x 17.02 (9.2 x 3 x 6.7")	24.38 x 8.67 x 19.56 (9.6 x 3.4 x 7.7")	6 Tubes	EAS	Biohazard Bag	56879	10029-408	50
27.94 x 26.67 x 10.16 (11 x 10.5 x 4")	28.58 x 27.31 x 10.8 (11.25 x 10.75 x 4.25")	6 x 9" Bag (x10)	EHS	Biohazard Bag	56875	10029-400	6
17.78 x 10.16 x 10.16 (7 x 4 x 4")	20.32 x 10.8 x 10.8 (8 x 4.25 x 4.25")	6 Tubes	EHS	Biohazard Bag	56874	10029-398	16
13.34 x 10.16 x 9.53 (5.25 x 4 x 3.75")	15.24 x 10.8 x 10.16 (6 x 4.25 x 4")	50 mL	EHS	3-Wall Speciman Bag	56872	10029-394	36
13.34 x 10.16 x 9.53 (5.25 x 4 x 3.75")	15.24 x 10.8 x 10.16 (6 x 4.25 x 4")	4 Tubes	EHS	Biohazard Bag	56873	10029-396	36
23.37 x 7.6 x 17.02 (9.2 x 3 x 6.7")	24.38 x 8.64 x 19.56 (9.6 x 3.4 x 7.7")	6 Tubes	EHS	Biohazard Bag	56871	10029-392	50
		10 mL Tubes (x10)					
10.41 x 15.24 (4.1 x 6")	16.51 x 16.51 x 20.32 (6.5 x 6.5 x 8")	2 mL Tubes (x57)	Cat. A	95kPa Canister	563208	76203-876	12

The image shows a white box of Therapak and a biohazard specimen bag. The box is labeled 'THERAPAK' and has a blue band. The specimen bag is orange and white, featuring a biohazard symbol and the text 'EXEMPT ANIMAL SPECIMEN' and 'BIOHAZARDOUS MATERIALS'.

COMPONENTS USED IN AMBIENT SYSTEMS

Avantor Services | vwr.com/avantor_services 7

Combo shipping systems

When shipping multiple temperature range samples in one shipment, Therapak's Combo Shipping Systems can simplify the process by giving you everything you need to ship Frozen/Refrigerated or Frozen/Ambient sample combinations in a single and convenient box.

Each Shipping System includes a shipping box, Expanded Polystyrene (EPS) coolers, absorbent material, and labels and inserts necessary for compliant shipping.



Cat. No. 10029-422



Cat. No. 10029-386



Cat. No. 10029-406



Cat. No. 10029-276

ID, cm (in.)	OD, cm (in.)	Capacity	Combination	Compliance	Leakproof Feature	Therapak No.	VWR Cat. No.	Unit
20.32 x 15.24 x 20.32 (8 x 6 x 8")	28.19 x 24.38 x 33.53 (11.1 x 9.6 x 13.2")	6 Tubes	Frozen/Ambient	Cat. B	95kPa Bag/Canister	56726	10029-386	Case/ 6
40.64 x 30.48 x 40.64 (8 x 6 x 8"); Two Chambers	48.26 x 27.69 x 27.18 (19 x 10.9 x 10.7")	6 Tubes	Frozen/Refrigerated	Cat. B	95kPa Canister	37810	10029-276	Each
20.32 x 15.24 x 20.32 (8 x 6 x 8")	27.94 x 22.86 x 27.94 (11 x 9 x 11")	4 Tubes	Frozen/Ambient	EAS	Biohazard Bag	56886	10029-422	Case/ 2
20.32 x 15.24 x 20.32 (8 x 6 x 8")	28.19 x 22.86 x 27.94 (11.1 x 9.6 x 13.2")	4 Tubes	Frozen/Ambient	EHS	Biohazard Bag	56878	10029-406	Case/ 2

COMPONENTS USED IN COMBO SYSTEMS

VWR Cat. No.	Therapak No.	Description	Details	Unit	VWR Cat. No. and Therapak No.			
					10029-304 37818	10029-272 37808	10029-274 37809	10029-302 37817
89170-946	15650G	95kPa Bag	17.78 x 27.94 cm (7 x 11")	Pack/ 100	•			
—	16022G	Zip Style Bag	15.24 x 15.24 cm (6 x 6")	Case/ 1,000	•			
10497-176	16019G	Biohazard Bag	15.24 x 22.86 cm (6 x 9")	Case/ 1,000			•	•
89170-958	23072G	95kPa Medium Canister	7.37 cm dia x 15.24 cm (3" dia. x 6")	Case/ 15	•	•		
—	44912G	Bubble Wrap Sheet	30.48 x 30.48 cm (12 x 12")	Roll/ 750			•	•
89170-934	10314G	Aqui-Pak Absorbent Pouch	4 Tubes	Case/ 900			•	•
89170-936	10316G	Aqui-Pak Absorbent Pouch	6 Tubes	Case/ 900	•	•		
10029-352	54607G	Biohazard Symbol Label	7.62 x 5.08 cm (3 x 2")	Roll/ 1,000	•	•	•	•
10029-360	54752G	UN3373 Category B Label	8.89 x 10.16 cm (3.5 x 4")	Roll/ 500	•			
10029-364	54782G	UN3373 Category B Label	20.32 x 10.16 cm (8 x 4")	Roll/ 500		•		
10029-358	54751G	EAS Label	8.89 x 2.54 cm (3.5 x 1")	Roll/ 500			•	
10029-368	54814G	EHS Label	6.35 x 5.08 cm (2.5 x 2")	Roll/ 500				•
10029-356	54635G	Double Up Arrow Label	7.62 x 10.16 cm (3 x 4")	Roll/ 500		•	•	
10029-362	54765G	List of Contents Card	10.8 x 7 cm (4 1/4 x 2 3/4")	Pack/ 100	•	•		
10029-350	54530G	Dry Ice Label	13.97 x 13.97 cm (5 1/2 x 5 1/2")	Roll/ 500	•	•	•	•
10497-252	56514G	Small EPS Cooler	27.94 x 25.4 x 27.94 cm (11 x 10 x 11")	Case/ 12	•	•	•	•
—	23347G	Box or Vertical Combo Shipper	27.94 x 24.13 x 33.02 cm (11 x 9 1/2 x 13")	Pack/ 25	•		•	•
—	23800G	Small Corrugated Insert	24.13 x 13.97 cm (9 1/2 x 5 1/2")	Pack/ 25	•	•		•
—	56553G	Box for x2 #56514	48.26 x 27.94 x 27.31 cm (19 x 11 x 10 3/4")	Pack/ 25		•		
—	23681G	Insert for #56514	1 Piece	Each			•	
89171-108	56424G	Insulated Envelope	22.86 x 27.94 cm (9 x 11")	Case/ 32	•		•	•
89171-106	56422G	Ambient Gel Wrap	46.99 x 15.24 cm (18.5 x 6")	Case/ 16	•		•	•
10029-372	56400G	Refrigerant Pack	2.27 k (8 oz.)	Case/ 72		•		