

Material Safety Data Sheet

Cary Blair Vials

ACC# 89659

Section 1 - Chemical Product and Company Identification

MSDS Name: Cary Blair Vials**Catalog Numbers:** 2300547, 61008**Synonyms:** None.**Company Identification:**

Fisher Diagnostics
 Fisher Scientific Company, LLC
 8365 Valley Pike
 Middletown, VA 22645-0307

For information, call: 800-524-0294**Emergency Number:** 800-524-0294**For CHEMTREC assistance, call:** 800-424-9300**For International CHEMTREC assistance, call:** 703-527-3887

Section 2 - Composition, Information on Ingredients

CAS#	Chemical Name	Percent	EINECS/ELINCS
0-00-0	Agar	<1.0	unlisted
143-74-8	Phenol Red	<1.0	205-609-7
367-51-1	Sodium thioglycollate	<1.0	206-696-4
7558-79-4	Sodium phosphate dibasic	<1.0	231-448-7
7647-14-5	Sodium chloride	<1.0	231-598-3
7778-77-0	Dihydrogen potassium phosphate	<1.0	231-913-4
7732-18-5	Water	Balance	231-791-2

Section 3 - Hazards Identification

EMERGENCY OVERVIEW

Appearance: reddish liquid. **Caution!** May cause irritation.**Target Organs:** None.**Potential Health Effects****Eye:** May cause eye irritation.**Skin:** May cause skin irritation.**Ingestion:** May cause irritation of the digestive tract.**Inhalation:** May cause respiratory tract irritation.**Chronic:** No information found.

Section 4 - First Aid Measures