

**Catalog Number: BDH28559.400**

**Catalog Description: Tetrahydrofuran (Inhibitor Free)**

**Grade: HiPerSolv CHROMANORM for HPLC**

**Date Current Specification Set: 2-16-16**

*All batches of this product manufactured after this date are tested to this specification. A new certificate will be published on the web in the unlikely event of specification changes.*


## Characteristics

## Specifications

Assay (by GC-FID) (Corrected for Water)	Min 99.70 %
Water (by KF coulometry)	Max 0.1000 %
Residue after evaporation	Max 0.0005 %
Acidity	Max 0.0005 meq/g
Transmittance at	
230 nm	Min 40 %
240 nm	Min 60 %
250 nm	Min 70 %
260 nm	Min 80 %
280 nm	Min 96 %
300 nm	Min 98 %

Signed on behalf of VWR:

Quality Control Manager



### Additional Information (optional):