



## VWR® BASIC PROTECTION SMS LAB COATS

As per the manufacturer, the below product meets the following criteria:

<b>VWR Catalog Number:</b>	414004-356
<b>Lot Number:</b>	A1191345
<b>Description:</b>	VWR LABCT BP SMS XL BLU CS30
<b>Date of Manufacture:</b>	December 13, 2019
<b>Country of Origin:</b>	CHINA

### Quality System Compliance

Products are manufactured under the ISO 9001:2015 & ISO13485:2003 standard. Products are Inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP. Inspection records are reviewed and signed off by qualified personnel for product release.

### QC Testing

Inspected according to AQL level II 4.0 sampling plan per MIL-STD-105E

<b>-Break Tensile properties(N):</b>	ASTM D5034standard,	MD=91.09	CD=85.19
<b>-Elongation rate (%):</b>	ASTM D5034standard,	MD=79.27	CD=88.78

### Product Specifications

<b>Material:</b>	SMS 50g/m <sup>2</sup> Light Blue (Spunbonded Polypropylene fabric+Melt bown+Spunbonded Polypropylene fabric)
<b>BSE/TSE Statement:</b>	The above product contains no ingredients of animal origin and no material derived from, or exposed to animals affected by or under quarantine for Transmitting, Animal Spongiform Encephalopathy (TSE)/Bovine Spongiform Encephalopathy (BSE)
<b>Non-Pyrogenic:</b>	NA
<b>DNase &amp; RNase Free:</b>	NA
<b>Cytotoxicity:</b>	NA
<b>Latex Free Statement:</b>	The above product is latex free

Signed:



Jamie Ethier  
VP Global Quality  
VWR, Part of Avantor

**Date:** January 04, 2020