



**Lithium Chloride,
Granular
BAKER ANALYZED ACS
Reagent**

Product No. 2370
Specifications current as of: Jun 8
2011

| TEST | SPECIFICATION |
|-------------------------------|---------------|
| Meets A.C.S. Specifications | |
| Assay (LiCl) (by Ag titrn) | min. 99.0 % |
| Insoluble Matter | max. 0.01 % |
| Loss on Drying | max. 1.0 % |
| Titrable Base (meq/g) | max. 0.008 |
| Sulfate (SO ₄) | max. 0.01 % |
| Barium (Ba) | max. 0.002 % |
| Heavy Metals (as Pb) | max. 0.002 % |
| Iron (Fe) | max. 0.001 % |
| Calcium (Ca) | max. 0.01 % |
| Potassium (K) | max. 0.01 % |
| Sodium (Na) | max. 0.20 % |

For Laboratory, Research or Manufacturing Use