


# The Multi-Task Genie™ Family



**THE ORIGINAL**



**TIMER**



**DIGITAL**

U.S. Pat. 4,781,487 / 5,707,861

Vortex-Genie Model	Vortex-Genie® 2	Vortex-Genie® 2T	Digital Vortex-Genie® 2
Supplied	3-inch Platform, Pop-off Cup	3-inch Platform, Pop-off Cup	3-inch Platform, Pop-off Cup
Speed (RPM)	600-3200	600-3200	500-3000
Timer - Touch On: Hands-Free:	N/A	1-60 seconds 1-60 minutes or continuous	1-99 seconds 1-99 minutes or continuous
Dimensions (base) (D x W x H) / Weight	6.5 x 4.8 x 6.5in / 8.8 lb	6.5 x 4.8 x 6.5in / 8.8 lb	6.5 x 4.8 x 6.5in / 8.8 lb
CATALOG NUMBERS			ABI Model
120V (North America)	58815-234	14216-188	89037-464 97060-142

## Multi-Task Genie™ Family Accessories



Fig.	Cat. No.	Description	Fig.	Cat. No.	Description
2,3	58815-230	Model H301 Multiple Sample Starter Set	13	89202-322	TurboMix™ Attachment
--	58815-212	6-inch Platform	14	89202-324	Universal Microtube Holder - (Requires 89202-322)
2	58815-206	Microtube Foam Insert (2) - (Requires 58815-212)	15	13502-388	Large Ampule/Tube Attachment
3	58815-185	Microplate Foam Insert (2) - (Requires 58815-212)	16	13502-390	Small Ampule/Tube Attachment
5-10	58815-224	Model H302 Large Sample Set	17	97060-174	Horizontal 50ml Tube Holder (6 tubes)
5	58815-222	Recessed Platform without Bands	18	97060-176	Horizontal 15ml Tube Holder (12 tubes)
6	58815-188	9-16mm Tube Foam Inserts (2) - (Requires 58815-222)	19	97060-178	Horizontal Microtube Holder (24 tubes)
7	58815-186	14-29mm Tube Foam Inserts (2) - (Requires 58815-222)	20	89401-416	Vertical High Speed 50ml Tube Holder (3 tubes)
8	58815-184	29-37mm Tube Foam Inserts (2) - (Requires 58815-222)	21	89401-418	Vertical High Speed 15ml Tube Holder (6 tubes)
9	13502-386	Blank Foam Inserts (2) - (Requires 58815-222)	22	89401-420	Vertical High Speed Microtube Holder (9 tubes)
6-9	58815-220	Assorted Foam Inserts (4) - (Requires 58815-222)	23	97060-180	Vertical 50ml Tube Holder (6 tubes)
10	58815-218	Elastic Bands (2) - (Requires 58815-222)	24	97060-182	Vertical 15ml Tube Holder (12 tubes)
11	58815-210	Pop-off Cup	25	97060-184	Vertical Microtube Holder (24 tubes)
12	14221-972	3-inch Platform with Rubber Cover	26	97060-186	Vertical (6) 50ml, (6) 15ml, (12) Microtube Holder
--	58815-180	Rubber Cover only for 3-inch Platform	--	89202-338	Bottom Closure with Screws Kit