


| | | | | |
|---------------------------|--|-------------------------|--|---|
| PDS No. 662102/ 662160 | PRODUCT DATA SHEET | | | Page 1 of 1 |
| Revision 04 | 24 Well Cell Culture Multiwell Plate, CELLSTAR® | | |  greiner bio-one |
| | Greiner Item-No. 662102/662160 | | | |
| Valid for Item-No.: | 662102 (sterile) | 662160 (sterile) | | |

| 1. | Description / Specification | |
|-----|-----------------------------|--|
| 1.1 | Description | PS Plate, 24 well, each well with a flattened, raised ring to reduce cross contamination, chimney well, flat bottom, alphanumeric well coding, single position lid with condensation rings, sterile 662102: hydrophobic surface (without treatment) 662160: physical surface treatment |
| 1.2 | Dimensions | <u>Plate</u> : length: 127,76 mm (+/- 0,2 mm) width: 85,48 mm (+/- 0,2 mm) curvature: ≤ 200µm <u>Lid</u> : length: 127,4 mm (+/- 0,2 mm) width: 85,1 mm (+/- 0,2 mm) curvature: +/-300 µm |
| 1.3 | Volume per well | Total volume: 3,3 ml (mathematically calculated) 662160: - working volume: 0,5-1,5 ml - growth area: 1,9 cm ² |
| 1.4 | Material / Resin | <u>Plate</u> : PS (Polystyrene), free of heavy metal <u>Lid</u> : PS (Polystyrene), free of heavy metal |
| 1.5 | Colour | <u>Plate</u> : clear <u>Lid</u> : clear |
| 1.6 | Sterilization | SAL 10 ⁻³ |
| 1.7 | Quality Control | - <u>Raw Material-Control</u> : physical testing - <u>Product-Control</u> : testing of attributive and variable characteristics in accordance with the valid specification |
| 1.8 | Other Information | For single use only |

| 2. | Features | |
|-----|--------------------------|---|
| 2.1 | Basic features | Free of detectable DNase/RNase, human DNA and pyrogens. Contents non-cytotoxic |
| 2.2 | Temperature range | -20°C to +60°C |
| 2.3 | Autoclavability | No |
| 2.4 | Centrifugation, max. RCF | 4800 x g: swinging-bucket rotor |
| 2.5 | Chemical Resistance | See homepage: www.gbo.com/bioscience →Products →Literature →Technical Information →Chemical Resistance of Resins |
| 2.6 | Shelf life | 662102: 5 years after month of production 662160: 4 years after month of production |
| 2.7 | Other Information | - |

| 3. | Packaging | |
|-----|-------------------|--|
| 3.1 | Pieces / Bag | 1 |
| 3.2 | Pieces / Box | 100 |
| 3.3 | Lot-No. | E YY MM XXX (manufacturing facility, year, month, consecutive SAP-No.) |
| 3.4 | Other Information | Certificate of Quality |

| 4. | Other Information |
|----|-------------------|
| | - |

Data Sheet subject to change without notice!

| Prior Issue | Drawn | Approved | Released | CONFIDENTIAL: Information contained in this document or drawing is confidential and proprietary to Greiner Bio-One GmbH. This document may not be reproduced for any reason without written permission from Greiner Bio-One GmbH. All rights of design, invention, and copyright are reserved. |
|--------------------|--------------------------|--------------------------|--------------------------|---|
| Revision 03 | Date 11 December 2009 | Date 14 December 2009 | Date 14 December 2009 | |
| Date 30.11.2005 | Name S. Kaelberer | Name Dr. R. Heller | Name A. Schulz | |