



Item Number M1194
Item Mineral Oil, White, Medium, FCC
CAS Number 8042-47-5

Molecular Formula

Molecular Weight

Carbon number at 5% distillation point..... 17- 25
Specific Gravity.....0.848-0.867
Ultraviolet Absorbance (Polynuclear Hydrocarbons)..... To pass test
Viscosity @ 100°C.....3-11 cSt

MAXIMUM LIMITS

Lead..... 1 mg/kg
Readily Carbonizable Substances..... To pass test

Spectrum Chemicals & Laboratory Products

Corporate Headquarters: East Coast Facility:

14422 S. San Pedro St. 755 Jersey Ave.

Gardena, CA 90248 New Brunswick, NJ 08901

(800) 813-1514 • (310) 516-8000 • Fax (310) 516-9843