



Chloroform AR (ACS)

Product No. 4440
Specifications current as of: Jun 8 2011

TEST	SPECIFICATION
Meets A.C.S. Specifications	
Appearance (clear, colorless liquid)	Passes Test
Assay (CHCl ₃) (by GC)	min. 99.8 %
Acetone, Aldehydes (as (CH ₃) ₂ CO)	Passes Test
Acid and Chloride	Passes Test
Color (APHA)	max. 10
Free Chlorine	Passes Test
Lead (Pb)	max. 0.05 ppm
Residue after Evaporation	max. 0.001 %
Substances Darkened by H ₂ SO ₄	Passes Test
Suitability for Use in Dithizone Tests	Passes Test
For Laboratory, Research or Manufacturing Use	