



Material	38332-700
Material description	Ethyl alcohol, anhydrous, denatured
CAS Number	64-17-5

## Characteristics

## Specifications

Assay (Ethanol) (V/V)	Min. 84.46 %
Color (Pt-Co)	Max. 5 APHA
Density (20°C)	Max. 0.7922 g/ml
Acidity (as CH <sub>3</sub> COOH) (W/V)	Max. 0.003 %
Ethyl acetate (V/V)	0.75 --> 0.95 %
Methanol (V/V)	12.95 --> 14.45 %
Water (V/V)	Max. 0.14 %

## Signature

We certify that this batch conforms to the specifications listed above.

This document has been produced electronically and is valid without a signature.

Michel Boisvert, Head of the Laboratory - Lachine  
VWR Anachemia Canada, 255 rue Norman, Lachine, H8R 1A3,  
QC, Canada

For Professional use in Laboratory or Manufacturing. Not for use as an Active Pharmaceutical Ingredient or Food or Animal Feed. Suitability and intended use of the product remains the responsibility of the user.