



Item Number SO168
Item Sodium Gluconate, FCC
CAS Number 527-07-1
Molecular Formula $C_6H_{11}NaO_7$
Molecular Weight 218.14
Assay ($C_6H_{11}NaO_7$)..... 98.0-102.0%

MAXIMUM LIMITS

Lead..... 2 mg/kg
Reducing Substances..... 0.5%

Spectrum Chemicals & Laboratory Products

Corporate Headquarters: **East Coast Facility:**
14422 S. San Pedro St. 755 Jersey Ave.
Gardena, CA 90248 New Brunswick, NJ 08901
(800) 813-1514 • (310) 516-8000 • Fax (310) 516-9843