



Item Number A1293
Item Ammonia Solution, FCC
CAS Number 1336-21-6
Molecular Formula NH₄OH
Molecular Weight 35.05
Assay (NH₃; by Wt.)..... 27.0-30.0%

MAXIMUM LIMITS

Lead..... 0.5 mg/kg
Nonvolatile Residue..... 0.02%
Readily Oxidizable Substances..... To pass test

Spectrum Chemicals & Laboratory Products

Corporate Headquarters: **East Coast Facility:**
14422 S. San Pedro St. 755 Jersey Ave.
Gardena, CA 90248 New Brunswick, NJ 08901
(800) 813-1514 • (310) 516-8000 • Fax (310) 516-9843