

Product Name:	VWR-HC-36S
Description:	<p>VWR Catalog No: 76579-118  VWR Plus Series Solid Door Laboratory Refrigerator  36 Cu. Ft. capacity  Two swing solid doors, self-closing, right hinged, not reversible  Microprocessor Temperature Controller  Adjustable operating temperature range: 1°-10°C  Digital Temperature Display and Alarm Module  Sample simulated digital temperature display; .1° Resolution  Adjustable temperature display in Celsius or Fahrenheit  Audible and visual high/low temperature alarms  2 Probes; 1 in air and 1 in sample bottle  Sensor error alarm  Power failure alarm  Resettable 1-point Min/Max History  Remote alarm contacts  Alarm validation (audible, visual, and remote alarm)  Alarm mute with ring back  Eight adjustable shelves (adjustable in ½" increments)  Keyed door lock  Forced draft air circulation  Magnetic door gasket for positive seal  Heavy duty door handle  LED interior lights are safety shielded and switch controlled  Swivel Casters (3 ½")  Probe access port (3/8")  White powder-coated steel exterior and interior  Hydrocarbon, natural refrigerant (R290)  Two year parts and labor warranty, plus an additional three compressor parts warranty  Exterior dimensions: 39 5/8" W x 34 3/4" D x 81 3/4" H with casters; 78 1/4" H without casters  Cycle defrost  115V, 60 Hz, 4.5 Amps, 1/4 HP  ETL listed  Shipping weight: 360 lbs.</p>