

Certificate of Quality



Part Number	MFLX31201-33
Lot Number	8022105
Part Description	ADAPTR 1.5IN CLAMP 1/4IN ID
Expiration Date	02AUG2024
Country of Origin	USA

Inspection & Testing

- It is the responsibility of the user to conduct tests that are deemed necessary to determine the suitability of this product for any particular use.

Product Specifications & Compliance

- Compound: Polypropylene with fused TPE gasket
- This product is Animal-Derived Component Free.
- Material meets USP Class VI requirements.
- Date of expiration (DOE) is set from 5 years from the date of manufacturing (DOM).
- It is the responsibility of the user to conduct tests that are deemed necessary to determine the suitability of this product for any particular use.
- This product is not approved for or intended for, and should not be used for medical, clinical, surgical, or other patient oriented applications.

Storage

- Storage in original packaging, at ambient conditions and away from direct sunlight, is recommended.

Certificate Approval

- This certificate was issued and approved by the Avantor Fluid Handling Quality and Regulatory Departments and is valid without signature.