



Item # _____

Job _____



Super Erecta Shelf® LAB ANIMAL RESEARCH (LAR) CARTS

LAR Feed Carts are designed to store and transport animal feed. Stainless steel finish provides superior corrosion resistance, while compact footprint makes it easy to maneuver in and around tight lab spaces. Clear totes provide easy identification of tote contents.

Lab Security Cart provides secure storage for animal cage, lab instrument, or hazardous chemical transport. Solid stainless steel shelf with raised ship's edge can be used as a work surface and can contain unwanted spills.

- **Super Erecta Shelf® Construction** makes units easy to assemble with absolute rigidity, and makes relocation of shelves quick and simple.
- **Wire Shelves:** Bright, modern, sanitary appearance. Open construction minimizes dust accumulation, maximizes visibility and air circulation.
- **Casters:** Light duty, 4" (102mm) resilient rubber casters are included with Catalog Number LAR1X. While 5" (127mm) polyurethane tread casters are included with Catalog Number LAR2X and SECMLAB.
- **Clear Totes:** Each cart comes with either one (LAR1) or two (LAR2) clear totes for clean and secure animal feed storage. A snap-on cover is included with each 12" (305mm) deep tote.
- **Cover Holder:** A convenient tote cover holder comes standard with each cart. Attaches easily to wire shelf and takes up minimal space.
- **Donut Bumpers** included with casters.
- **Height:** 38" (965mm) floor to handle with 4" (102mm) casters; 39" (990mm) floor to handle with 5" (127mm) casters.
- **Accessories:** A wide variety of Super Erecta-based accessories to choose from.



LAR Feed Cart
(One Tote Capacity)
Cat. No. LAR1



LAR Feed Cart
(Two Tote Capacity)
Cat. No. LAR2



Lab Security Cart
(Shown with
Optional Accessories)
Cat. No. SECMLAB



InterMetro Industries Corporation
North Washington Street
Wilkes-Barre, PA 18705
www.metro.com





Super Erecta Shelf® LAB ANIMAL RESEARCH (LAR) CARTS

Specifications

- All cart components, including wire shelves and handles, are stainless steel, unless otherwise noted. Split sleeves and donut bumpers are plastic. Caster and tote information is noted below.

Dimensions

Cat. No.	Width		Length		Height		Cart Load Rating ¹		Description	Approx. Pkd. Wt.	
	(in.)	(mm)	(in.)	(mm)	(in.)	(mm)	(lbs.)	(kg)		(lbs.)	(kg)
LAR1	18	457	24	610	38	965	375	171	Light-duty feed cart One tote capacity/shelf 4" (102mm) resilient casters ²	42	18.7
LAR2	24	610	36	914	39	990	600	273	Heavy-duty feed cart Two tote capacity/shelf 5" (127mm) polyurethane casters ²	65	29.2
SECMLAB ³	24	610	30	760	39	990	600	273	Stainless steel security cart with flip-up door 5" (127mm) polyurethane casters	85	38.6

¹Cart load rating listed is based on an evenly distributed load.

²Cat. No. LAR1 features four L4D series swivel casters with resilient rubber tread. Cat. No. LAR2 features four 5MP series swivel casters with polyurethane tread.

³LAB Security Cart features two stainless steel wire shelves, one stainless steel security module, and one stainless steel solid shelf. Cart comes standard with two brake casters (Cat. No. 5MPB) and two swivel casters (Cat. No. 5MP).

Lab Animal Research Accessories

Durable, clear totes provide clean and secure storage for animal feed. Totes are injection-molded from clear, polypropylene material, are stackable, and can be combined with clear snap-on covers and cardholders for security and identification of contents.



Clear Tote Boxes
(Shown with Snap-on Cover and Cardholder Accessories)

Clear Tote Box

Cat. No.	Clear Polypropylene	Outside Dimensions				Inside Dimensions	
		Length x Width x Height				Bottom	
		Top (in.)	(mm)	Bottom (in.)	(mm)	(in.)	(mm)
TB93120CLR		22.5x17.5x12	572x445x305	21x16	533x406	20.2x15.1x11.5	513x384x292

Note 1: All outside dimensions have an overall tolerance of (+/-) .050". All inside dimensions have an overall tolerance of (+/-) .100".
Note 2: Tote not recommended for freezer environments (temperatures lower than -20°F [-29°C]).

Clear Snap-on Cover

Cat. No.	Length		Outside Dimensions		Height	
	(in.)	(mm)	(in.)	(mm)	(in.)	(mm)
CO93000CLR	22.5	572	17.5	445	1.25	32

Tote Box Cardholders

Cat. No.	Material	Length		Inside Dimensions		Height	
		(in.)	(mm)	(in.)	(mm)	(in.)	(mm)
OP2500CLR	Vinyl	1	25.4	8	203	8	203
OP2501CLR	Vinyl	4.5	114.3	8	203	8	203
OP2535CLR	Vinyl	2 ¹¹ / ₁₆	60.6	5	127	5	127

Cover Holder

Stainless steel wire cover holder provides a convenient storage option for tote box covers when not in use. Accessory easily mounts on any Super Erecta® based wire shelf.



Cover Holder
(Shown mounted on Cart with one Tote Cover)

Cat. No.	Width		Length		Depth		Approx. Carton Wt.	
	(in.)	(mm)	(in.)	(mm)	(in.)	(mm)	(lbs.)	(kg)
LAR-CH	18 ¹ / ₂	470	18 ¹ / ₄	464	4	101	6	2.7

All Metro Catalog Sheets are available on our Web Site: www.metro.com



InterMetro Industries Corporation

North Washington Street, Wilkes-Barre, PA 18705
Phone: 570-825-2741 • Fax: 570-825-2852

For Product Information Call: 1-800-433-2232

L03-239

Printed in U.S.A. 3/02

Information and specifications are subject to change without notice. Please confirm at time of order.

Copyright © 2002 InterMetro Industries Corp.

