
Product Specifications

VWR Part Number: IR2430
ImmunoReagents, Inc. Part Number: SP-004-VX10

Product: NORMAL GOAT SERUM 10ML

Concentration: 60.0 mg/ml (Bradford, IgG standards)

Amount: 10 ml

Conjugate: Unconjugated

Form: Lyophilized

Purification: Lipid extracted and dialyzed against 10 mM Sodium Phosphate, 0.15 M Sodium Chloride, pH 7.2

Purity: N/A

Host: Normal

Immunogen: N/A

Buffer: 10 mM Sodium Phosphate, 0.15 M Sodium Chloride, pH 7.2

Storage: Store lyophilized material at 2-8 °C. For long term storage after reconstitution, prepare small aliquots and store at -20 °C. For storage at 2-8 °C, add a preservative to prevent growth of bacteria.

Shelf Life: 1 year from date of receipt. Prepare working dilution prior to use and then discard.

Specificity: Non-immune

Cross Reactivity: N/A

Country of Origin: Normal Goat Serum was obtained from healthy animals of US origin and under the care of a registered veterinarian.

Applications: This product is used as a blocking reagent or control for most immunoassay applications.

Disclaimer: For in vitro Laboratory Use Only. Not for diagnostic or therapeutic use. Not for human or animal consumption. Suggested applications of our products are not recommendations to use our products in violation of any patent or as a license under any patent of ImmunoReagents, Inc.

Product may not be resold or modified for resale without prior written approval of ImmunoReagents, Inc.