



VWR ANALYTICAL

Product Specification

Catalog Number: BDH1131-4LP, BDH1131-4LP, BDH1131-5GL, BDH1131-204L, BDH1131-4LP, BDH1131-4LP, BDH1131-5GL, BDH1131-204L,

Catalog Description: Isopropyl Alcohol 70% Water 30% (v/v)

Grade:

Date Current Specification Set: June 26 2008

All batches of this product manufactured after this date are tested to this specification. A new certificate will be published on the web in the unlikely event of specification changes.

Characteristics	Specifications
Appearance:	Clear, colorless liquid
Isopropyl Alcohol:	68.0-72.0%
Water:	28.0-32.0%
Specific gravity @ 20°C:	0.863-0.883
Color (APHA):	Less than 10

Signed on behalf of VWR:

Neal Fox
Quality Control Supervisor