



Item Number B1196
Item Butylated Hydroxytoluene, Granular, NF
CAS Number 128-37-0
Molecular Formula C₁₅H₂₄O
Molecular Weight 220.36

Assay..... Min. 99.0%
Congealing Temperature..... Min. 69.2°C

MAXIMUM LIMITS

Residue on Ignition..... 0.002%
Heavy Metals (as Pb)..... 0.001%
Related compounds..... To pass test
Residual Solvents..... To pass test

Spectrum Chemicals & Laboratory Products

Corporate Headquarters: **East Coast Facility:**
14422 S. San Pedro St. 755 Jersey Ave.
Gardena, CA 90248 New Brunswick, NJ 08901
(800) 813-1514 • (310) 516-8000 • Fax (310) 516-9843