



VWR ANALYTICAL

# Product Specification

Material	BDH7999-4
Material description	Isopropyl alcohol 70% v/v
CAS Number	67-63-0
Molecular formula	C <sub>3</sub> H <sub>8</sub> O
Molecular mass	60,10
Additional information	

## Characteristics

## Specifications

Assay (V/V)	69,0 --> 71,0 %
Appearance	Clear liquid
Specific gravity	0,873 --> 0,883
Prepared to formulation on file	Confirmed

## Signature

We certify that this batch conforms to the specifications listed above.

This document has been produced electronically and is valid without a signature.

Alicia Bourdeau, Head of the Laboratory - Rouses Point  
VWR Chemicals, LLC, 3 Lincoln Blvd, Rouses Point, NY  
12979-1045, USA

## Additional information

The weights and/or volumes used to prepare this product are NIST traceable. All balances used in the preparation of the product are calibrated weekly against NIST traceable weights. The balances are maintained and serviced on a regular basis by an outside certified company. All volumetric glassware used is NIST traceable and certified as meeting Class A specifications. Unless otherwise indicated above, all chemicals used in the preparation of this product are Reagent ACS grade.