



# SAFETY DATA SHEET

according to the (US) Hazard Communication Standard (29 CFR 1910.1200)

Date of issue: 04/28/2014

Version 1.0

## SECTION 1. Identification

### Product identifier

Product number PX1840  
Product name *iso*-Propyl Alcohol <br/>in Water 70% (v/v)

### Relevant identified uses of the substance or mixture and uses advised against

Identified uses Reagent for analysis

### Details of the supplier of the safety data sheet

Company EMD Millipore Corporation | 290 Concord Road, Billerica, MA 01821,  
United States of America | General Inquiries: +1-978-715-4321 |  
Monday to Friday, 9:00 AM to 4:00 PM Eastern Time (GMT-5)

Emergency telephone 800-424-9300 CHEMTREC (USA)  
+1-703-527-3887 CHEMTREC (International)  
24 Hours/day; 7 Days/week

## SECTION 2. Hazards identification

### GHS Classification

Flammable liquid, Category 2, H225  
Eye irritation, Category 2, H319  
Specific target organ systemic toxicity - single exposure, Category 3, H336  
For the full text of the H-Statements mentioned in this Section, see Section 16.

### GHS-Labeling

#### Hazard pictograms



#### Signal Word

Danger

#### Hazard Statements

H225 Highly flammable liquid and vapor.  
H319 Causes serious eye irritation.  
H336 May cause drowsiness or dizziness.

# SAFETY DATA SHEET

according to the (US) Hazard Communication Standard (29 CFR 1910.1200)

Product number

PX1840

Version 1.0

Product name

*iso*-Propyl Alcohol <br/>in Water 70% (v/v)

---

## *Precautionary Statements*

P210 Keep away from heat/sparks/open flames/hot surfaces. - No smoking.

P233 Keep container tightly closed.

P305 + P351 + P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

## OSHA Hazards

This material is considered hazardous by the OSHA Hazard Communication Standard (29 CFR 1910.1200). This information is based on 29 CFR 1910.1200 criteria prior to adoption of the GHS and may deviate from the GHS information.

## Other hazards

None known.

---

## SECTION 3. Composition/information on ingredients

Chemical nature

Aqueous solution

### Hazardous ingredients

*Chemical Name (Concentration)*

CAS-No.

*2-Propanol (>= 70 % - < 90 % )*

67-63-0

Exact percentages are being withheld as a trade secret.

---

## SECTION 4. First aid measures

### Description of first-aid measures

#### *Inhalation*

After inhalation: fresh air. Consult doctor if feeling unwell.

#### *Skin contact*

After skin contact: wash off with plenty of water. Remove contaminated clothing.

#### *Eye contact*

After eye contact: rinse out with plenty of water. Immediately call in ophthalmologist.

#### *Ingestion*

After swallowing: caution if victim vomits. Risk of aspiration! Keep airways free. Pulmonary failure possible after aspiration of vomit. Call a physician immediately.

Never give anything by mouth to an unconscious person.

### Most important symptoms and effects, both acute and delayed

irritant effects, respiratory paralysis, Drowsiness, Dizziness, Unconsciousness, narcosis, inebriation, Headache, drowsiness, Coma

### Indication of any immediate medical attention and special treatment needed

No information available.

---

## SECTION 5. Fire-fighting measures

### Extinguishing media

---

## SAFETY DATA SHEET

according to the (US) Hazard Communication Standard (29 CFR 1910.1200)

Product number PX1840  
Product name *iso*-Propyl Alcohol <br/>in Water 70% (v/v)

Version 1.0

### *Suitable extinguishing media*

Carbon dioxide (CO<sub>2</sub>), Foam, Dry powder

### *Unsuitable extinguishing media*

For this substance/mixture no limitations of extinguishing agents are given.

### **Special hazards arising from the substance or mixture**

Combustible material

Vapors are heavier than air and may spread along floors.

Forms explosive mixtures with air at ambient temperatures.

Pay attention to flashback.

Development of hazardous combustion gases or vapors possible in the event of fire.

### **Advice for firefighters**

#### *Special protective equipment for fire-fighters*

In the event of fire, wear self-contained breathing apparatus.

#### *Further information*

Cool closed containers exposed to fire with water spray. Suppress (knock down) gases/vapors/mists with a water spray jet. Prevent fire extinguishing water from contaminating surface water or the ground water system.

---

## **SECTION 6. Accidental release measures**

### **Personal precautions, protective equipment and emergency procedures**

Advice for non-emergency personnel: Avoid substance contact. Do not breathe vapors, aerosols. Ensure adequate ventilation. Evacuate the danger area, observe emergency procedures, consult an expert.

Advice for emergency responders: Protective equipment see section 8.

### **Environmental precautions**

Do not empty into drains. Risk of explosion.

### **Methods and materials for containment and cleaning up**

Cover drains. Collect, bind, and pump off spills.

Observe possible material restrictions (see sections 7 and 10).

Take up with liquid-absorbent material (e.g. Chemisorb®). Dispose of properly. Clean up affected area.

---

## **SECTION 7. Handling and storage**

### **Precautions for safe handling**

Work under hood. Do not inhale substance/mixture. Avoid generation of vapors/aerosols.

Observe label precautions.

#### *Advice on protection against fire and explosion*

Keep away from open flames, hot surfaces and sources of ignition. Take precautionary measures against static discharge.

### **Conditions for safe storage, including any incompatibilities**

Requirements for storage areas and containers Keep away from heat and sources of ignition.

Keep container tightly closed in a dry and well-ventilated place.

# SAFETY DATA SHEET

according to the (US) Hazard Communication Standard (29 CFR 1910.1200)

Product number PX1840  
Product name *iso*-Propyl Alcohol <br/>in Water 70% (v/v)

Version 1.0

No declaration

## SECTION 8. Exposure controls/personal protection

### Exposure limit(s)

#### Ingredients

Basis	Value	Threshold limits	Remarks
<i>2-Propanol 67-63-0</i>			
ACGIH	Time Weighted Average (TWA):	200 ppm	
	Short Term Exposure Limit (STEL):	400 ppm	
NIOSH/GUIDE	Recommended exposure limit (REL):	400 ppm 980 mg/m <sup>3</sup>	
	Short Term Exposure Limit (STEL):	500 ppm 1,225 mg/m <sup>3</sup>	
OSHA_TRANS	PEL:	400 ppm 980 mg/m <sup>3</sup>	
Z1A	Short Term Exposure Limit (STEL):	500 ppm 1,225 mg/m <sup>3</sup>	
	Time Weighted Average (TWA):	400 ppm 980 mg/m <sup>3</sup>	

### Engineering measures

Technical measures and appropriate working operations should be given priority over the use of personal protective equipment.

### Individual protection measures

Protective clothing should be selected specifically for the workplace, depending on concentration and quantity of the hazardous substances handled. The chemical resistance of the protective equipment should be inquired at the respective supplier.

### Hygiene measures

Immediately change contaminated clothing. Apply skin- protective barrier cream. Wash hands and face after working with substance.

### Eye/face protection

Safety glasses

### Hand protection

Chemical-resistant, impervious gloves complying with an approved standard should be worn at all times when handling chemical products if a risk assessment indicates this is necessary.

### Other protective equipment:

Flame retardant antistatic protective clothing

### Respiratory protection

required when vapors/aerosols are generated.

Use a properly fitted, air-purifying or air-fed respirator complying with an approved standard if a risk assessment indicates this is necessary. Respirator selection must be based on known or anticipated exposure levels, the hazards of the product and the safe working limits of the selected respirator.

# SAFETY DATA SHEET

according to the (US) Hazard Communication Standard (29 CFR 1910.1200)

Product number

PX1840

Version 1.0

Product name

*iso*-Propyl Alcohol <br/>in Water 70% (v/v)

---

## SECTION 9. Physical and chemical properties

Physical state	liquid
Color	colorless
Odor	characteristic
Odor Threshold	No information available.
pH	at 68 °F (20 °C) neutral
Melting point	No information available.
Boiling point	No information available.
Flash point	ca. 64 °F (18 °C) Method: DIN 51755 Part 1
Evaporation rate	No information available.
Flammability (solid, gas)	No information available.
Lower explosion limit	2 %(V) (2-Propanol)
Upper explosion limit	13.4 %(V) (2-Propanol)
Vapor pressure	ca.42 hPa at 68 °F (20 °C)
Relative vapor density	No information available.
Density	0.85 g/cm <sup>3</sup> at 68 °F (20 °C)
Relative density	No information available.
Water solubility	at 68 °F (20 °C) soluble
Partition coefficient: n-octanol/water	log Pow: 0.05 OECD Test Guideline 107 (2-Propanol) Bioaccumulation is not expected.
Autoignition temperature	No information available.
Decomposition temperature	No information available.

## SAFETY DATA SHEET

according to the (US) Hazard Communication Standard (29 CFR 1910.1200)

Product number

PX1840

Version 1.0

Product name

*iso*-Propyl Alcohol <br/>in Water 70% (v/v)

---

Viscosity, dynamic No information available.

Explosive properties No information available.

Oxidizing properties No information available.

Ignition temperature ca. 797 °F (425 °C)

---

### SECTION 10. Stability and reactivity

#### Reactivity

Vapors may form explosive mixture with air.

#### Chemical stability

The product is chemically stable under standard ambient conditions (room temperature) .

#### Possibility of hazardous reactions

Risk of ignition or formation of inflammable gases or vapors with:

Alkali metals, Alkaline earth metals, powdered aluminum

Exothermic reaction with:

Oxidizing agents, Aldehydes, Amines, fuming sulfuric acid, Iron, Nitric acid

Risk of explosion with:

chlorates, Phosgene, organic nitro compounds, hydrogen peroxide, nitrogen oxides

#### Conditions to avoid

Warming.

A range from approx. 15 Kelvin below the flash point is to be rated as critical.

#### Incompatible materials

various plastics, rubber

#### Hazardous decomposition products

no information available

---

### SECTION 11. Toxicological information

#### Information on toxicological effects

##### *Likely route of exposure*

Inhalation, Eye contact, Skin contact

##### *Target Organs*

Eyes

Skin

Respiratory system

##### *Skin irritation*

Repeated or prolonged exposure may cause skin irritation and dermatitis, due to degreasing properties of the product.

##### *Eye irritation*

Mixture causes serious eye irritation.

---

# SAFETY DATA SHEET

according to the (US) Hazard Communication Standard (29 CFR 1910.1200)

Product number

PX1840

Version 1.0

Product name

*iso*-Propyl Alcohol <br/>in Water 70% (v/v)

---

### *Specific target organ systemic toxicity - single exposure*

Mixture may cause drowsiness or dizziness.

### *Specific target organ systemic toxicity - repeated exposure*

The substance or mixture is not classified as specific target organ toxicant, repeated exposure.

### *Aspiration hazard*

Regarding the available data the classification criteria are not fulfilled.

## **Carcinogenicity**

IARC	No ingredient of this product present at levels greater than or equal to 0.1% is identified as probable, possible or confirmed human carcinogen by IARC.
OSHA	No ingredient of this product present at levels greater than or equal to 0.1% is identified as a carcinogen or potential carcinogen by OSHA.
NTP	No ingredient of this product present at levels greater than or equal to 0.1% is identified as a known or anticipated carcinogen by NTP.
ACGIH	No ingredient of this product present at levels greater than or equal to 0.1% is identified as a carcinogen or potential carcinogen by ACGIH.

## **Further information**

Systemic effects:

After absorption:

Headache, Dizziness, inebriation, Unconsciousness, narcosis

After uptake of large quantities:

respiratory paralysis, Coma

Further data:

Other dangerous properties can not be excluded.

Handle in accordance with good industrial hygiene and safety practice.

## **Ingredients**

### *2-Propanol*

#### *Acute oral toxicity*

LDLO human: 3,570 mg/kg (RTECS)

LD50 rat: 5,045 mg/kg (RTECS)

#### *Acute inhalation toxicity*

LC50 rat: 46.5 mg/l; 4 h (External MSDS)

#### *Acute dermal toxicity*

LD50 rabbit: 12,800 mg/kg (RTECS)

#### *Eye irritation*

rabbit

Result: Eye irritation  
(RTECS)

# SAFETY DATA SHEET

according to the (US) Hazard Communication Standard (29 CFR 1910.1200)

Product number

PX1840

Version 1.0

Product name

*iso*-Propyl Alcohol <br/>in Water 70% (v/v)

---

#### *Sensitization*

Sensitization test: guinea pig

Result: negative

(IUCLID)

#### *Germ cell mutagenicity*

##### *Genotoxicity in vivo*

Mutagenicity (mammal cell test): micronucleus.

Result: negative

(IUCLID)

#### *Genotoxicity in vitro*

Ames test

Result: negative

(IUCLID)

#### *Carcinogenicity*

Did not show carcinogenic effects in animal experiments. (IUCLID)

#### *Reproductive toxicity*

No impairment of reproductive performance in animal experiments. (IUCLID)

#### *Teratogenicity*

Did not show teratogenic effects in animal experiments. (IUCLID)

---

## SECTION 12. Ecological information

### Ecotoxicity

No information available.

### Persistence and degradability

No information available.

### Bioaccumulative potential

*Partition coefficient: n-octanol/water*

log Pow: 0.05

OECD Test Guideline 107

(2-Propanol) Bioaccumulation is not expected.

### Mobility in soil

No information available.

## Ingredients

### *2-Propanol*

#### *Toxicity to fish*

LC50 *Lepomis macrochirus* (Bluegill sunfish): 1,400 mg/l; 96 h (ECOTOX Database)

#### *Toxicity to daphnia and other aquatic invertebrates*

EC5 *E.sulcatum*: 4,930 mg/l; 72 h (maximum permissible toxic concentration) (Lit.)

EC50 *Daphnia magna* (Water flea): 13,299 mg/l; 48 h (IUCLID)

#### *Toxicity to algae*

IC50 *Desmodesmus subspicatus* (green algae): > 1,000 mg/l; 72 h (IUCLID)

#### *Toxicity to bacteria*

EC5 *Pseudomonas putida*: 1,050 mg/l; 16 h (Lit.)

---



# SAFETY DATA SHEET

according to the (US) Hazard Communication Standard (29 CFR 1910.1200)

Product number

PX1840

Version 1.0

Product name

*iso*-Propyl Alcohol <br/>in Water 70% (v/v)

---

*Biodegradability*

95 %; 21 d

OECD Test Guideline 301E

Readily biodegradable.

*Theoretical oxygen demand (ThOD)*

2,400 mg/g

(Lit.)

*Ratio BOD/ThBOD*

BOD5 49 %

(IUCLID)

*Ratio COD/ThBOD*

96 %

(Lit.)

Substance does not meet the criteria for PBT or vPvB according to Regulation (EC) No 1907/2006, Annex XIII.

---

## SECTION 13. Disposal considerations

The information presented only applies to the material as supplied. The identification based on characteristic(s) or listing may not apply if the material has been used or otherwise contaminated. It is the responsibility of the waste generator to determine the toxicity and physical properties of the material generated to determine the proper waste identification and disposal methods in compliance with applicable regulations. Disposal should be in accordance with applicable regional, national and local laws and regulations.

---

## SECTION 14. Transport information

### Land transport (DOT)

UN number	UN 1219
Proper shipping name	ISOPROPANOL SOLUTION
Class	3
Packing group	II
Environmentally hazardous	--

### Air transport (IATA)

UN number	UN 1219
Proper shipping name	ISOPROPANOL SOLUTION
Class	3
Packing group	II
Environmentally hazardous	--
Special precautions for user	no

### Sea transport (IMDG)

UN number	UN 1219
Proper shipping name	ISOPROPANOL SOLUTION
Class	3
Packing group	II
Environmentally hazardous	--

SAFETY DATA SHEET

according to the (US) Hazard Communication Standard (29 CFR 1910.1200)

Product number

PX1840

Version 1.0

Product name

*iso*-Propyl Alcohol <br/>in Water 70% (v/v)

---

**Special precautions for user**

yes

EmS

F-E S-D

---

**SECTION 15. Regulatory information**

**United States of America**

**OSHA Hazards**

Flammable Liquid

Target organ effects

Eye irritant

This information is based on 29 CFR 1910.1200 criteria prior to adoption of the GHS, and may deviate from the GHS information on the label and in section 2.

**SARA 311/312 Hazards**

Fire Hazard

Acute Health Hazard

Chronic Health Hazard

**SARA 313**

The following components are subject to reporting levels established by SARA Title III, Section 313:

*Ingredients*

2-Propanol

67-63-0

70 %

**SARA 302**

SARA 302: No chemicals in this material are subject to the reporting requirements of SARA Title III, Section 302.

**Clean Water Act**

This product does not contain any Hazardous Substances listed under the U.S. CleanWater Act, Section 311, Table 116.4A.

This product does not contain any Hazardous Chemicals listed under the U.S. CleanWater Act, Section 311, Table 117.3.

**DEA List I**

Not listed

**DEA List II**

Not listed

**US State Regulations**

**Massachusetts Right To Know**

*Ingredients*

2-Propanol

**Pennsylvania Right To Know**

*Ingredients*

2-Propanol

water

**New Jersey Right To Know**

*Ingredients*

# SAFETY DATA SHEET

according to the (US) Hazard Communication Standard (29 CFR 1910.1200)

Product number

PX1840

Version 1.0

Product name

*iso*-Propyl Alcohol <br/>in Water 70% (v/v)

---

2-Propanol  
water

### California Prop 65 Components

This product does not contain any chemicals known to the State of California to cause cancer, birth, or any other reproductive defects.

### Notification status

TSCA: All components of the product are listed in the TSCA-inventory.

DSL: All components of this product are on the Canadian DSL.

---

## SECTION 16. Other information

### Training advice

Provide adequate information, instruction and training for operators.

### Full text of H-Statements referred to under sections 2 and 3.

H225 Highly flammable liquid and vapor.  
H319 Causes serious eye irritation.  
H336 May cause drowsiness or dizziness.

### Key or legend to abbreviations and acronyms used in the safety data sheet

Used abbreviations and acronyms can be looked up at [www.wikipedia.org](http://www.wikipedia.org).

Date of issue: 04/28/2014

---

The information contained herein is based on the present state of our knowledge. It characterizes the product with regard to appropriate safety precautions. It does not represent a warranty of any product properties and we assume no liability for any loss or injury which may result from the use of this information. Users should conduct their own investigations to determine the suitability of the information.

*All rights reserved. Millipore and the "M" Mark are registered trademarks of Merck KGaA, Darmstadt, Germany.*