



Material	67482-540
Material description	Petroleum ether
Grade	ACS
CAS Number	8032-32-4
Molecular formula	—

Characteristics

Specifications

Acidity	Passes test
Boiling range (35-60°C)	Passes test
Color (Pt-Co)	Max. 10 APHA
Residue after evaporation	Max. 0,001 %
Meets MIL-PRF-680 Type 1 specifications	Conforms

Signature

We certify that this batch conforms to the specifications listed above.

This document has been produced electronically and is valid without a signature.

Michel Boisvert, Head of the Laboratory - Lachine
VWR Anachemia Canada, 255 rue Norman, Lachine, H8R 1A3,
QC, Canada