

Product Specification

| Test | Specification |
|--|---------------|
| Assay (NaCl) (by Ag titrn) | ≥ 99.0 % |
| pH of 5% Solution at 25°C | 5.0 – 9.0 |
| ACS – Insoluble Matter | ≤ 0.005 % |
| Iodide (I) | ≤ 0.002 % |
| Bromide (Br) | ≤ 0.01 % |
| Chlorate and Nitrate (as NO ₃) | ≤ 0.003 % |
| ACS – Phosphate (PO ₄) | ≤ 5 ppm |
| Sulfate (SO ₄) | ≤ 0.004 % |
| Barium (Ba) | Passes Test |
| ACS – Heavy Metals (as Pb) | ≤ 5 ppm |
| Iron (Fe) | ≤ 2 ppm |
| Calcium (Ca) | ≤ 0.002 % |
| Magnesium (Mg) | ≤ 0.001 % |
| Potassium (K) | ≤ 0.005 % |

For Laboratory, Research, or Manufacturing Use
Meets Reagent Specifications for testing USP/NF monographs

Packaging Site: Paris Mfg Ctr & DC