

Product Specification

Test	Specification
Appearance (Clear, colorless liquid)	Passes Test
Normality	9.9 – 10.1
Sodium Carbonate (Na ₂ CO ₃)	≤ 2.0 %
Chloride (Cl)	≤ 5 ppm
ACS – Heavy Metals (as Pb)	≤ 1 ppm
Trace Impurities – Iron (Fe)	≤ 1 ppm

For Laboratory, Research, or Manufacturing Use
Standardization at 25°C traceable to NIST Standard Reference Material.

SRM No

IMPORTANT: Material will freeze if stored below 41°F (5°C).

Storage Condition: IMPORTANT: Material will freeze if stored below 41 °F (5°C). Protect from air to avoid absorption of carbon dioxide.

Packaging Site: Phillipsburg Mfg Ctr & DC