



order on
VWR.COM

EcoTek[®] FINALLY. BIODEGRADABILITY WITHOUT PERFORMANCE LOSS.



Powerform S6

AQL
1.5

Specifications

Material:	Nitrile
Color:	Teal
Inner Surface:	Powder-Free
Outer Surface:	Fully Textured
Length:	9.5"
Finger Thickness:	5.8 mil
Palm Thickness:	5.0 mil
Dispenser Quantity:	100 Gloves
Case Quantity:	1000 Gloves

Mfr. No.	Size
76051-116	S: 6.5-7
76051-118	M: 7.5-8
76051-120	L: 8.5-9
76051-122	XL: 9.5-10
76051-124	XXL: 10.5-11



Stellar S6

AQL
1.5

Specifications

Material:	Nitrile
Color:	Black
Inner Surface:	Powder-Free
Outer Surface:	Fully Textured
Length:	9.5"
Finger Thickness:	5.6 mil
Palm Thickness:	5.0 mil
Dispenser Quantity:	100 Gloves
Case Quantity:	1000 Gloves

Mfr. No.	Size
76051-106	S: 6.5-7
76051-108	M: 7.5-8
76051-110	L: 8.5-9
76051-112	XL: 9.5-10
76051-114	XXL: 10.5-11



Supercede S6

AQL
1.5

Specifications

Material:	Nitrile
Color:	Light Blue
Inner Surface:	Powder-Free
Outer Surface:	Fully Textured
Length:	9.5"
Finger Thickness:	6.1 mil
Palm Thickness:	5.4 mil
Dispenser Quantity:	100 Gloves
Case Quantity:	1000 Gloves

Mfr. No.	Size
76205-196	S: 6.5-7
76205-198	M: 7.5-8
76205-200	L: 8.5-9
76205-166	XL: 9.5-10
76205-168	XXL: 10.5-11

SW nitrile gloves with EcoTek biodegrade by 92.6% in only 2 years 6 months under test method ASTM D5526-12, while competing products show virtually no degradation in the same time period.

Learn More: www.VWR.COM/SWSAFETY



100% AIR AND WATER TESTED. THAT'S THE SW® DIFFERENCE.



Aloeform X6

AQL
0.65

Specifications

Material:	Nitrile
Color:	Light Green
Inner Surface:	Powder-Free/ACTIValoe
Outer Surface:	Fully Textured
Length:	9.5"
Finger Thickness:	5.9 mil
Palm Thickness:	5.1 mil
Dispenser Quantity:	100 Gloves
Case Quantity:	1000 Gloves

Mfr. No.	Size
76053-592	S: 6.5-7
76053-086	M: 7.5-8
76053-088	L: 8.5-9
76053-090	XL: 9.5-10
76053-092	XXL: 10.5-11



Color-Q pH

AQL
0.65

Specifications

Material:	Nitrile
Color:	Blue / White
Inner Surface:	Powder-Free/pH Natural
Outer Surface:	Fully Textured
Length:	9.5"
Finger Thickness:	5.1 mil
Palm Thickness:	4.3 mil
Dispenser Quantity:	100 Gloves
Case Quantity:	1000 Gloves

Mfr. No.	Size
76050-976	XS: 5.5-6
76050-978	S: 6.5-7
76050-980	M: 7.5-8
76050-982	L: 8.5-9
76050-984	XL: 9.5-10
76050-986	XXL: 10.5-11



Hydrex

AQL
0.65

Specifications

Material:	Nitrile
Color:	Blue
Inner Surface:	Powder-Free / DriTek
Outer Surface:	Fully Textured
Length:	9.5"
Finger Thickness:	6.9 mil
Palm Thickness:	6.1 mil
Dispenser Quantity:	50 Gloves
Case Quantity:	500 Gloves

Cat. No.	Size
76050-998	XS: 5.5-6
76051-000	S: 6.5-7
76051-002	M: 7.5-8
76051-004	L: 8.5-9
76051-006	XL: 9.5-10
76051-008	XXL: 10.5-11



800 932 5000 | VWR.COM

Prices and product details are current when published and subject to change without notice. | Certain products may be limited by federal, state, provincial, or local regulations. | VWR, part of Avantor, makes no claims or warranties concerning sustainable/green products. Any claims concerning sustainable/green products are the sole claims of the manufacturer and not those of VWR International, LLC and/or Avantor, Inc. or affiliates. All prices are in US dollars unless otherwise noted. Offers valid in US and Canada, void where prohibited by law or company policy, while supplies last. | Trademarks are owned by Avantor, Inc. or its affiliates, unless otherwise noted. | Visit vwr.com to view our privacy policy, trademark owners, and additional disclaimers. © 2019 Avantor, Inc. All rights reserved.