

Glassware washer buying guide

Undercounter & freestanding units



START HERE:

Do you need a single unit (6 ft³ volume with 2 racks) or are you interested in two, vertically stacked units (12 ft³ volume with 4 racks)?

Two Vertically Stacked Units:
See Vertical SpaceSaver™ decision tree (reverse side)

Single Unit:
Do you need a glassware washer to fit under the counter of a lab cabinet or do you prefer a freestanding/mobile unit?

Undercounter
Do you need the glassware washer to hold narrow neck glassware such as volumetric or Erlenmeyer flasks?

Freestanding/ Mobile
Do you need the glassware washer to hold narrow neck glassware such as volumetric or Erlenmeyer flasks?

YES
You need a washer that is spindle rack ready. All Spindle Rack Ready units include a DI rinse function. Is your incoming DI water pressure 18 psi or greater? (Spindle racks are sold separately)



Spindle Rack
Cat. No. 10002-950 (Set)

NO
Do you need a DI rinse function?

YES
You need a washer that is spindle rack ready. All Spindle Rack Ready units include a DI rinse function. Is your incoming DI water pressure 18 psi or greater? (Spindle racks are sold separately)



Spindle Rack
Cat. No. 10002-950 (Set)

NO
Do you need a DI rinse function?

YES
Recommended VWR glassware washer: **Cat. No. 10002-858** (Spindle racks are sold separately)

NO
For incoming DI water pressure of less than 18 psi, you will need a DI booster pump. Recommended VWR glassware washer: **Cat. No. 10002-874** (Spindle racks are sold separately)

YES
Is the incoming DI water pressure 18 psi or greater?

YES
Recommended VWR glassware washer: **Cat. No. 10002-854**

NO
Recommended VWR glassware washer: **Cat. No. 10002-850**

NO
For incoming DI water pressure of less than 18 psi, you will need a washer with a DI booster pump. Recommended VWR glassware washer: **Cat. No. 10002-870**

YES
Recommended VWR glassware washer: **Cat. No. 10002-860** (Spindle racks are sold separately)

NO
For incoming DI water pressure of less than 18 psi, you will need a booster pump. Recommended VWR glassware washer: **Cat. No. 10002-876** (Spindle racks are sold separately)

YES
Recommended VWR glassware washer: **Cat. No. 10002-856**

YES
Is the incoming DI water pressure 18 psi or greater?

NO
For incoming DI water pressure of less than 18 psi, you will need a washer with a DI booster pump. Recommended VWR glassware washer: **Cat. No. 10002-872**

NO
Recommended VWR glassware washer: **Cat. No. 10002-852**

Glassware washer buying guide

Vertical SpacesaverTM units



START HERE:

Do you need a single unit (6 ft³ volume with 2 racks) or are you interested in two, vertically stacked units (12 ft³ volume with 4 racks)?

Single Unit:
undercounter & freestanding units
decision tree (reverse side)

Two Vertically Stacked Units:
Vertical SpaceSaverTM units feature two, independent washers unit 1 (top) and unit 2 (bottom) securely stacked on top of each other in a stainless steel housing.

For Unit 2: Do you need this washer to hold narrow neck glassware such as volumetric or Erlenmeyer flasks? (If you select yes, you will later have the choice to include this option for unit 1.)

YES
You need a washer that is Spindle Rack Ready. (Spindle racks are sold separately.)

All Spindle Rack Ready units include a DI rinse function. Is your incoming DI water pressure 18 psi or greater?

NO
Do you need a DI rinse?

YES
Do you need the same Spindle Rack Ready features for unit 2?

NO
If your incoming DI water pressure is less than 18 psi, you will need a booster pump. Booster pumps are purchased separately: **Cat. No. 82021-092**

YES
All bottom units with DI rinse are also Spindle Rack Ready, but may be used with the included, standard racks. Is the incoming DI water pressure 18 psi or greater?

NO
If unit 2 is a standard tap water washer, unit 1 must be the same. Recommended VWR SpaceSaverTM washer stack: **Cat. No. 10002-866**

YES
Recommended VWR SpaceSaverTM washer stack:
Cat. No. 10002-864
(Spindle racks are sold separately)

Spindle Rack
Cat. No. 10002-950 (Set)



NO
Unit 2 only will be Spindle Rack Ready. Recommended VWR SpaceSaverTM washer stack: **Cat. No. 10002-862**
(Spindle racks are sold separately)

Spindle Rack
Cat. No. 10002-950 (Set)



YES
When unit 2 has a DI rinse function, unit 1 must also have a DI rinse function but does not have to be Spindle Rack Ready. Recommended VWR SpaceSaverTM washer stack is: **Cat. No. 10002-862**

NO
Recommended VWR SpaceSaverTM washer stack: **Cat. No. 10002-862**. However you will need to purchase a DI booster pump **Cat. No. 82021-092** separately if your incoming DI water pressure is less than 18 psi.

Vertical Stack Washer pictured is representational only. Spindle and standard racks cannot be used in the same unit.



part of avantor

800 932 5000 | VWR.COM

Prices and product details are current when published and subject to change without notice. | Certain products may be limited by federal, state, provincial, or local regulations. | VWR, part of Avantor, makes no claims or warranties concerning sustainable/green products. Any claims concerning sustainable/green products are the sole claims of the manufacturer and not those of VWR International, LLC and/or Avantor, Inc. or affiliates. All prices are in US dollars unless otherwise noted. Offers valid in US and Canada, void where prohibited by law or company policy, while supplies last. | Trademarks are owned by Avantor, Inc. or its affiliates, unless otherwise noted. | Visit vwr.com to view our privacy policy, trademark owners, and additional disclaimers. © 2019 Avantor, Inc. All rights reserved.