



Item Number S1454
Item Sodium Sulfate, Anhydrous, FCC
CAS Number 7757-82-6
Molecular Formula Na₂SO₄
Molecular Weight 142.04

Assay (Na₂SO₄; Dried Basis)..... 99.0-100.5%

MAXIMUM LIMITS

Lead..... 2 mg/kg
Loss on Drying..... 1%
Selenium..... 0.003%

Spectrum Chemicals & Laboratory Products

Corporate Headquarters: **East Coast Facility:**
14422 S. San Pedro St. 755 Jersey Ave.
Gardena, CA 90248 New Brunswick, NJ 08901
(800) 813-1514 • (310) 516-8000 • Fax (310) 516-9843