

# Certificate of Quality

<b>VWR Catalogue Number</b>	See List
<b>Description</b>	VWR® Cleanroom Sponges
<b>Country of Origin</b>	Manufactured in USA
<b>Date of Issue (yyyy-mm-dd)</b>	2021-09-13

VWR Catalogue Number	Description	Sterility
89495-464	PolyMesh Polyester/Polyester Abrasive Mesh Covered Sponge	Non Sterile
89495-460	PolyTex Sponge	Non Sterile
89495-462	PolyTex Covered Sponge, 100% Polyester Textured Cover, Irradiated	Sterile
89495-456	Microfiber Sponge, 80% Polyester, 20% Polyimide	Non Sterile
89495-458	Microfiber Sponge, 80% Polyester, 20% Polyimide, Irradiated	Sterile

## Quality System Compliance

Products are manufactured under the [ISO 9001:2015](#) standard. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

## QC Testing

Representative products samples are collected and inspected in accordance with current applicable product specifications.

## Product Specifications

<b>Material</b>	100% Textured Polyester Cover, Urethane Foam.
<b>Sterilization</b>	Products labelled as sterile are irradiated using VDmax25 sterilization method with a specified dose range of 25 kGy – 40 kGy.
<b>ATP Assay</b>	Not applicable.
<b>DNase &amp; RNase Free</b>	Not applicable.

<b>BSE/TSE</b>	No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).
<b>Cytotoxicity</b>	Not applicable.
<b>Non-Pyrogenic Statement</b>	Not applicable.
<b>Latex Statement</b>	Not applicable.
<b>BPA Statement</b>	Not applicable.
<b>DEHP Statement</b>	Not Applicable.
<b>RoHS</b>	Not applicable.
<b>REACH Statement</b>	Not applicable.
<b>Storage Conditions</b>	Store away from direct sunlight. Store in space that is protected from extreme temperature and humidity fluctuations and away from Ultraviolet light sources.
<b>Shelf Life</b>	Products labelled sterile has a shelf life of about 5 years whereas shelf life is not applicable for non-sterile products.

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.  
This document has been produced electronically and is valid without a signature.