


TECHNICAL DATA SHEET

Description

HALYARD* **PUREZERO*** MARIN* (9.5") Nitrile Gloves are designed for use in laboratories, research environments and clean areas with applications in pharmaceutical, medical device manufacturing, biotechnology and food processing/handling. HALYARD* **PUREZERO*** MARIN* Nitrile Gloves are tested for use against 57 chemicals, 14 chemotherapy drugs and Fentanyl. The beaded cuff allows for ease of donning, and textured fingertips provide tactile agility.



Physical Properties


AQL	1.0
Non-Sterile	✓
Ambidextrous	✓
Textured Fingertips	✓
Not Made With Natural Rubber Latex	✓
Powder-Free	✓
Tensile Strength ¹	42 MPa (Target)
Ultimate Elongation ¹	600%
Shelf Life	3 Years



Manufactured at
OUR OWN Safeskin
 Facility in Thailand.

SMARTPULL* Dispenser

PUREZERO* MARIN* gloves have a SMARTPULL* Dispenser, which incorporates two separate openings on the box. The first, smaller opening is used when the box is full to control dispensing and lessen waste. When the box is half empty, the second, larger opening allows easier access to the gloves.



Glove Dimensions

	X-SMALL	SMALL	MEDIUM	LARGE	X-LARGE
MARIN* Glove Length (inch/cm)	9.5"/24	9.5"/24	9.5"/24	9.5"/24	9.5"/24
Middle Finger Length (mm)	72	78	82	86	95
Width of Palm (mm)	70	80	95	110	120
Fingertip Thickness	0.15 mm (6 mil)	0.15 mm (6 mil)	0.15 mm (6 mil)	0.15 mm (6 mil)	0.15 mm (6 mil)
Palm Thickness	0.12 mm (4.7 mil)	0.12 mm (4.7 mil)	0.12 mm (4.7 mil)	0.12 mm (4.7 mil)	0.12 mm (4.7 mil)
Cuff Thickness	0.10 mm (4 mil)	0.10 mm (4 mil)	0.10 mm (4 mil)	0.10 mm (4 mil)	0.10 mm (4 mil)



HALYARD* PUREZERO* MARIN* Nitrile Gloves For Industrial/Lab/Food Handling Use

TECHNICAL DATA SHEET

Ordering Information

Product Name	VWR Code	Code	Size	Box Qty.	Case Qty.
HALYARD* PUREZERO* MARIN* Nitrile Gloves - 9.5"/24cm Length	76553-996	LFS1117XS	XS	100	1000
	76553-992	LFS1117SM	S	100	1000
	76554-028	LFS1117MD	M	100	1000
	76554-040	LFS1117LG	L	100	1000
	76554-022	LFS1117XL	XL	90	900

Quality & Regulatory Standards

Compliant to these regulatory standards:

- ISO 9001
- ISO 13485: 2016
- ISO 10993
- ISO 14001: 2015 Environmental Management System

Compliant to these food handling regulatory standards:

- FDA 21 CFR 177-2600
- Commission Regulation (EU) No 10/2011
- EN ISO 374-5:2016 Virus Protection
- EN ISO 374-1:2016/Type C K-Low Chemical Protection
- Compliant with the REACH regulation

US Standards


- ASTM F1671
- ASTM F739-12
- ASTM D6978-05

Brazil PPE

- ISO 21420:2020+A1:2022
- ABNT NBR ISO 374-1:2018
- ISO 374-2:2019
- BS EN 16523-1:2015+A1:2018
- ISO 374-4:2019

EU PPE

- EN ISO 21420:2020
- EN ISO 374-1:2016+A1:2018
- EN ISO 374-2:2019
- EN 16523-1:2015+A1:2018
- EN ISO 374-4:2019
- ISO-16604:2004



For additional information
or samples, contact your
local distributor.

1 Tested per ASTM D412

This fact sheet has been created using the most recent information. In the interest of continuous improvement, the characteristics of the products may change without prior notice.



HALYARD* is an Owens & Minor brand