

MSDS

m-CRESOL PURPLE INDICATORS, <= 0.1% (w/v) Aqueous Solutions

Material Safety Data Sheet

Section 1: Chemical Product and Company Identification

Catalog Number: 2400, 2405, C-215, C-220, C027000, SC027000	
Product Identity: m-CRESOL PURPLE INDICATORS, <= 0.1% (w/v) Aqueous Solutions	
Manufacturer's Name: RICCA CHEMICAL COMPANY LLC	Emergency Contact(24 hr) -- CHEMTREC® Domestic: 800-424-9300 International: 703-527-3887
CAGE Code: 0V553	
Address: 448 West Fork Dr Arlington, TX 76012	Telephone Number For Information: 817-461-5601
Date Prepared: 3/31/99	Revision: 6 Last Revised: 03/14/2007 Date Printed: 06/15/2007 4:53:01 pm

Section 2. Composition/Information on Ingredients

Component	CAS Registry #	Concentration	ACGIH TLV	OSHA PEL
m-Cresol Purple, Sodium Salt	62625-31-4	<= 0.1	Not Available Not Available	Not Available Not Available
Water, Deionized	7732-18-5	Balance	Not Available Not Available	Not Available Not Available

Section 3: Hazard Identification

Emergency Overview: Non-flammable, non-toxic, non-corrosive. Does not present any significant health hazards at this low concentration. Wash areas of contact with water.

Target Organs: eyes, skin.

Eye Contact: May cause slight irritation.

Inhalation: Not likely to be hazardous by inhalation.

Skin Contact: May cause slight irritation and temporary discoloration of skin.

Ingestion: Large doses may cause stomach upset.

Chronic Effects/Carcinogenicity: None

IARC - No.

NTP - No.

OSHA - No.

Reproductive Information: Not Applicable.

Teratology (Birth Defect) Information: Not Applicable.

MSDS

m-CRESOL PURPLE INDICATORS, <= 0.1% (w/v) Aqueous Solutions

Section 4: First Aid Measures - In all cases, seek qualified evaluation.

Eye Contact: Irrigate immediately with large quantity of water for at least 15 minutes. Call a physician if irritation develops.

Inhalation: Not expected to require first aid measures.

Skin Contact: Wash areas of contact with soap and water for at least 15 minutes. Call a physician if irritation develops.

Ingestion: Dilute with water or milk. Call a physician if necessary.

Section 5: Fire Fighting Measures

Flash Point: Not Available.

Method Used: Not Available.

LFL: Not Available.

UFL: Not Available.

Extinguishing Media: Use any means suitable for extinguishing surrounding fire.

Fire & Explosion Hazards: Not considered to be a fire or explosion hazard.

Fire Fighting Instructions: Use normal procedures/instructions.

Fire Fighting Equipment: Use protective clothing and breathing equipment appropriate for the surrounding fire.

Section 6: Accidental Release Measures

Absorb with suitable material and treat as normal refuse. Small amounts of the liquid may be flushed to the drain with excess water. Always dispose of in accordance with local regulations.

Section 7. Handling and Storage

As with all chemicals, wash hands thoroughly after handling. Avoid contact with eyes and skin. Protect from freezing and physical damage.

Safety Storage Code: General

Section 8: Exposure Control/Personal Protection

Engineering Controls: No specific controls are needed. Normal room ventilation is adequate.

Respiratory Protection: Normal room ventilation is adequate.

Skin Protection: Chemical resistant gloves.

Eye Protection: Safety glasses or goggles.

Section 9: Physical and Chemical Properties

Appearance: Yellow to Purple liquid

Odor: Odorless

Solubility in Water: Infinite

Specific Gravity: Approximately 1

pH: 7 - 9

Boiling Point(°C): Approximately 100

Melting Point(°C): Approximately 0

Vapor Pressure: Not Applicable.

Section 10: Stability and Reactivity

Chemical Stability: Stable under normal conditions of use and storage.

Incompatibility: Strong oxidizing agents

Hazardous Decomposition Products: Oxides of Sulfur, Carbon Monoxide and Carbon Dioxide may form when heated to decomposition.

Hazardous Polymerization: Will not occur.

Section 11. Toxicological Information

No LD50 information found.

Section 12. Ecological Information

MSDS

m-CRESOL PURPLE INDICATORS, <= 0.1% (w/v) Aqueous Solutions

Ecotoxicological Information: No information found.

Chemical Fate Information: No information found.

Section 13. Disposal Considerations

Dilute with water and flush to the sewer if local regulations allow. If not allowed, containerize for proper disposal at an approved waste disposal facility. Always dispose of in accordance with local, state and federal regulations.

Section 14. Transport Information

Part Numbers:

This product is not regulated.

Section 15. Regulatory Information (Not meant to be all inclusive - selected regulation represented)

OSHA Status: The above items either do not contain any specifically hazardous material or the potentially hazardous material is present in such low concentration that the items do not present any immediate threat to health and safety. These items do not meet the OSHA Hazard Communication Standard (29 CFR 1910.1200) definition of a hazardous material.

TSCA Status: All components of this solution are listed on the TSCA Inventory or are mixtures (hydrates) of items listed on the TSCA Inventory.

Sara Title III:

Section 302 Extremely Hazardous Substances: Not Applicable.

Section 311/312 Hazardous Categories: No

Section 313 Toxic Chemicals: Not Applicable.

California: None Reported.

Pennsylvania: None Reported.

RCRA Status: Not Applicable.

CERCLA Reportable Quantity: None Reported.

WHMIS: Not Applicable.

Section 16. Other Information

NFPA Ratings:

Health: 0

Flammability: 0

Reactivity: 0

Special Notice Key: None

HMIS Ratings:

Health: 0

Flammability: 0

Reactivity: 0

Protective Equipment: B (Protective Eyewear, Gloves)



MSDS

For RICCA, SpectroPure, Red Bird, and Solutions Plus Brands
Emergency Contact (24 hr) -- CHEMTREC®
Domestic: 800-424-9300
International: 703-527-3887

m-CRESOL PURPLE INDICATORS, <= 0.1% (w/v) Aqueous Solutions

Rev 1, 9-20-2000: Reformatted from WordPerfect® to Microsoft Word®; (Section 1) Revised emergency telephone number to CHEMTREC® 800-424-9300; (Section 3) added skin to target organs; (Section 7) added storage code.
Rev 2, 10-09-2001: Reformatted to electronic data format.
Rev 3, 10-31-2002: (Section 1) added Solutions Plus catalog number C027000.
Rev 4, 02-02-2006: (Section 1) added Red Bird catalog number C-215.
Rev 5, 07-11-2006: (Section 9) revised color from burgundy, revised pH from approximately 7.
Rev 6, 03-14-2007: (Section 1) added Red Bird catalog number C-220.

When handled properly by qualified personnel, the product described herein does not present a significant health or safety hazard. Alteration of its characteristics by concentration, evaporation, addition of other substances, or other means may present hazards not specifically addressed herein and which must be evaluated by the user. The information furnished herein is believed to be accurate and represents the best data currently available to us. No warranty, expressed or implied, is made and RICCA CHEMICAL COMPANY assumes no legal responsibility or liability whatsoever resulting from its use.