



Material	23510-540
Material description	Chloroform
Grade	ACS
CAS Number	67-66-3
Molecular formula	CHCl <sub>3</sub>
Molecular mass	119.38

## Characteristics

## Specifications

Assay (additives included)	Min. 99.8 %
Acetone and aldehyde	Passes test
Acid and chloride	Passes test
Free chlorine	Passes test
Presence of preservative (Ethanol)	Passes test
Substances darkened by sulphuric acid	Passes test
Color (Pt-Co)	Max. 10 APHA
Residue after evaporation	Max. 0.001 %
Pb (Lead)	Max. 0.05 ppm

## Signature

We certify that this batch conforms to the specifications listed above.

This document has been produced electronically and is valid without a signature.

Michel Boisvert, Head of the Laboratory - Lachine  
VWR Anachemia Canada, 255 rue Norman, Lachine, H8R 1A3,  
QC, Canada

For Professional use in Laboratory or Manufacturing. Not for use as an Active Pharmaceutical Ingredient or Food or Animal Feed. Suitability and intended use of the product remains the responsibility of the user.