



# CERTIFICATE OF CONFORMANCE

Product Part Number: ITCAA-9IR  
Lot Number: 010380-000  
Irradiation Reference Number: 0002085-00  
Lot Quantity: 60 CASE(S)  
Date of manufacture: 01/20/2019  
Shelf Life: N/A  
Description: ISOLATOR TOOL COVER, MICROFIBER, IRR'D (50/CS)  
Size -  
Width: 2.5"  
Length: 9"  
Height: N/A

*This certification is provided as full assurance that the following product code and lot number was manufactured in accordance with prescribed procedures and specifications.*

01/23/2019

Authorized Signature  
Quality Control Department

Date

CC: With Product  
Customer File  
C of C File

3431 WEST LOMITA BOULEVARD • TORRANCE, CA 90505-5010 Tel. (310) 784-6990 • FAX (310) 784-6980

STERIS (S) 01/19/2019 09:39:31 AM

STERIS (S) 01/21/2019 10:15 AM

Prepared for MICRONOVA MED INC

Gamma Process Run ID 113218A

Product Code	Product Lot Number	Quantity	UOM
HCAA 9IR	010380 000	60	CS

Other Information:

0002085 00

Processing Run Start Date/Time: 20-Jan-2019 09:39:31 am      Approx. Downtime (hours): 0.00  
 Processing Run End Date/Time: 20-Jan-2019 02:33:29 pm

Minimum Specified Dose (kGy):	40.0	Minimum Delivered Dose (kGy):	42.7
Maximum Specified Dose (kGy):	60.0	Maximum Delivered Dose (kGy):	55.8
Product meets Customer specifications; zero nonconformities occurred during this irradiation run.			

Signature Manifest

Reviewed and E-Signed By  
*Rony Mejia (QS/RC Analyst)*  
 Document Content Revision: 1

Signed On 1/21/2019 at 10:15 AM  
 UTC / GMT Offset (hh:mm) -8:00

Processing Location:  
 STERIS  
 1000 South Sarah Place  
 Ontario, CA 91761  
 Phone: 909-390-9942  
 Fax: 909-390-7854

Operating facilities are in compliance with applicable state and federal regulations (FDA, NRC, EPA, and OSHA) and provide services under a quality system which meets the requirements of FDA QSR, EN/ISO 13485, and in alignment with EN ANSIA/MIL/ISO 11137. STERIS certifies that these processed items received the indicated doses within the precision and accuracy of the dosimetry system used.