

Name of the Product : **M-Endo Broth MF (MF Endo Medium) (M-Coliform Broth)**
 Code No. : **M1103**
 Section 1 : **Chemical Identification**
 Code No. : M1103
 Name of the Product : M-Endo Broth MF (MF Endo Medium) (M-Coliform Broth)
 Produced by : HiMedia Laboratories Pvt. Ltd.
 Address : 23, Vadhani Indl. Estate, LBS Marg, Mumbai 400 086, India.
 Tel. No. : 2500 0970, 2500 1607 Fax No. 022 2500 2468

Section 2 : **Composition**

Ingredients	Grams/Litre	Ingredients	Grams/Litre
Tryptose	10.00	Casein enzymic hydrolysate	5.00
Peptone, special	5.00	Yeast extract	1.50
Lactose	12.50	*Sodium deoxycholate	0.10
Dipotassium phosphate	4.375	Monopotassium phosphate	1.375
Sodium chloride	5.00	*Sodium lauryl sulphate	0.05
*Sodium sulphite	2.10	*Basic fuchsin	1.05
Final pH (at 25°C) 7.2±0.2			

* Harmful material



Sodium deoxycholate
 CAS No. : 302-95-4
 R : 22-37 S : 22-26-36
 RTECS : FZ 2250000



Sodium lauryl sulphate
 CAS No. : 151-21-3
 R : 22-36/38 S : 26-36
 RTECS : WT 1050000



Sodium sulphite
 CAS No. : 7757-83-7
 R : 22-31-36/37/38 S : 26-36/37
 RTECS : WE 2150000



Basic fuchsin
 CAS No. : 632-99-5
 R : 40 S : 36/37
 RTECS : CX 9850000

Section 3 : **Hazards Identification**



Hazard : Inhalation and ingestion of, or skin penetration by these substances is harmful to one's health. Non-recurring, recurring or lengthy exposure to these substances may result in irreversible damage.
 Caution : Avoid contact with the human body, including inhalation of the vapours and in cases of malaise consult a doctor.

Section 4 : **First - Aid Measures**

In case of swallowing, wash out mouth with water provided person is conscious. In case of skin contact, flush with copious amounts of water for at least 15 minutes. Remove contaminated clothing and shoes. In case of inhalation, remove to fresh air. If breathing becomes difficult, call physician. In case of contact with eyes, wash with copious amounts of water for at least 15 minutes. Assure adequate washing by separating the eyelids with fingers. In serious conditions, call a physician, show the container or label.

Section 5 : **Fire Fighting Measures**

Extinguishing Media. Water spray.

Carbon dioxide, dry chemical powder or appropriate foam. Special fire-fighting procedures. Wear self-contained breathing apparatus and protective clothing to prevent contact with skin and eyes. Unusual fire and explosions hazards Emits toxic fumes under fire conditions.

Section 6 : **Accidental Release Measures**

Wear respirator, chemical safety goggles, rubber boots and heavy rubber gloves. Sweep up, place in a bag and hold for waste disposal. Avoid raising dust. Ventilate area and wash spill site after material pickup is complete.

Section 7 : **Handling and Storage**

Handling - Refer to Section 8

Storage - Store below 30°C

Section 8 : **Exposure Controls/Personal Protection**

Wear appropriate NIOSH/MSHA-approved respirator, chemical-resistant gloves, safety goggles, other protective clothing. Mechanical exhaust required.

Section 9 : **Physical and Chemical Properties**

Appearance : Homogeneous powder

Colour : Light pink to purple

Section 10 : **Stability and Reactivity**

Stability : Product is stable if stored as per the conditions specified under storage of Section No. 7. Product loses its potency/performance above 45°C.

Conditions to avoid : Heat and light.

Decomposition products : if product is decomposed toxic fumes of carbon monoxide and carbon dioxide may evolve. Hazardous polymerization will not occur.

Section 11 : **Toxicological Information**

Acute Effects : May be harmful if swallowed.

Exposure can cause : Stomach pains, vomiting, diarrhoea. Prolonged or repeated exposure may cause allergic reactions in certain sensitive individuals.

Target organ(s) : G.I. System

Exposure can cause : Stomach pains, vomiting, diarrhoea. Prolonged or repeated exposure may cause allergic reactions in certain sensitive individuals.

RTECS No. : FZ 2250000 WT 1050000 WE 2150000 CX 9850000

For complete information RTECS may be referred.

Section 12 : **Ecological Information**

Data not available

Section 13 : **Disposal Considerations**

On completion of work all used or unusable preparations of this product and derivatives thereof are to be disposed off by autoclaving and/or by incineration. Dispose off waste in accordance with all applicable Federal, State and local laws.

Section 14 : **Transport Information**

UN No. : Not applicable.

Section 15 : **Regulatory Information**

Risk Phrases : Harmful if swallowed

Risk Phrases : Irritating to respiratory system

Risk Phrases : Irritating to eyes and skin

Risk Phrases : Contact with acids liberates toxic gas

Risk Phrases : Irritating to eyes, respiratory system and skin

Risk Phrases : Limited evidence of a carcinogenic effect

Safety Phrases : Do not breathe dust

Safety Phrases : In case of contact with eyes, rinse immediately with plenty of water and seek medical advice

Safety Phrases : Wear suitable protective clothing and gloves



Significance of signs : Inhalation and ingestion of, or skin penetration by these substances is harmful to one's health. Non-recurring, recurring or lengthy exposure to these substances may result in irreversible damage.

Section 16 : **Other Information**

The information contained in this data sheet represents the best information currently available to us. However, no warranty is made with respect to its completeness and we assume no liability resulting from its use. The information is offered solely for user's obligation to investigate and determine the suitability of the information for their particular purpose.

