Material Safety Data Sheet



Issuing Date 12/22/2010 Revision Number 0

1. PRODUCT AND COMPANY IDENTIFICATION

Product Name Thymol Blue Indicator

Product Code(s) 2213

Synonyms none

Recommended Use Laboratory chemicals. Industrial (not for food or food contact use). Test kit reagent.

Company LaMotte Company, Inc.

802 Washington Avenue

P.O. Box 329

Chestertown, MD 21620

USA

Emergency Telephone Number 24 Hour Emergency Number (CHEM-TEL):

USA, Canada, Puerto Rico 1-800-255-3924

Outside North American Continent (Call collect) 813-248-0585

2. HAZARDS IDENTIFICATION

Emergency Overview

May stain skin

The product contains no substances which at their given concentration are considered to be hazardous to health

Appearance Varies, Amber to blue Physical State Liquid Odor None

OSHA Regulatory Status This material is not considered hazardous by the OSHA Hazard Communication Standard

(29 CFR 1910.1200).

Potential Health Effects

Principle Routes of Exposure Skin contact

Acute Toxicity

Eves May cause slight irritation.

SkinSubstance may cause slight skin irritation.InhalationNot an expected route of exposure.IngestionNot an expected route of exposure.

Chronic Effects

3. COMPOSITION/INFORMATION ON INGREDIENTS

Chemical Name	CAS-No	Weight %
Thymol blue	76-61-9	<0.1
Water	7732-18-5	to 100%

Sodium hydroxide and/or hydrochloric acid added for pH adjustment.

4. FIRST AID MEASURES

Thymol Blue Indicator Product Code(s) 2213

General Advice Do not get in eyes, on skin, or on clothing.

Eye Contact Rinse immediately with plenty of water, also under the eyelids, for at least 15 minutes. If

irritation persists or develops, contact a physician.

Skin Contact Wash skin with soap and water. If irritation develops or persists, consult physician.

Inhalation Not expected to require first aid measures.

Ingestion Drink 1 or 2 glasses of water. If gastrointestinal distress occurs contact physician.

Notes to Physician Treat symptomatically.

Protection of First-aiders

Use personal protective equipment. Do not use mouth-to-mouth method if victim ingested or

inhaled the substance; induce artificial respiration with the aid of a pocket mask equipped

with a one-way valve or other proper respiratory medical device.

5. FIRE-FIGHTING MEASURES

Flammable Properties Not a fire hazard.

Flash Point Not applicable

Suitable Extinguishing Media Water spray, dry chemical, carbon dioxide (CO₂), or foam.

Explosion Data

NFPA Health Hazard 0 Flammability 0 Stability 0 Physical and Chemical Hazards -

6. ACCIDENTAL RELEASE MEASURES

Personal Precautions Refer to Section 8.

Methods for Containment Absorb spill with inert material (e.g. dry sand or earth), then place in a chemical waste

container.

Methods for Cleaning Up Soak up with inert absorbent material. After cleaning, flush away traces with water.

7. HANDLING AND STORAGE

Handling Handle in accordance with good industrial hygiene and safety practice. Prevent contact with

skin, eyes and clothing. Do not ingest. Do not eat, drink or smoke when using this product.

Storage Keep tightly closed in a dry and cool place. Keep out of the reach of children.

8. EXPOSURE CONTROLS / PERSONAL PROTECTION

Chemical Name	ACGIH TLV	OSHA PEL	NIOSH IDLH
Thymol blue	None Known	None Known	None Known
76-61-9			
Water	None Known	None Known	None Known
7732-18-5			

Personal Protective Equipment

Eye/Face Protection Safety glasses with side-shields.

Skin and Body Protection Protective gloves.

Respiratory Protection None required under normal usage.

Hygiene Measures Do not eat, drink or smoke when using this product.

Thymol Blue Indicator Product Code(s) 2213

9. PHYSICAL AND CHEMICAL PROPERTIES

AppearanceVaries Amber to blueOdorNonePhysical StateLiquidpH7.5 - 9.0Flash PointNot applicableBoiling Point/Rangeca. 100 °C

Vapor Pressure <17 mmHg @ 20°C Vapor Density <1(Air = 1)

10. STABILITY AND REACTIVITY

Stability Stable under normal conditions of use and storage.

Incompatible Products None known.

Conditions to Avoid Extremes of temperature and direct sunlight.

Hazardous Decomposition Products None under normal use.

Hazardous Polymerization Hazardous polymerization does not occur.

11. TOXICOLOGICAL INFORMATION

Acute Toxicity

Chemical Name	LD50 Oral	LD50 Dermal	LC50 Inhalation	
Thymol blue	None Known	None Known	None Known	
Water	90 mL/kg (Rat)	None Known	None Known	

Chronic Toxicity

Carcinogenicity Contains no ingredient listed as a carcinogen.

Chemical Name	ACGIH	IARC	NTP	OSHA
Thymol blue	None Known	None Known	None Known	None Known
Water	None Known	None Known	None Known	None Known

Chemical Name	EU - Endocrine Disrupters Candidate List	EU - Endocrine Disruptors - Evaluated Substances	Japan - Endocrine Disruptor Information
Thymol blue	None Known	None Known	None Known
Water	None Known	None Known	None Known

12. ECOLOGICAL INFORMATION

Ecotoxicity

No information available.

Chemical Name	Toxicity to Algae	Toxicity to Fish	Microtox	Daphnia Magna (Water Flea)
Thymol blue	None Known	None Known	None Known	None Known
Water	None Known	None Known	None Known	None Known

Persistence and Degradability No information available.

Bioaccumulation/AccumulationNo information available.

Chemical Name	Log Pow
Thymol blue	None Known
Water	None Known

13. DISPOSAL CONSIDERATIONS

Waste Disposal Method

Dispose of in accordance with local regulations. This material, as supplied, is not a hazardous waste according to Federal regulations (40 CFR 261). This material could become a hazardous waste if it is mixed with or otherwise comes in contact with a hazardous waste, if chemical additions are made to this material, or if the material is processed or otherwise altered. Consult 40 CFR 261 to determine whether the altered material is a hazardous waste. Consult the appropriate state, regional, or local regulations for additional requirements.

Chemical Name
Thymol blue - 76-61-9
Water - 7732-18-5

Chemical Name	RCRA - Halogenated Organic Compounds	RCRA - P Series Wastes	RCRA - F Series Wastes	RCRA - K Series Wastes
Thymol blue - 76-61-9	None Known	None Known	None Known	None Known
Water - 7732-18-5	None Known	None Known	None Known	None Known

14. TRANSPORT INFORMATION

DOT Not regulated

IATA Not regulated

IMDG/IMO Not regulated

15. REGULATORY INFORMATION

International Inventories

Component	TSCA	DSL	EINECS/ELIN CS	ENCS	IECSC	KECL	PICCS	AICS
Thymol blue 76-61-9 (<0.1)	Present	Х	X	ENCS	X	KECL	Х	X
Water 7732-18-5 (to 100%)	Present	Х	X	ENCS	X	KE-35400	Х	Х

U.S. Federal Regulations

SARA 313

Section 313 of Title III of the Superfund Amendments and Reauthorization Act of 1986 (SARA). This product does not contain any chemicals which are subject to the reporting requirements of the Act and Title 40 of the Code of Federal Regulations, Part 372.

Chemical Name	CAS-No	Weight %	SARA 313 - Threshold	
			Values %	
Thymol blue	76-61-9	<0.1	None Known	
Water	7732-18-5	to 100%	None Known	

SARA 311/312 Hazard Categories Acute Health Hazard

No

Published Date: 22-Dec-2010 Page 4/6

Chronic Health Hazard No
Fire Hazard No
Sudden Release of Pressure Hazard No
Reactive Hazard No

Clean Water Act

This product does not contain any substances regulated as pollutants pursuant to the Clean Water Act (40 CFR 122.21 and 40 CFR 122.42).

Component	CWA - Reportable	CWA - Toxic Pollutants	CWA - Priority Pollutants	
	Quantities			Substances
Thymol blue	None Known	None Known	None Known	None Known
76-61-9 (<0.1)				
Water	None Known	None Known	None Known	None Known
7732-18-5 (to 100%)				

Clean Air Act, Section 112 Hazardous Air Pollutants (HAPs) (see 40 CFR 61)

This product does not contain any substances regulated as hazardous air pollutants (HAPS) under Section 112 of the Clean Air Act Amendments of 1990.

Chemical Name	CAS-No	Weight %	HAPS data	VOC Chemicals	Class 1 Ozone Depletors	Class 2 Ozone Depletors
Thymol blue	76-61-9	<0.1	None Known	None Known	None Known	None Known
Water	7732-18-5	to 100%	None Known	None Known	None Known	None Known

CERCLA

Chemical Name	Hazardous Substances RQs	Extremely Hazardous Substances RQs
Thymol blue	None Known	None Known
Water	None Known	None Known

U.S. State Regulations

California Proposition 65

This product does not contain any Proposition 65 chemicals

Chemical Name	CAS-No	California Prop. 65
Thymol blue	76-61-9	None Known
Water	7732-18-5	None Known

Chemical Name	Massachusetts	New Jersey	Pennsylvania	Illinois	Rhode Island
Thymol blue	None Known	None Known	None Known	None Known	None Known
Water	None Known	None Known	None Known	None Known	None Known

International Regulations

Mexico - Grade

Chemical Name	Carcinogen Status	Exposure Limits
Thymol blue	None Known	None Known
Water	None Known	None Known

Canada


WHMIS Hazard Class

Non-controlled

Published Date: 22-Dec-2010 Page 5/6

Thymol Blue Indicator Product Code(s) 2213

16. OTHER INFORMATION



Prepared By Regulatory Affairs Department

Issuing Date 12/22/2010

Revision Date

Revision Note Initial Release.

Disclaimer

The information provided on this MSDS is correct to the best of our knowledge, information and belief at the date of its publication. The information given is designed only as a guide for safe handling, use, processing, storage, transportation, disposal and release and is not to be considered as a warranty or quality specification. The information relates only to the specific material designated and may not be valid for such material used in combination with any other material or in any process, unless specified in the text.

End of MSDS