



Item Number P1456
Item Propylene Glycol, FCC
CAS Number 57-55-6
Molecular Formula $C_3H_8O_2$
Molecular Weight 76.10

Assay ($C_3H_8O_2$; by Wt.)..... Min. 99.5%
Acidity..... To pass test
Distillation Range..... 185°-189°C
Specific Gravity @ 25°C..... 1.035-1.037

MAXIMUM LIMITS

Lead..... 1 mg/kg
Residue on Ignition..... 0.007%
Water..... 0.2%

Spectrum Chemicals & Laboratory Products

Corporate Headquarters: East Coast Facility:
14422 S. San Pedro St. 755 Jersey Ave.
Gardena, CA 90248 New Brunswick, NJ 08901
(800) 813-1514 • (310) 516-8000 • Fax (310) 516-9843