

Material Safety Data Sheet

HARLECO® Gill Modified EA



Section 1. Product and Company Identification

Product name : HARLECO® Gill Modified EA
Product code : 65068
Synonym : Dye solution
Material uses : Other non-specified industry: Laboratory Reagent
Manufacturer : EMD Chemicals Inc.
P.O. Box 70
480 Democrat Road
Gibbstown, NJ 08027
856-423-6300 Technical Service
Monday - Friday: 8:00 - 5:00 PM
Validation date : **9/19/2007.**
Print date : 9/19/2007.
In case of emergency : 800-424-9300 CHEMTREC (USA)
613-996-6666 CANUTEC (Canada)
24 Hours/Day: 7 Days/Week

Section 2. Hazards Identification

Physical state : Liquid.
Odor : Alcohol-like.
OSHA/HCS status : This material is considered hazardous by the OSHA Hazard Communication Standard (29 CFR 1910.1200).
Emergency overview : DANGER!
POISON!
FLAMMABLE LIQUID AND VAPOR.
VAPOR MAY CAUSE FLASH FIRE.
VAPOR HARMFUL.
MAY BE FATAL OR CAUSE BLINDNESS IF SWALLOWED.
CANNOT BE MADE NONPOISONOUS.
HARMFUL IF INHALED OR ABSORBED THROUGH SKIN.
CAUSES EYE AND SKIN IRRITATION.
BIRTH DEFECT HAZARD.
CONTAINS MATERIAL WHICH CAN CAUSE BIRTH DEFECT.
CONTAINS MATERIAL WHICH MAY CAUSE DAMAGE TO THE FOLLOWING
ORGANS: BLOOD, REPRODUCTIVE SYSTEM, LIVER, GASTROINTESTINAL TRACT,
RESPIRATORY TRACT, SKIN, CENTRAL NERVOUS SYSTEM, EYE, LENS OR
CORNEA, TEETH
Do not ingest. Avoid contact with eyes, skin and clothing. Avoid breathing vapor or mist.
Keep away from heat, sparks and flame. Keep container closed. Use only with
adequate ventilation. Wash thoroughly after handling. Avoid exposure during
pregnancy.
Routes of entry : Dermal contact. Eye contact. Inhalation. Ingestion.
Potential acute health effects
Eyes : Irritating to eyes.
Skin : Toxic in contact with skin. Irritating to skin.
Inhalation : Toxic by inhalation.
Ingestion : Toxic if swallowed.
Carcinogenic effects : No known significant effects or critical hazards.
Mutagenic effects : No known significant effects or critical hazards.
**Teratogenicity /
Reproductive toxicity** : Contains material which can cause birth defects.

Continued on Next Page

Section 2. Hazards Identification

Medical conditions aggravated by over-exposure : Repeated skin exposure can produce local skin destruction or dermatitis. Repeated or prolonged contact with spray or mist may produce chronic eye irritation and severe skin irritation. Repeated or prolonged exposure to the substance can produce target organs damage.

See toxicological information (section 11)

Section 3. Composition/Information on Ingredients

United States

<u>Name</u>	<u>CAS number</u>	<u>% by Weight</u>
Ethanol	64-17-5	>67
Methanol	67-56-1	25
Light Green SF, Yellowish	5141-20-8	0.1 - 0.9
Eosin Y	17372-87-1	0.1 - 0.9
Acetic Acid	64-19-7	2
Phosphotungstic Acid	12501-23-4	0.1 - 0.9
Water	7732-18-5	3

Section 4. First Aid Measures

- Eye contact** : Get medical attention immediately. Immediately flush eyes with plenty of water for at least 15 minutes, occasionally lifting the upper and lower eyelids. Check for and remove any contact lenses.
- Skin contact** : Get medical attention immediately. Flush contaminated skin with plenty of water. Continue to rinse for at least 10 minutes. Remove contaminated clothing and shoes. Wash contaminated clothing thoroughly with water before removing or wear gloves. Wash clothing before reuse. Clean shoes thoroughly before reuse.
- Inhalation** : Get medical attention immediately. Move exposed person to fresh air. If it is suspected that fumes are still present, the rescuer should wear an appropriate mask or self-contained breathing apparatus. Keep person warm and at rest. If not breathing, if breathing is irregular or if respiratory arrest occurs, provide artificial respiration or oxygen by trained personnel. It may be dangerous to the person providing aid to give mouth-to-mouth resuscitation. If unconscious, place in recovery position and get medical attention immediately. Maintain an open airway. Loosen tight clothing such as a collar, tie, belt or waistband.
- Ingestion** : Get medical attention immediately. Wash out mouth with water. Remove dentures if any. Move exposed person to fresh air. Keep person warm and at rest. If material has been swallowed and the exposed person is conscious, give small quantities of water to drink. Stop if the exposed person feels sick as vomiting may be dangerous. Do not induce vomiting unless directed to do so by medical personnel. If vomiting occurs, the head should be kept low so that vomit does not enter the lungs. Never give anything by mouth to an unconscious person. If unconscious, place in recovery position and get medical attention immediately. Maintain an open airway. Loosen tight clothing such as a collar, tie, belt or waistband.
- Protection of first-aiders** : No action shall be taken involving any personal risk or without suitable training. If it is suspected that fumes are still present, the rescuer should wear an appropriate mask or self-contained breathing apparatus. It may be dangerous to the person providing aid to give mouth-to-mouth resuscitation.

Section 5. Fire Fighting Measures

Flammability of the product : Flammable liquid and vapor. Vapor may cause flash fire. Vapors may accumulate in low or confined areas or travel a considerable distance to a source of ignition and flash back. Runoff to sewer may create fire or explosion hazard.

Products of combustion : These products are carbon oxides (CO, CO₂).

Extinguishing media

Suitable : Use dry chemical, CO₂, water spray (fog) or foam.

Not suitable : Do not use water jet.

Section 5. Fire Fighting Measures

- Special exposure hazards** : Not available.
- Special protective equipment for fire-fighters** : Fire-fighters should wear appropriate protective equipment and self-contained breathing apparatus (SCBA) with a full face-piece operated in positive pressure mode.
- Special remarks on fire hazards** : Dangerous fire and explosion risk. Container explosion may occur under fire conditions or when heated. Vapor may travel a considerable distance to source of ignition and flash back. (Methanol)
- Special remarks on explosion hazards** : Vapor may travel a considerable distance to source of ignition and flash back.

Section 6. Accidental Release Measures

- Personal precautions** : Immediately contact emergency personnel. Eliminate all ignition sources. Keep unnecessary personnel away. Use suitable protective equipment. Do not touch or walk through spilled material.
- Environmental precautions** : Avoid dispersal of spilled material and runoff and contact with soil, waterways, drains and sewers.
- Methods for cleaning up** : If emergency personnel are unavailable, contain spilled material. For small spills, add absorbent (soil may be used in the absence of other suitable materials) and use a non-sparking or explosion-proof means to transfer material to a sealable, appropriate container for disposal. For large spills, dike spilled material or otherwise contain material to ensure runoff does not reach a waterway. Place spilled material in an appropriate container for disposal.

Section 7. Handling and Storage

- Handling** : Do not ingest. Avoid contact with eyes, skin and clothing. Keep container closed. Use only with adequate ventilation. Avoid breathing vapor or mist. Keep away from heat, sparks and flame. To avoid fire or explosion, dissipate static electricity during transfer by grounding and bonding containers and equipment before transferring material. Use explosion-proof electrical (ventilating, lighting and material handling) equipment. Wash thoroughly after handling.
- Storage** : Store in a segregated and approved area. Keep container in a cool, well-ventilated area. Keep container tightly closed and sealed until ready for use. Avoid all possible sources of ignition (spark or flame).

Section 8. Exposure Controls/Personal Protection

Product name

Exposure limits

United States

Ethanol

ACGIH TLV (United States, 1/2005). Notes: 1996 Adoption Refers to Appendix A -- Carcinogens.

TWA: 1880 mg/m³ 8 hour/hours. Form: All forms

TWA: 1000 ppm 8 hour/hours. Form: All forms

NIOSH REL (United States, 12/2001).

TWA: 1900 mg/m³ 10 hour/hours. Form: All forms

TWA: 1000 ppm 10 hour/hours. Form: All forms

OSHA PEL (United States, 8/1997).

TWA: 1900 mg/m³ 8 hour/hours. Form: All forms

TWA: 1000 ppm 8 hour/hours. Form: All forms

OSHA PEL 1989 (United States, 3/1989).

TWA: 1900 mg/m³ 8 hour/hours. Form: All forms

TWA: 1000 ppm 8 hour/hours. Form: All forms

Methanol

ACGIH (United States, 1994). Skin

TWA: 262 mg/m³

STEL: 328 mg/m³

OSHA (United States, 1989). Skin

TWA: 260 mg/m³

STEL: 325 mg/m³

NIOSH REL (United States, 12/2001). Skin

Section 8. Exposure Controls/Personal Protection

STEL: 325 mg/m³ 15 minute/minutes. Form: All forms

STEL: 250 ppm 15 minute/minutes. Form: All forms

TWA: 260 mg/m³ 10 hour/hours. Form: All forms

TWA: 200 ppm 10 hour/hours. Form: All forms

OSHA PEL (United States, 8/1997).

TWA: 260 mg/m³ 8 hour/hours. Form: All forms

TWA: 200 ppm 8 hour/hours. Form: All forms

OSHA PEL 1989 (United States, 3/1989). Skin

STEL: 325 mg/m³ 15 minute/minutes. Form: All forms

STEL: 250 ppm 15 minute/minutes. Form: All forms

TWA: 260 mg/m³ 8 hour/hours. Form: All forms

TWA: 200 ppm 8 hour/hours. Form: All forms

ACGIH TLV (United States, 1/2005). Skin Notes: Substances for which there is a Biological Exposure Index or Indices

STEL: 328 mg/m³ 15 minute/minutes. Form: All forms

STEL: 250 ppm 15 minute/minutes. Form: All forms

TWA: 262 mg/m³ 8 hour/hours. Form: All forms

TWA: 200 ppm 8 hour/hours. Form: All forms

ACGIH TLV (United States, 1/2006).

STEL: 37 mg/m³ 15 minute/minutes. Form: All forms

STEL: 15 ppm 15 minute/minutes. Form: All forms

TWA: 25 mg/m³ 8 hour/hours. Form: All forms

TWA: 10 ppm 8 hour/hours. Form: All forms

NIOSH REL (United States, 12/2001).

STEL: 37 mg/m³ 15 minute/minutes. Form: All forms

STEL: 15 ppm 15 minute/minutes. Form: All forms

TWA: 25 mg/m³ 10 hour/hours. Form: All forms

TWA: 10 ppm 10 hour/hours. Form: All forms

OSHA PEL (United States, 8/1997).

TWA: 25 mg/m³ 8 hour/hours. Form: All forms

TWA: 10 ppm 8 hour/hours. Form: All forms

OSHA PEL 1989 (United States, 3/1989).

TWA: 25 mg/m³ 8 hour/hours. Form: All forms

TWA: 10 ppm 8 hour/hours. Form: All forms

Acetic Acid

Consult local authorities for acceptable exposure limits.

Engineering measures : Use only with adequate ventilation. If user operations generate dust, fumes, vapor or mist, use process enclosures, local exhaust ventilation or other engineering controls to keep worker exposure to airborne contaminants below any recommended or statutory limits. The engineering controls also need to keep gas, vapor or dust concentrations below any lower explosive limits. Use explosion-proof ventilation equipment.

Personal protection

- Eyes** : Safety eyewear complying with an approved standard should be used when a risk assessment indicates this is necessary to avoid exposure to liquid splashes, mists, gases or dusts.
Recommended: face shield
- Skin** : Personal protective equipment for the body should be selected based on the task being performed and the risks involved and should be approved by a specialist before handling this product.
Body: Recommended: lab coat and gloves
- Respiratory** : Use a properly fitted, air-purifying or air-fed respirator complying with an approved standard if a risk assessment indicates this is necessary. Respirator selection must be based on known or anticipated exposure levels, the hazards of the product and the safe working limits of the selected respirator.
Recommended: Vapor respirator. Be sure to use an approved/certified respirator or equivalent. Vapor respirator or self-contained breathing apparatus (SCBA). Wear appropriate respirator when ventilation is inadequate.
- Hands** : Chemical-resistant, impervious gloves complying with an approved standard should be worn at all times when handling chemical products if a risk assessment indicates this is necessary.

Section 8. Exposure Controls/Personal Protection

Hygiene measures : Wash hands, forearms and face thoroughly after handling chemical products, before eating, smoking and using the lavatory and at the end of the working period. Appropriate techniques should be used to remove potentially contaminated clothing. Wash contaminated clothing before reusing. Ensure that eyewash stations and safety showers are close to the workstation location.

Section 9. Physical and Chemical Properties

Physical state : Liquid.

Flash point : Closed cup: 25.556°C (78°F).

Auto-ignition temperature : The lowest known value is 398.85°C (749.9°F) (Ethanol).

Flammable limits : Lower: 3.3% Upper: 19%

Color : Clear.

Odor : Alcohol-like.

Boiling/condensation point : 78°C (172.4°F)

Melting/freezing point : May start to solidify at 16.67°C (62°F) based on data for: Acetic Acid . Weighted average: -103.6°C (-154.5°F)

Relative density : Weighted average: 0.79 (Water = 1)

Vapor pressure : 2.5 kPa (19 mm Hg) (at 20°C)

Vapor density : The highest known value is 2.1 (Air = 1) (Acetic Acid). Weighted average: 1.48 (Air = 1)

Odor threshold : The lowest known value is 1 ppm (Acetic Acid) Weighted average: 92.67 ppm

Evaporation rate : The highest known value is 1.7 (Ethanol) Weighted average: 1.63 compared with (n-Butyl Acetate =1)

Section 10. Stability and Reactivity

Stability and reactivity : The product is stable.

Conditions of instability : Avoid exposure to heat, spark or open flame.

Incompatibility with various substances : Highly reactive or incompatible with the following materials: oxidizing materials. Reactive or incompatible with the following materials: metals, acids and alkalis. Avoid contact with strong oxidizers, excessive heat, sparks or open flame. Reacts with heavy metals and metallic salts.

Hazardous polymerization : Will not occur.

Conditions of reactivity : Highly flammable in the presence of the following materials or conditions: open flames, sparks and static discharge, heat, shocks and mechanical impacts and oxidizing materials.
Highly explosive in the presence of the following materials or conditions: open flames, sparks and static discharge, heat, shocks and mechanical impacts and oxidizing materials.
Vapor may travel a considerable distance to source of ignition and flash back.

Section 11. Toxicological Information

Toxicity data

United States

<u>Product/ingredient name</u>	<u>Test</u>	<u>Result</u>	<u>Route</u>	<u>Species</u>
Ethanol	LD50	7060 mg/kg	Oral	Rat
	LD50	6300 mg/kg	Oral	Rabbit
	LD50	3450 mg/kg	Oral	Mouse
	LDLo	1400 mg/kg	Oral	human
	LDLo	5500 mg/kg	Oral	Dog
Methanol	LD50	5628 mg/kg	Oral	Rat
	LD50	14200 mg/kg	Oral	Rabbit
	LD50	7300 mg/kg	Oral	Mouse
	LD50	15800 mg/kg	Dermal	Rabbit

Continued on Next Page

Section 11. Toxicological Information

	LDLo	143 mg/kg	Oral	human
	LDLo	428 mg/kg	Oral	human
	LDLo	6422 mg/kg	Oral	man
	LDLo	393 mg/kg	Dermal	Monkey.
	LC50	64000 ppm (4 hour/hours)	Inhalation	Rat
Acetic Acid	LD50	3310 mg/kg	Oral	Rat
	LD50	4960 mg/kg	Oral	Mammal
	LD50	1060 mg/kg	Dermal	Mammal
	LDLo	600 mg/kg	Oral	Rabbit
	LDLo	600 mg/kg	Oral	Rabbit
	LC50	5620 ppm (1 hour/hours)	Inhalation	Muskrat

Chronic effects on humans : **CARCINOGENIC EFFECTS** Classified A4 (Not classifiable for humans or animals.) by ACGIH [Ethanol]. Classified 3 (Not classifiable for humans.) by IARC [Eosin Y]. Classified 3 (Not classifiable for humans.) by IARC [Light Green SF, Yellowish]. Contains material which causes damage to the following organs: blood, the reproductive system, liver, gastrointestinal tract, upper respiratory tract, skin, central nervous system (CNS), eye, lens or cornea, teeth.

Other toxic effects on humans : Very hazardous in case of skin contact (irritant, permeator), of eye contact (irritant), of inhalation (lung irritant). Hazardous in case of ingestion.

Special remarks on other toxic effects on humans : Irritations of mucous membranes in the mouth, pharynx, oesophagus and gastrointestinal tract.

Specific effects

- Carcinogenic effects** : No known significant effects or critical hazards.
- Mutagenic effects** : No known significant effects or critical hazards.
- Teratogenicity / Reproductive toxicity** : Contains material which can cause birth defects.

Sensitization

- Ingestion** : No known significant effects or critical hazards.
- Inhalation** : No known significant effects or critical hazards.
- Eyes** : Irritating to eyes.
- Skin** : Irritating to skin.

Section 12. Ecological Information

Ecotoxicity data

United States

Product/ingredient name

<u>Product/ingredient name</u>	<u>Species</u>	<u>Period</u>	<u>Result</u>
Ethanol	Daphnia magna (EC50)	48 hour/hours	2 mg/l
	Daphnia magna (EC50)	48 hour/hours	9.3 mg/l
	Daphnia magna (EC50)	48 hour/hours	>100 mg/l
	Daphnia magna (LC50)	96 hour/hours	>100 mg/l
	Pimephales promelas (LC50)	96 hour/hours	>100 mg/l
	Oncorhynchus mykiss (LC50)	96 hour/hours	13200 mg/l
Methanol	Daphnia magna (EC50)	48 hour/hours	>10000 mg/l
	Oncorhynchus mykiss (EC50)	48 hour/hours	13200 mg/l
	Lepomis macrochirus (EC50)	48 hour/hours	16000 mg/l
	Daphnia magna (LC50)	96 hour/hours	>100 mg/l
	Pimephales promelas (LC50)	96 hour/hours	>100 mg/l
	Lepomis macrochirus (LC50)	96 hour/hours	15400 mg/l
Acetic Acid	Daphnia magna (EC50)	48 hour/hours	65 mg/l
	Lepomis macrochirus (LC50)	96 hour/hours	75 mg/l
	Pimephales promelas (LC50)	96 hour/hours	79 mg/l
	Pimephales promelas (LC50)	96 hour/hours	88 mg/l

Continued on Next Page

Section 12. Ecological Information

- Environmental precautions** : No known significant effects or critical hazards.
- Products of degradation** : These products are carbon oxides (CO, CO₂) and water.
- Toxicity of the products of biodegradation** : The products of degradation are as toxic as the product itself.

Section 13. Disposal Considerations


- Waste disposal** : The generation of waste should be avoided or minimized wherever possible. Avoid dispersal of spilled material and runoff and contact with soil, waterways, drains and sewers. Disposal of this product, solutions and any by-products should at all times comply with the requirements of environmental protection and waste disposal legislation and any regional local authority requirements.
- RCRA classification** : **Code:** (1)

Disposal should be in accordance with applicable regional, national and local laws and regulations. Local regulations may be more stringent than regional or national requirements.

The information presented below only applies to the material as supplied. The identification based on characteristic(s) or listing may not apply if the material has been used or otherwise contaminated. It is the responsibility of the waste generator to determine the toxicity and physical properties of the material generated to determine the proper waste identification and disposal methods in compliance with applicable regulations.

Refer to Section 7: HANDLING AND STORAGE and Section 8: EXPOSURE CONTROLS/PERSONAL PROTECTION for additional handling information and protection of employees.

Section 14. Transport Information

Regulatory information	UN number	Proper shipping name	Class	PG*	Label	Additional information
DOT Classification	UN1987	Alcohols, N.O.S. (Ethanol, Methanol)	3	II		Not available.

PG* : Packing group

Section 15. Regulatory Information

United States

- HCS Classification** : Flammable liquid
Toxic material
Irritating material
Target organ effects
- U.S. Federal regulations** : TSCA 8(b) inventory: Listed

Section 15. Regulatory Information

SARA 302/304/311/312 extremely hazardous substances: No products were found.
 SARA 302/304 emergency planning and notification: No products were found.
 SARA 302/304/311/312 hazardous chemicals: Ethanol; Methanol; Acetic Acid
 SARA 311/312 MSDS distribution - chemical inventory - hazard identification: Ethanol:
 Fire hazard, Immediate (acute) health hazard, Delayed (chronic) health hazard;
 Methanol: Fire hazard, Immediate (acute) health hazard, Delayed (chronic) health
 hazard; Acetic Acid : Fire hazard, Immediate (acute) health hazard, Delayed (chronic)
 health hazard
 Clean Water Act (CWA) 307: No products were found.
 Clean Water Act (CWA) 311: Acetic Acid
 Clean Air Act (CAA) 112 accidental release prevention: No products were found.
 Clean Air Act (CAA) 112 regulated flammable substances: No products were found.
 Clean Air Act (CAA) 112 regulated toxic substances: No products were found.

SARA 313

	<u>Product name</u>	<u>CAS number</u>	<u>Concentration</u>
Form R - Reporting requirements	: Methanol	67-56-1	25
Supplier notification	: Methanol	67-56-1	25

SARA 313 notifications must not be detached from the MSDS and any copying and redistribution of the MSDS shall include copying and redistribution of the notice attached to copies of the MSDS subsequently redistributed.

State regulations : Pennsylvania RTK: Ethanol: (generic environmental hazard); Methanol: (environmental hazard, generic environmental hazard); Acetic Acid : (environmental hazard, generic environmental hazard)
 Massachusetts RTK: Ethanol; Methanol; Acetic Acid
 New Jersey: Gill Modified EA

Canada

WHMIS (Canada) : Class B-2: Flammable liquid
 Class D-1A: Material causing immediate and serious toxic effects (Very toxic).
 Class D-2A: Material causing other toxic effects (Very toxic).

CEPA DSL/CEPA NDSL : CEPA DSL: Ethanol; Methanol; Acetic Acid ; Water

This product has been classified according to the hazard criteria of the CPR and the MSDS contains all the information required by the CPR.

EU regulations

Hazard symbol/symbols :



Risk phrases : R10- Flammable.
 R23/24/25- Toxic by inhalation, in contact with skin and if swallowed.
 R39/23/24/25- Toxic: danger of very serious irreversible effects through inhalation, in contact with skin and if swallowed.

Safety phrases : S1/2- Keep locked up and out of the reach of children.
 S36/37- Wear suitable protective clothing and gloves.
 S45- In case of accident or if you feel unwell, seek medical advice immediately (show the label where possible).
 S63- In case of accident by inhalation: remove casualty to fresh air and keep at rest.

International regulations

International lists : Australia (NICNAS): Ethanol; Methanol; Acetic Acid ; Phosphotungstic Acid; Eosin Y; Light Green SF, Yellowish ; Water

China: Ethanol; Methanol; Acetic Acid ; Eosin Y; Light Green SF, Yellowish

Germany water class: Ethanol; Methanol; Acetic Acid

Section 15. Regulatory Information

Japan (METI): Ethanol; Methanol; Acetic Acid ; Phosphotungstic Acid; Eosin Y; Light Green SF, Yellowish ; Water

Korea (TCCL): Ethanol; Methanol; Acetic Acid ; Eosin Y; Water

Philippines (RA6969): Ethanol; Methanol; Acetic Acid ; Eosin Y; Water

Section 16. Other Information

Label requirements

: DANGER!
POISON!

FLAMMABLE LIQUID AND VAPOR.

VAPOR MAY CAUSE FLASH FIRE.

VAPOR HARMFUL.

MAY BE FATAL OR CAUSE BLINDNESS IF SWALLOWED.

CANNOT BE MADE NONPOISONOUS.

HARMFUL IF INHALED OR ABSORBED THROUGH SKIN.

CAUSES EYE AND SKIN IRRITATION.

BIRTH DEFECT HAZARD.

CONTAINS MATERIAL WHICH CAN CAUSE BIRTH DEFECT.

CONTAINS MATERIAL WHICH MAY CAUSE DAMAGE TO THE FOLLOWING


ORGANS: BLOOD, REPRODUCTIVE SYSTEM, LIVER, GASTROINTESTINAL TRACT,

RESPIRATORY TRACT, SKIN, CENTRAL NERVOUS SYSTEM, EYE, LENS OR

CORNEA, TEETH

National Fire Protection Association (U.S.A.)

:



Notice to reader

The statements contained herein are based upon technical data that EMD Chemicals Inc. believes to be reliable, are offered for information purposes only and as a guide to the appropriate precautionary and emergency handling of the material by a properly trained person having the necessary technical skills. Users should consider these data only as a supplement to other information gathered by them and must make independent determinations of suitability and completeness of information from all sources to assure proper use, storage and disposal of these materials and the safety and health of employees and customers and the protection of the environment. EMD CHEMICALS INC. MAKES NO REPRESENTATION OR WARRANTY OF ANY KIND, EXPRESS OR IMPLIED, INCLUDING MERCHANTABILITY OR FITNESS FOR A PARTICULAR USE, WITH RESPECT TO THE INFORMATION HEREIN OR THE PRODUCT TO WHICH THE INFORMATION REFERS.