

# Certificate of Quality



|                   |  |
|-------------------|--|
| Part Number       | <b>MFLX06417-31</b>                    |
| Lot Number        | <b>0424149811</b>                      |
| Part Description  | <b>TUBING PTFE 0.032X0.056IN 100FT</b> |
| Expiration Date   | <b>02JUN2025</b>                       |
| Country of Origin | <b>USA</b>                             |

## Inspection & Testing

- It is the responsibility of the user to conduct tests that are deemed necessary to determine the suitability of this product for any particular use.

## Product Specifications & Compliance

- Compound: PTFE
- All ingredients used in these products are approved under the FDA Code of Federal Regulations 21, subchapter B Part 177.1550, for use with food and food products.
- The tubing itself has not been tested for approval by the FDA or any other testing body.
- Material is manufactured in an ISO 9001-2015 approved facility.
- All PTFE tubing is manufactured with virgin resin.
- The proper cleaning and certification of this product for food applications is the responsibility of the end-user.
- Material can be sterilized in autoclave.
- Material is RoHS Complaint.

## Storage

- Storage in original packaging, at ambient conditions and away from direct sunlight, is recommended.

## Certificate Approval

- This certificate was issued and approved by the Avantor Fluid Handling Quality and Regulatory Departments and is valid without signature.