

Monoethanolamine
BAKER ANALYZED® A.C.S. Reagent
(2-aminoethanol)



Material No.: 9314-01

Revision No.: 0

Product Specification

Meets ACS Reagent Chemical Requirements

Test	Specification
Assay (HOCH ₂ CH ₂ NH ₂) (by GC, corrected for water)	≥ 99.0 %
Color (APHA)	≤ 15
Water (H ₂ O)	≤ 0.30 %
Heavy Metals (by ICP-OES)	≤ 5 ppm
Trace Impurities - Iron (Fe)	≤ 5 ppm

For Laboratory, Research, or Manufacturing Use
Meets Reagent Specifications for testing USP/NF monographs Product
Information (not specifications): Density (g/mL) @ 25°C (typical)=1.012

Packaging Site: Phillipsburg Mfg Ctr & DC


Jamie Ethier
Vice President Global Quality

For questions on this Product Specification please contact Technical Services at 855.282.6867 or +1.610.386.1700

Avantor Performance Materials, LLC

100 Matsonford Rd, Suite 200, Radnor, PA 19087. U.S.A. Phone 610.386.1700