

Specifications for Physician Seating



	92 OPERATOR STOOL	96 OPERATOR STOOL WITH BACK	101 PNUEMATIC STOOL	106 PNUEMATIC STOOL WITH BACK
Warranty:	4 year	4 year	4 year	4 year
Base Construction:	fully welded, tubular steel base	fully welded, tubular steel base	cast aluminum	cast aluminum
Base Finish:	durable powder coated finish	durable powder coated finish	polished aluminum	polished aluminum
Number of Legs:	5	5	5	5
Foot Ring:	yes	yes	no*	no*
Caster Style:	2" ball bearing, hard surface caster	2" ball bearing, hard surface caster	2" carpet caster	2" carpet caster
Adjustment System:	threaded stem	threaded stem	pneumatic gas cylinder	pneumatic gas cylinder
Adjustment By:	rotation of seat	rotation of seat	hand lever under seat	hand lever under seat
Height:	18.5" to 25"	18.5" to 25"	20" to 30"	20" to 30"
Back Support:	no	yes	no	yes
Back Adjustment:	N/A	3" vertical	N/A	3" vertical
Seat Dimensions:	14"W x 2.5" H	14"W x 2.5" H	14"W x 2.5" H	14"W x 2.5" H
Weight Capacity:	350 lbs.	225 lbs.	350 lbs.	225 lbs.
Stool Weight:	16 lbs.	21 lbs.	13 lbs.	18 lbs.
Shipping:	assembly required	assembly required	assembly required	assembly required
Color Selection:	Seat: Black, Cambridge Blue, Desert Sand, Earthen Brown, English Brown, Emerald, Fawn, Ninja, Ocean Jade, Smoke, Windsor Wine			
<i>*available as an add on accessory</i>				