



# SAFETY DATA SHEET

according to the (US) Hazard Communication Standard (29 CFR 1910.1200)

Date of issue: 01/19/2015

Version 1.0

## SECTION 1. Identification

### Product identifier

Product number	WX0001
Product name	Water LC-MS Grade For Liquid Chromatography -  Mass Spectrometry OmniSolv®
CAS-No.	7732-18-5

### Relevant identified uses of the substance or mixture and uses advised against

Identified uses	Reagent for analysis
-----------------	----------------------

### Details of the supplier of the safety data sheet

Company	EMD Millipore Corporation   290 Concord Road, Billerica, MA 01821, United States of America   General Inquiries: +1-978-715-4321   Monday to Friday, 9:00 AM to 4:00 PM Eastern Time (GMT-5)
---------	--

Emergency telephone	800-424-9300 CHEMTREC (USA) +1-703-527-3887 CHEMTREC (International) 24 Hours/day; 7 Days/week
---------------------	--

## SECTION 2. Hazards identification

### GHS-Labeling

Not a dangerous substance according to GHS.

### Other hazards

None known.

## SECTION 3. Composition/information on ingredients

Formula	H <sub>2</sub> O	H <sub>2</sub> O (Hill)
Molar mass	18.02 g/mol	
Remarks	No hazardous ingredients according to the OSHA Hazard Communication Standard 29 CFR 1910.1200.	

## SECTION 4. First aid measures

### Description of first-aid measures

# SAFETY DATA SHEET

according to the (US) Hazard Communication Standard (29 CFR 1910.1200)

Product number

WX0001

Version 1.0

Product name

Water LC-MS Grade For Liquid Chromatography - <br/>Mass Spectrometry  
OmniSolv®

---

## *General advice*

No hazards which require special first aid measures.

Never give anything by mouth to an unconscious person.

## **Most important symptoms and effects, both acute and delayed**

We have no description of any toxic symptoms.

## **Indication of any immediate medical attention and special treatment needed**

No information available.

---

## **SECTION 5. Fire-fighting measures**

### **Extinguishing media**

#### *Suitable extinguishing media*

Use extinguishing measures that are appropriate to local circumstances and the surrounding environment.

#### *Unsuitable extinguishing media*

For this substance/mixture no limitations of extinguishing agents are given.

### **Special hazards arising from the substance or mixture**

Not combustible.

### **Advice for firefighters**

#### *Special protective equipment for fire-fighters*

none

#### *Further information*

none

---

## **SECTION 6. Accidental release measures**

### **Personal precautions, protective equipment and emergency procedures**

none

### **Environmental precautions**

No special precautionary measures necessary.

### **Methods and materials for containment and cleaning up**

Observe possible material restrictions (see sections 7 and 10).

Pour into drain.

---

## **SECTION 7. Handling and storage**

### **Precautions for safe handling**

Observe label precautions.

### **Conditions for safe storage, including any incompatibilities**

Tightly closed.

Store at room temperature.

---

# SAFETY DATA SHEET

according to the (US) Hazard Communication Standard (29 CFR 1910.1200)

Product number

WX0001

Version 1.0

Product name

Water LC-MS Grade For Liquid Chromatography -   
Mass Spectrometry  
OmniSolv®

---

## SECTION 8. Exposure controls/personal protection

### Exposure limit(s)

Contains no substances with occupational exposure limit values.

### Engineering measures

Technical measures and appropriate working operations should be given priority over the use of personal protective equipment.

### Individual protection measures

Protective clothing should be selected specifically for the workplace, depending on concentration and quantity of the hazardous substances handled. The chemical resistance of the protective equipment should be inquired at the respective supplier.

### *Hygiene measures*

not required

### *Eye/face protection*

not required

### *Hand protection*

Chemical-resistant, impervious gloves complying with an approved standard should be worn at all times when handling chemical products if a risk assessment indicates this is necessary.

### *Respiratory protection*

not required

Use a properly fitted, air-purifying or air-fed respirator complying with an approved standard if a risk assessment indicates this is necessary. Respirator selection must be based on known or anticipated exposure levels, the hazards of the product and the safe working limits of the selected respirator.

---

## SECTION 9. Physical and chemical properties

Physical state	liquid
Color	colorless
Odor	odorless
Odor Threshold	Not applicable
pH	at 68 °F (20 °C) neutral
Melting point	0 °C
Boiling point/boiling range	212 °F (100 °C) at 1,013 hPa
Flash point	Not applicable
Evaporation rate	No information available.
Flammability (solid, gas)	No information available.

---

# SAFETY DATA SHEET

according to the (US) Hazard Communication Standard (29 CFR 1910.1200)

Product number

WX0001

Version 1.0

Product name

Water LC-MS Grade For Liquid Chromatography - <br/>Mass Spectrometry  
OmniSolv®

---

Lower explosion limit	Not applicable
Upper explosion limit	Not applicable
Vapor pressure	23 hPa at 68 °F (20 °C)
Relative vapor density	No information available.
Density	1.00 g/cm <sup>3</sup> at 68 °F (20 °C)
Relative density	No information available.
Water solubility	completely soluble
Partition coefficient: n-octanol/water	Not applicable
Autoignition temperature	Not applicable
Decomposition temperature	Distillable in an undecomposed state at normal pressure.
Viscosity, dynamic	0.952 mPa.s at 68 °F (20 °C)
Explosive properties	Not classified as explosive.
Oxidizing properties	none
Ignition temperature	Not applicable
Minimum ignition energy	Not applicable

---

## SECTION 10. Stability and reactivity

### Reactivity

See below

### Chemical stability

The product is chemically stable under standard ambient conditions (room temperature) .

### Possibility of hazardous reactions

Violent reactions possible with:

The generally known reaction partners of water.

### Conditions to avoid

none

### Incompatible materials

# SAFETY DATA SHEET

according to the (US) Hazard Communication Standard (29 CFR 1910.1200)

Product number

WX0001

Version 1.0

Product name

Water LC-MS Grade For Liquid Chromatography - <br/>Mass Spectrometry  
OmniSolv®

---

no information available

## Hazardous decomposition products

none

---

## SECTION 11. Toxicological information

### Information on toxicological effects

#### *Likely route of exposure*

Inhalation, Eye contact, Skin contact

#### *Specific target organ systemic toxicity - single exposure*

The substance or mixture is not classified as specific target organ toxicant, single exposure.

#### *Specific target organ systemic toxicity - repeated exposure*

The substance or mixture is not classified as specific target organ toxicant, repeated exposure.

#### *Aspiration hazard*

Regarding the available data the classification criteria are not fulfilled.

### Carcinogenicity

IARC

No ingredient of this product present at levels greater than or equal to 0.1% is identified as probable, possible or confirmed human carcinogen by IARC.

OSHA

No ingredient of this product present at levels greater than or equal to 0.1% is identified as a carcinogen or potential carcinogen by OSHA.

NTP

No ingredient of this product present at levels greater than or equal to 0.1% is identified as a known or anticipated carcinogen by NTP.

ACGIH

No ingredient of this product present at levels greater than or equal to 0.1% is identified as a carcinogen or potential carcinogen by ACGIH.

### Further information

No toxic effects are to be expected when the product is handled appropriately.

---

## SECTION 12. Ecological information

### Ecotoxicity

No information available.

### Persistence and degradability

No information available.

### Bioaccumulative potential

*Partition coefficient: n-octanol/water*

Not applicable

### Mobility in soil

---

# SAFETY DATA SHEET

according to the (US) Hazard Communication Standard (29 CFR 1910.1200)

Product number

WX0001

Version 1.0

Product name

Water LC-MS Grade For Liquid Chromatography -   
Mass Spectrometry  
OmniSolv®

---

No information available.

### *Additional ecological information*

No ecological problems are to be expected when the product is handled and used with due care and attention.

---

## **SECTION 13. Disposal considerations**

The information presented only applies to the material as supplied. The identification based on characteristic(s) or listing may not apply if the material has been used or otherwise contaminated. It is the responsibility of the waste generator to determine the toxicity and physical properties of the material generated to determine the proper waste identification and disposal methods in compliance with applicable regulations. Disposal should be in accordance with applicable regional, national and local laws and regulations.

---

## **SECTION 14. Transport information**

### **Land transport (DOT)**

Not classified as dangerous in the meaning of transport regulations.

### **Air transport (IATA)**

Not classified as dangerous in the meaning of transport regulations.

### **Sea transport (IMDG)**

Not classified as dangerous in the meaning of transport regulations.

---

## **SECTION 15. Regulatory information**

### **United States of America**

#### **SARA 313**

This material does not contain any chemical components with known CAS numbers that exceed the threshold (De Minimis) reporting levels established by SARA Title III, Section 313.

#### **SARA 302**

No chemicals in this material are subject to the reporting requirements of SARA Title III, Section 302.

#### **Clean Water Act**

This product does not contain any Hazardous Substances listed under the U.S. CleanWater Act, Section 311, Table 116.4A.

This product does not contain any Hazardous Chemicals listed under the U.S. CleanWater Act, Section 311, Table 117.3.

#### **DEA List I**

Not listed

#### **DEA List II**

Not listed

### **US State Regulations**

#### **Massachusetts Right To Know**

---

# SAFETY DATA SHEET

according to the (US) Hazard Communication Standard (29 CFR 1910.1200)

Product number

WX0001

Version 1.0

Product name

Water LC-MS Grade For Liquid Chromatography -   
Mass Spectrometry  
OmniSolv®

---

## Remarks

No components are subject to the Massachusetts Right to Know Act.

### California Prop 65 Components

This product does not contain any chemicals known to the State of California to cause cancer, birth, or any other reproductive defects.

### Notification status

TSCA:

All components of the product are listed in the TSCA-inventory.

DSL:

All components of this product are on the Canadian DSL.

---

## SECTION 16. Other information

### Training advice

Provide adequate information, instruction and training for operators.

### Key or legend to abbreviations and acronyms used in the safety data sheet

Used abbreviations and acronyms can be looked up at [www.wikipedia.org](http://www.wikipedia.org).

Date of issue:01/19/2015

---

The information contained herein is based on the present state of our knowledge. It characterizes the product with regard to appropriate safety precautions. It does not represent a warranty of any product properties and we assume no liability for any loss or injury which may result from the use of this information. Users should conduct their own investigations to determine the suitability of the information.

*All rights reserved. Millipore and the "M" Mark are registered trademarks of Merck KGaA, Darmstadt, Germany.*