

DURAN®

YOUTILITY BOTTLES

Glass laboratory bottles for use in a wide and diverse range of scientific research applications. The innovation does not stop at the improved ergonomic bottle shape, but extends to dedicated accessories that incorporate many new innovative features to improve handling, sample identification, and ease of use. The new YOUTILITY bottle system helps to make laboratory work easier, safer, more economical and fun.

Bottle Features:

- Borosilicate 3.3, a pharmacopoeial Type I neutral glass
- Offered in amber color. *(Fully meets the requirements of the USP<660> and Ph. Eur. for UV light protection)*
- Integrated Retrace Code for full product traceability of the borosilicate glass bottle
- Flat bottoms for use with magnetic stirrers and stir bars
- Gripping zones permit easier and safer handling
- Very good chemical resistance and high temperature resistance
- The volume is clearly and indicated in durable white ceramic
- Slimmer bottle shape allows a more optimal use of limited space in autoclaves and laboratory refrigerators

Cap Features:

- Shortened thread permits faster opening and closing, yet it is fully compatible with standard GL45 threads
- Ergonomic shape has comfortable grooves and ridges for a more efficient and easier tightening or removal, especially with smaller sized or gloved hands
- Cap plug sealing system ensures a liquid tight seal
- Manufactured from a food grade polypropylene that is fully autoclavable and totally reusable

Also Available...

Contact us for other bottles and accessories.

5563 SERIES



VWR™
part of avantor

order on
VWR.COM

 **ACE GLASS**
INCORPORATED

DURAN®

YOUTILITY BOTTLES

The accurate and reliable labelling of scientific samples is an integral part of the success of your laboratory. Labels carry critical information such as the description of the contents, composition, concentrations, pH, make-up date, lot number, expiration date, storage conditions, user information, and without them, your research could be unreliable and possibly invalid. Labels are the foundation of an effective laboratory bottle management process.

Bottle Tag:

- Easily and securely attached around the neck of the YOUTILITY bottle, even with the screw cap in place
- Tab that makes them easy to remove with bare or gloved hands
- Available in eight colors: red, orange, yellow, green, blue, purple, black and white
- Made of silicone: chemically robust, heat resistant, fully autoclavable and totally re-usable

Bottle & Cap Labels:

- Stay on when you need them to do so, and are easy to remove when no longer required
- Manufactured from a tear resistant, matt white, polyester material coated with a high performance adhesive
- Chemically resistant to typical laboratory chemicals such as disinfectants and many solvents
- -40 °C to +150 °C temperature performance range, making them suitable for use in fridges, freezers, autoclaves, incubators and water baths without the risk of the labels drying out or falling off

Also Available...

Contact us for other bottles, caps, and accessories.

LABELLING



vwr[™]
part of avantor

1.800.932.5000 | VWR.COM

Prices and product details are current when published and subject to change without notice. | Certain products may be limited by federal, state, provincial, or local regulations. | VWR, part of Avantor, makes no claims or warranties concerning sustainable/green products. Any claims concerning sustainable/green products are the sole claims of the manufacturer and not those of VWR International, LLC and/or Avantor, Inc. or affiliates. All prices are in US dollars unless otherwise noted. Offers valid in US, void where prohibited by law or company policy, while supplies last. | Trademarks are owned by Avantor, Inc. or its affiliates, unless otherwise noted. | Visit vwr.com to view our privacy policy, trademark owners, and additional disclaimers. © 2019 Avantor, Inc. All rights reserved.