

Application comparison table

VWR® LIFE SCIENCE POWER SUPPLIES



VWR Life Science Power Supply	VWR Mini	VWR 500V	VWR 300V HC
Description	Basic and economical unit perfect for small nucleic acid gel tanks.	Ideal for basic electrophoresis applications such as nucleic acid, protein separation and small Wet blot Western blotting applications.	
Specifications	100V or 200V	500V / 400mA	300V / 2000mA High Current
Cat. No.	76196-450	76196-454	76196-458
Horizontal Mini (up to 7 x 10cm)	●	●	●
Horizontal Midi Gels (up to 10 x 10cm)	●	●	●
Horizontal Maxi Gels (up to 20 x 20cm)		●	●
Mini PAGE Gels		●	●
Mini PAGE Gels, two or more tanks		●	●
Mini Wet Blotting 8 x 10cm, 10 x 10cm			●
Mini Wet Blotting, two or more tanks			●
Maxi Wet Blotting 20 x 20cm			●

115 VAC units contain a standard US plug and cord.
230 VAC units contain an international plug and cord.