

Certificate of Quality

VWR Catalogue Number	See List
Description	VWR® Red Vacuum Rubber Tubing
Country of Origin	Manufactured in USA
Date of Issue (yyyy-mm-dd)	2022-01-21

VWR Catalogue Number	Coil length	Inner Dimension x Outer Dimension	Wall Thickness
62995-335	9 ³ / ₄ '	¹ / ₄ x ⁵ / ₈ "	³ / ₁₆ "
62995-368	9 ³ / ₄ '	¹ / ₂ x 1 ¹ / ₈ "	⁵ / ₁₆ "
62995-390	9 ³ / ₄ '	⁵ / ₈ x 1 ³ / ₈ "	³ / ₈ "
62995-426	9 ³ / ₄ '	³ / ₄ x 1 ¹ / ₂ "	³ / ₈ "
62995-059	50'	¹ / ₄ x ¹ / ₂ "	¹ / ₈ "
62995-062	50'	¹ / ₄ x ⁵ / ₈ "	³ / ₁₆ "
62995-095	50'	¹ / ₂ x 1 ¹ / ₈ "	⁵ / ₁₆ "
62995-120	50'	⁵ / ₈ x 1 ³ / ₈ "	³ / ₈ "
62995-153	50'	³ / ₄ x 1 ¹ / ₂ "	³ / ₈ "

Quality System Compliance

Products are manufactured under the [ISO-13485:2016](#) standard. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable product specifications.

Product Specifications

Material	Dry Natural Rubber.
Sterilization	Not applicable.
ATP Assay	Not applicable.
DNase & RNase Free	Not applicable.

BSE/TSE	Not applicable.
Cytotoxicity	Not applicable.
Non-Pyrogenic Statement	Not applicable.
Latex Statement	This product contains dry natural rubber.
BPA Statement	Not applicable.
DEHP Statement	Not applicable.
RoHS	Not applicable.
REACH Statement	Not applicable.
Storage Conditions	Not applicable.
Shelf Life	Not applicable.

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.
This document has been produced electronically and is valid without a signature.