



**Item Number** PR130  
**Item** Propylene Glycol, USP  
**CAS Number** 57-55-6  
**Molecular Formula** C<sub>3</sub>H<sub>8</sub>O<sub>2</sub>  
**Molecular Weight** 76.10

Assay..... Min. 99.5%  
Specific Gravity @ 25°C..... 1.035-1.037

**MAXIMUM LIMITS**

Acidity..... To pass test  
Water..... 0.2%  
Residue on Ignition..... 0.007%  
Chloride (Cl)..... 0.007%  
Sulfate (SO<sub>4</sub>)..... 0.006%  
Heavy Metals (as Pb)..... 5 ppm  
Residual solvents..... To pass test

***Spectrum Chemicals & Laboratory Products***

***Corporate Headquarters: East Coast Facility:***

14422 S. San Pedro St. 755 Jersey Ave.

Gardena, CA 90248 New Brunswick, NJ 08901

(800) 813-1514 • (310) 516-8000 • Fax (310) 516-9843