



**Item Number** CH123  
**Item** Chlorobutanol, Anhydrous, NF  
**CAS Number** 57-15-8  
**Molecular Formula** C<sub>4</sub>H<sub>7</sub>Cl<sub>3</sub>O  
**Molecular Weight** 177.46

Assay (Anhydrous Basis).....98.0-100.5%  
Reaction (to Litmus)..... Neutral

MAXIMUM LIMITS

Water..... 1.0%  
Chloride..... 0.07%  
Residual Solvents..... To pass test

***Spectrum Chemicals & Laboratory Products***

---

**Corporate Headquarters:**      **East Coast Facility:**  
14422 S. San Pedro St.      755 Jersey Ave.  
Gardena, CA 90248      New Brunswick, NJ 08901  
(800) 813-1514 • (310) 516-8000 • Fax (310) 516-9843